

Beckman Coulter Systems Line


New

SARS-CoV-2 SEROLOGICAL ASSAY

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
----------	-----------	---------	--------	--------------	-----------	----------------------

Please contact us for more information

SARS-CoV-2 SEROLOGICAL ASSAY based on PECIA Colloidal Gold



ALBUMIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	ALTUR-B00/AU	1x50 ml	1x15 ml	417	36 months	Serum & plasma	1,096 - 120.44 (g/l)
	ALTUR-H00/AU	2x50 ml	2x15 ml	833			
	ALTUR-L00/AU	12x50 ml	12x15 ml	5000			
KIT for AU®400-640-2700-5400	ALTUR-B00/AU	1x50 ml	1x15 ml	222	36 months	Serum & plasma	1,096 - 120.44 (g/l)
	ALTUR-H00/AU	2x50 ml	2x15 ml	444			
	ALTUR-L00/AU	12x50 ml	12x15 ml	2667			
KIT for AU®5800	ALTUR-B00/AU	1x50 ml	1x15 ml	510	36 months	Serum & plasma	1,096 - 120.44 (g/l)
	ALTUR-H00/AU	2x50 ml	2x15 ml	1020			
	ALTUR-L00/AU	12x50 ml	12x15 ml	6122			
KIT for DxC 700 AU®	ALTUR-B00/AU	1x50 ml	1x15 ml	510	36 months	Serum & plasma	1,096 - 120.44 (g/l)
	ALTUR-H00/AU	2x50 ml	2x15 ml	1020			
	ALTUR-L00/AU	12x50 ml	12x15 ml	6122			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALBUMIN CSF							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	MATUR-B00/AU	1x50 ml	1x10 ml	417	36 months	CSF	10,86 - 494,33 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	833			
	MATUR-L00/AU	12x50 ml	12x10 ml	5000			
KIT for AU®400-640-2700-5400	MATUR-B00/AU	1x50 ml	1x10 ml	333	36 months	CSF	10,86 - 494,33 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	667			
	MATUR-L00/AU	12x50 ml	12x10 ml	4000			
KIT for AU®5800	MATUR-B00/AU	1x50 ml	1x10 ml	556	36 months	CSF	10,86 - 494,33 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	1111			
	MATUR-L00/AU	12x50 ml	12x10 ml	6667			
KIT for DxC 700 AU®	MATUR-B00/AU	1x50 ml	1x10 ml	417	36 months	CSF	10,86 - 494,33 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	833			
	MATUR-L00/AU	12x50 ml	12x10 ml	5000			

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MACOZ-002	1x2 ml	24 months
H	MACOX-002	1x2 ml	24 months
M	MACON-002	1x2 ml	24 months
L	MACOS-002	1x2 ml	24 months

ALPHA-1 GLYCOPROTEIN ACIDE/OROSOMUCOID

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	GPTUR-B00/AU	1x60 ml	1x12 ml	300	36 months	Serum	0,143 - 2,42 (g/l)
	GPTUR-H00/AU	2x60 ml	2x12 ml	600			
KIT for AU®400-640-2700-5400	GPTUR-B00/AU	1x60 ml	1x12 ml	240	36 months	Serum	0,143 - 2,42 (g/l)
	GPTUR-H00/AU	2x60 ml	2x12 ml	480			
KIT for AU®5800	GPTUR-B00/AU	1x60 ml	1x12 ml	400	36 months	Serum	0,143 - 2,42 (g/l)
	GPTUR-H00/AU	2x60 ml	2x12 ml	800			
KIT for DxC 700 AU®	GPTUR-B00/AU	1x60 ml	1x12 ml	400	36 months	Serum	0,143 - 2,42 (g/l)
	GPTUR-H00/AU	2x60 ml	2x12 ml	800			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA-1 ANTITRYPSIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	ATTUR-H00/AU	1x60 ml	1x12 ml	300	36 months	Serum & plasma	0.007 - 3.50 (g/l)
KIT for AU®400-640-2700-5400	ATTUR-H00/AU	1x60 ml	1x12 ml	240	36 months	Serum & plasma	0.007 - 3.50 (g/l)
KIT for AU®5800	ATTUR-H00/AU	1x60 ml	1x12 ml	400	36 months	Serum & plasma	0.007 - 3.50 (g/l)
KIT for DxC 700 AU®	ATTUR-H00/AU	1x60 ml	1x12 ml	400	36 months	Serum & plasma	0.007 - 3.50 (g/l)

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA-1 MICROGLOBULIN (URINE)

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®400-480-640-680-2700-5400-5800	A1DGX-B00/AU	1 x 15 ml	1 x 15 ml	250	17 months	Serum & plasma	1,29 - 166,6 (mg/l)
	A1DGX-H00/AU	2 x 20 ml	2 x 20 ml	667			
	A1DGX-L00/AU	12 x 20 ml	12 x 20 ml	4000			
KIT for DxC 700 AU®	A1DGX-B00/AU	1 x 15 ml	1 x 15 ml	200	17 months	Serum & plasma	1,29 - 166,6 (mg/l)
	A1DGX-H00/AU	2 x 20 ml	2 x 20 ml	400			

CALIBRATORS

KIT 5 LEV	A1REK-000	5x1 ml	14 months
-----------	-----------	--------	-----------

CONTROLS

M	A1CON-002	1x2 ml	14 months
L	A1COS-002	1x2 ml	14 months

ALPHA-2 MACROGLOBULIN - ALPHA-2 MACROGLOBULIN URINARY

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	MGTUR-H00/AU	1x60 ml	1x12 ml	300	36 months	Serum & plasma or urine	0,026 - 4,96 (g/l) S/P 3,36 - 199 (mg/l) urine
	MGTUR-L00/AU	12x60 ml	12x12 ml	3600			
KIT for AU®400-640-2700-5400	MGTUR-H00/AU	1x60 ml	1x12 ml	240	36 months	Serum & plasma or urine	0,026 - 4,96 (g/l) S/P 3,36 - 199 (mg/l) urine
	MGTUR-L00/AU	12x60 ml	12x12 ml	2880			
KIT for AU®5800	MGTUR-H00/AU	1x60 ml	1x12 ml	400	36 months	Serum & plasma or urine	0,026 - 4,96 (g/l) S/P 3,36 - 199 (mg/l) urine
	MGTUR-L00/AU	12x60 ml	12x12 ml	4800			
KIT for DxC 700 AU®	MGTUR-H00/AU	1x60 ml	1x12 ml	400	36 months	Serum & plasma or urine	0,026 - 4,96 (g/l) S/P 3,36 - 199 (mg/l) urine
	MGTUR-L00/AU	12x60 ml	12x12 ml	4800			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months
URINARY LOW	MGCU1-002	1x2 ml	24 months
URINARY HIGH	MGCU2-002	1x2 ml	24 months

ANTI STREPTOLYSIN O (ASO)

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	ASLTX-B00/AU	1x60 ml	1x6.5 ml	295	20 months	Serum & plasma	19.99 to 864.85 (IU/ml)
	ASLTX-H00/AU	2x60 ml	2x6.5 ml	591			
KIT for AU®400-640-2700-5400	ASLTX-B00/AU	1x60 ml	1x6.5 ml	241	20 months	Serum & plasma	19.99 to 864.85 (IU/ml)
	ASLTX-H00/AU	2x60 ml	2x6.5 ml	481			
KIT for AU®5800	ASLTX-B00/AU	1x60 ml	1x6.5 ml	406	20 months	Serum & plasma	19.99 to 864.85 (IU/ml)
	ASLTX-H00/AU	2x60 ml	2x6.5 ml	813			
KIT for DxC 700 AU®	ASLTX-B00/AU	1x60 ml	1x6.5 ml	406	20 months	Serum & plasma	19.99 to 864.85 (IU/ml)
	ASLTX-H00/AU	2x60 ml	2x6.5 ml	813			

CALIBRATORS

H	ASREK-000	1x1 or (1x5) ml	18 months
---	-----------	-----------------	-----------

CONTROLS

M	ASCON-002	1x2 ml	20 months
---	-----------	--------	-----------

APOLIPOPROTEIN A1

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	AATUR-B00/AU	1x50 ml	1x13 ml	245	36 months	Serum & plasma	0,029 - 4 (g/l)
	AATUR-H00/AU	2x50 ml	2x13 ml	491			
KIT for AU®400-640-2700-5400	AATUR-B00/AU	1x50 ml	1x13 ml	197	36 months	Serum & plasma	0,029 - 4 (g/l)
	AATUR-H00/AU	2x50 ml	2x13 ml	394			
KIT for AU®5800	AATUR-B00/AU	1x50 ml	1x13 ml	325	36 months	Serum & plasma	0,029 - 4 (g/l)
	AATUR-H00/AU	2x50 ml	2x13 ml	650			
KIT for DxC 700 AU®	AATUR-B00/AU	1x50 ml	1x13 ml	325	36 months	Serum & plasma	0,029 - 4 (g/l)
	AATUR-H00/AU	2x50 ml	2x13 ml	650			

CALIBRATORS

KIT 5 LEV	ABREK-000	5x0,5 ml	18 months
-----------	-----------	----------	-----------

CONTROLS

H	ABCOX-002	1x2 ml	18 months
M	ABCON-002	1x2 ml	18 months

APOLIPOPROTEIN B

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	ABTUR-B00/AU	1x50 ml	1x10 ml	286	36 months	Serum & plasma	0.019 - 4 (g/l)
	ABTUR-H00/AU	2x50 ml	2x10 ml	571			
KIT for AU®400-640-2700-5400	ABTUR-B00/AU	1x50 ml	1x10 ml	200	36 months	Serum & plasma	0.019 - 4 (g/l)
	ABTUR-H00/AU	2x50 ml	2x10 ml	400			
KIT for AU®5800	ABTUR-B00/AU	1x50 ml	1x10 ml	333	36 months	Serum & plasma	0.019 - 4 (g/l)
	ABTUR-H00/AU	2x50 ml	2x10 ml	667			
KIT for DxC 700 AU®	ABTUR-B00/AU	1x50 ml	1x10 ml	333	36 months	Serum & plasma	0.019 - 4 (g/l)
	ABTUR-H00/AU	2x50 ml	2x10 ml	667			

CALIBRATORS

KIT 5 LEV	ABREK-000	5x0,5 ml	18 months
-----------	-----------	----------	-----------

CONTROLS

H	ABCOX-002	1x2 ml	18 months
M	ABCON-002	1x2 ml	18 months

BETA-2 MICROGLOBULIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	B2DGX-B00/AU	1x18 ml	1x8 ml	125	30 months	serum & plasma	0,132 to 95 (mg/l)
	B2DGX-H00/AU	2x44 ml	2x19 ml	594			
	B2DGX-L00/AU	8x44 ml	8x19 ml	2375			
KIT for AU®400-640-2700-5400	B2DGX-B00/AU	1x18 ml	1x8 ml	100	30 months	serum & plasma	0,132 to 95 (mg/l)
	B2DGX-H00/AU	2x44 ml	2x19 ml	475			
	B2DGX-L00/AU	8x44 ml	8x19 ml	1900			
KIT for AU®5800	B2DGX-B00/AU	1x18 ml	1x8 ml	167	30 months	serum & plasma	0,132 to 95 (mg/l)
	B2DGX-H00/AU	2x44 ml	2x19 ml	792			
	B2DGX-L00/AU	8x44 ml	8x19 ml	3167			
KIT for DxC 700 AU®	B2DGX-B00/AU	1x18 ml	1x8 ml	167	30 months	serum & plasma	0,132 to 95 (mg/l)
	B2DGX-H00/AU	2x44 ml	2x19 ml	792			
	B2DGX-L00/AU	8x44 ml	8x19 ml	3167			

CALIBRATORS

KIT 5 LEV	B2REK-000	5x1 ml	30 months
-----------	-----------	--------	-----------

CONTROLS

H	B2COX-001	1 x 1 ml	30 months
M	B2CON-003	1x3 ml	30 months
L	B2COS-003	1x3 ml	30 months

BETA-2 MICROGLOBULIN (URINE)

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	BUDGX-B00/AU	1 x 28 ml	1 x 12 ml	125	30 months	Urine	0,055 to 9.3 (mg/l)
	BUDGX-H00/AU	2 x 28 ml	2 x 12 ml	594			
KIT for AU®400-640-2700-5400	BUDGX-B00/AU	1 x 28 ml	1 x 12 ml	100	30 months	Urine	0,055 to 9.3 (mg/l)
	BUDGX-H00/AU	2 x 28 ml	2 x 12 ml	475			
KIT for AU®5800	BUDGX-B00/AU	1 x 28 ml	1 x 12 ml	167	30 months	Urine	0,055 to 9.3 (mg/l)
	BUDGX-H00/AU	2 x 28 ml	2 x 12 ml	792			
KIT for DxC 700 AU®	BUDGX-B00/AU	1 x 28 ml	1 x 12 ml	300	30 months	Urine	0,055 to 9.3 (mg/l)
	BUDGX-H00/AU	2 x 28 ml	2 x 12 ml	600			

CALIBRATORS

KIT 5 LEV	BUREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

M	B2CON-003	1x3 ml	30 months
L	B2COS-003	1x3 ml	30 months

CERULOPLASMIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	CETUR-B00/AU	1x30 ml	1x6 ml	250	36 months	Serum & plasma	9,48 - 1051,2 (mg/l)
	CETUR-H00/AU	1x60 ml	1x12 ml	500			
KIT for AU®400-640-2700-5400	CETUR-B00/AU	1x30 ml	1x6 ml	200	36 months	Serum & plasma	9,48 - 1051,2 (mg/l)
	CETUR-H00/AU	1x60 ml	1x12 ml	400			
KIT for AU®5800	CETUR-B00/AU	1x30 ml	1x6 ml	333	36 months	Serum & plasma	9,48 - 1051,2 (mg/l)
	CETUR-H00/AU	1x60 ml	1x12 ml	667			
KIT for DxC 700 AU®	CETUR-B00/AU	1x30 ml	1x6 ml	250	36 months	Serum & plasma	9,48 - 1051,2 (mg/l)
	CETUR-H00/AU	1x60 ml	1x12 ml	500			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

COMPLEMENT C3

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	C3TUR-B00/AU	1x50 ml	1x10 ml	400	36 months	Serum & plasma	0.015 - 4.31 (g/l)
	C3TUR-H00/AU	2x50 ml	2x10 ml	800			
KIT for AU®400-640-2700-5400	C3TUR-B00/AU	1x50 ml	1x10 ml	263	36 months	Serum & plasma	0.015 - 4.31 (g/l)
	C3TUR-H00/AU	2x50 ml	2x10 ml	526			
KIT for AU®5800	C3TUR-B00/AU	1x50 ml	1x10 ml	500	36 months	Serum & plasma	0.015 - 4.31 (g/l)
	C3TUR-H00/AU	2x50 ml	2x10 ml	1000			
KIT for DxC 700 AU®	C3TUR-B00/AU	1x50 ml	1x10 ml	400	36 months	Serum & plasma	0.015 - 4.31 (g/l)
	C3TUR-H00/AU	2x50 ml	2x10 ml	800			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

COMPLEMENT C4

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	C4TUR-B00/AU	1x50 ml	1x10 ml	400	36 months	Serum & plasma	0,002 - 0,825 (g/l)
	C4TUR-H00/AU	2x50 ml	2x10 ml	800			
KIT for AU®400-640-2700-5400	C4TUR-B00/AU	1x50 ml	1x10 ml	278	36 months	Serum & plasma	0,002 - 0,825 (g/l)
	C4TUR-H00/AU	2x50 ml	2x10 ml	556			
KIT for AU®5800	C4TUR-B00/AU	1x50 ml	1x10 ml	500	36 months	Serum & plasma	0,002 - 0,825 (g/l)
	C4TUR-H00/AU	2x50 ml	2x10 ml	1000			
KIT for DxC 700 AU®	C4TUR-B00/AU	1x50 ml	1x10 ml	400	36 months	Serum & plasma	0,002 - 0,825 (g/l)
	C4TUR-H00/AU	2x50 ml	2x10 ml	800			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CRP XL WIDE RANGE

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	CPTXL-B00/AU	1x60 ml	1x12 ml	120	36 months	Serum & plasma	0,607 - 425 (mg/l)
	CPTXL-H00/AU	2x60 ml	2x12 ml	923			
	CPTXL-L00/AU	12x60 ml	12x12 ml	5938			
KIT for AU®400-640-2700-5400	CPTXL-B00/AU	1x60 ml	1x12 ml	333	36 months	Serum & plasma	0,607 - 425 (mg/l)
	CPTXL-H00/AU	2x60 ml	2x12 ml	667			
	CPTXL-L00/AU	12x60 ml	12x12 ml	4000			
KIT for AU®5800	CPTXL-B00/AU	1x60 ml	1x12 ml	556	36 months	Serum & plasma	0,607 - 425 (mg/l)
	CPTXL-H00/AU	2x60 ml	2x12 ml	1111			
	CPTXL-L00/AU	12x60 ml	12x12 ml	6667			
KIT for DxC 700 AU®	CPTXL-B00/AU	1x60 ml	1x12 ml	462	36 months	Serum & plasma	0,607 - 425 (mg/l)
	CPTXL-H00/AU	2x60 ml	2x12 ml	923			
	CPTXL-L00/AU	12x60 ml	12x12 ml	5538			

CALIBRATORS

KIT 5 LEV	CPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	24 months
M	CPCON-002	1x2 ml	24 months
L	CXCON-002	1x2 ml	24 months

CRP GOLD LATEX

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	CPDGX-B00/AU	1x20 ml	1x20 ml	200	24 months	Serum & plasma	0,031 to 300 (mg/l)
	CPDGX-H00/AU	2x50 ml	2x50 ml	1000			
	CPDGX-L00/AU	8x50 ml	8x50 ml	4000			
KIT for AU®400-640-2700-5400	CPDGX-B00/AU	1x20 ml	1x20 ml	200	24 months	Serum & plasma	0,031 to 300 (mg/l)
	CPDGX-H00/AU	2x50 ml	2x50 ml	1000			
	CPDGX-L00/AU	8x50 ml	8x50 ml	4000			
KIT for AU®5800	CPDGX-B00/AU	1x20 ml	1x20 ml	333	24 months	Serum & plasma	0,031 to 300 (mg/l)
	CPDGX-H00/AU	2x50 ml	2x50 ml	1667			
	CPDGX-L00/AU	8x50 ml	8x50 ml	6667			
KIT for DxC 700 AU®	CPDGX-B00/AU	1x20 ml	1x20 ml	333	24 months	Serum & plasma	0,031 to 300 (mg/l)
	CPDGX-H00/AU	2x50 ml	2x50 ml	1667			
	CPDGX-L00/AU	8x50 ml	8x50 ml	6667			

CALIBRATORS

KIT 5 LEV	CPRLK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	24 months
M	CPCON-002	1x2 ml	24 months
L	CXCON-002	1x2 ml	24 months

CYSTATIN C - Colloidal Gold

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	CYCOL-B00/AU	1x20 ml	1x7.5 ml	125	24 months	Serum & plasma	0.15 to 9.355 (mg/l)
	CYCOL-H00/AU	2x50 ml	2x15 ml	500			
	CYCOL-L00/AU	12x50 ml	12x15 ml	3000			
KIT for AU®400-640-2700-5400	CYCOL-B00/AU	1x20 ml	1x7.5 ml	100	24 months	Serum & plasma	0.15 to 9.355 (mg/l)
	CYCOL-H00/AU	2x50 ml	2x15 ml	400			
	CYCOL-L00/AU	12x50 ml	12x15 ml	2400			
KIT for AU®5800	CYCOL-B00/AU	1x20 ml	1x7.5 ml	167	24 months	Serum & plasma	0.15 to 9.355 (mg/l)
	CYCOL-H00/AU	2x50 ml	2x15 ml	667			
	CYCOL-L00/AU	12x50 ml	12x15 ml	4000			
KIT for DxC 700 AU®	CYCOL-B00/AU	1x20 ml	1x7.5 ml	167	24 months	Serum & plasma	0.15 to 9.355 (mg/l)
	CYCOL-H00/AU	2x50 ml	2x15 ml	667			
	CYCOL-L00/AU	12x50 ml	12x15 ml	4000			

CALIBRATORS

KIT 6 LEV	CYREK-000	6x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	CYCOX-002	1x2 ml	24 months
M	CYCON-002	1x2 ml	24 months
L	CYCOS-002	1x2 ml	24 months

FERRITIN - Colloidal Gold

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	FNCOL-B00/AU	1x50 ml	1x9 ml	360	12 months	Serum	6,344 - 500 (µg/l)
	FNCOL-H00/AU	2x50 ml	2x9 ml	720			
KIT for AU®400-640-2700-5400	FNCOL-B00/AU	1x50 ml	1x9 ml	300	12 months	Serum	6,344 - 500 (µg/l)
	FNCOL-H00/AU	2x50 ml	2x9 ml	600			
KIT for AU®5800	FNCOL-B00/AU	1x50 ml	1x9 ml	450	12 months	Serum	6,344 - 500 (µg/l)
	FNCOL-H00/AU	2x50 ml	2x9 ml	900			
KIT for DxC 700 AU®	FNCOL-B00/AU	1x50 ml	1x9 ml	450	12 months	Serum	6,344 - 500 (µg/l)
	FNCOL-H00/AU	2x50 ml	2x9 ml	900			

CALIBRATORS

KIT 5 LEV	FNRGK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

M	FNCON-002	1x2 ml	18 months
---	-----------	--------	-----------

HAPTOGLOBIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	HATUR-B00/AU	1x30 ml	1x6 ml	150	36 months	Serum & plasma	0.002 - 3.86 (g/l)
	HATUR-H00/AU	1x60 ml	1x12 ml	300			
KIT for AU®400-640-2700-5400	HATUR-B00/AU	1x30 ml	1x6 ml	120	36 months	Serum & plasma	0.002 - 3.86 (g/l)
	HATUR-H00/AU	1x60 ml	1x12 ml	240			
KIT for AU®5800	HATUR-B00/AU	1x30 ml	1x6 ml	200	36 months	Serum & plasma	0.002 - 3.86 (g/l)
	HATUR-H00/AU	1x60 ml	1x12 ml	400			
KIT for DxC 700 AU®	HATUR-B00/AU	1x30 ml	1x6 ml	200	36 months	Serum & plasma	0.002 - 3.86 (g/l)
	HATUR-H00/AU	1x60 ml	1x12 ml	400			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS A

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	IATUR-B00/AU	1x50 ml	1x10 ml	370	36 months	Serum	0.02 - 8 (g/l)
	IATUR-H00/AU	2x50 ml	2x10 ml	741			
KIT for AU®400-640-2700-5400	IATUR-B00/AU	1x50 ml	1x10 ml	299	36 months	Serum	0.02 - 8 (g/l)
	IATUR-H00/AU	2x50 ml	2x10 ml	599			
KIT for AU®5800	IATUR-B00/AU	1x50 ml	1x10 ml	500	36 months	Serum	0.02 - 8 (g/l)
	IATUR-H00/AU	2x50 ml	2x10 ml	1000			
KIT for DxC 700 AU®	IATUR-B00/AU	1x50 ml	1x10 ml	370	36 months	Serum	0.02 - 8 (g/l)
	IATUR-H00/AU	2x50 ml	2x10 ml	741			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS E

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	IEDGX-B00/AU	1x20 ml	1x10 ml	217	20 months	Serum & plasma	21,0 - 1436 UI/ml
	IEDGX-H00/AU	2x30 ml	2x15 ml	652			
KIT for AU®400-640-2700-5400	IEDGX-B00/AU	1x20 ml	1x10 ml	167	20 months	Serum & plasma	21,0 - 1436 UI/ml
	IEDGX-H00/AU	2x30 ml	2x15 ml	500			
KIT for AU®5800	IEDGX-B00/AU	1x20 ml	1x10 ml	250	20 months	Serum & plasma	21,0 - 1436 UI/ml
	IEDGX-H00/AU	2x30 ml	2x15 ml	750			
KIT for DxC 700 AU®	IEDGX-B00/AU	1x20 ml	1x10 ml	250	20 months	Serum & plasma	21,0 - 1436 UI/ml
	IEDGX-H00/AU	2x30 ml	2x15 ml	750			

CALIBRATORS

KIT 5 LEV	IEREK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

M	IECON-003	1x3 ml	18 months
L	IECOS-003	1x3 ml	18 months

IMMUNOGLOBULINS G

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	IGTUR-B00/AU	1x50 ml	1x18 ml	353	36 months	Serum	0.125 - 35 (g/l)
	IGTUR-H00/AU	2x50 ml	2x18 ml	706			
KIT for AU®400-640-2700-5400	IGTUR-B00/AU	1x50 ml	1x18 ml	286	36 months	Serum	0.125 - 35 (g/l)
	IGTUR-H00/AU	2x50 ml	2x18 ml	571			
KIT for AU®5800	IGTUR-B00/AU	1x50 ml	1x18 ml	474	36 months	Serum	0.125 - 35 (g/l)
	IGTUR-H00/AU	2x50 ml	2x18 ml	947			
KIT for DxC 700 AU®	IGTUR-B00/AU	1x50 ml	1x18 ml	353	36 months	Serum	0.125 - 35 (g/l)
	IGTUR-H00/AU	2x50 ml	2x18 ml	706			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS G (CSF)

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	IGLCR-B00/AU	1 x 55 ml	1 x 10 ml	500	36 months	CSF	3,375 - 249,78 (mg/l)
	IGLCR-H00/AU	2 x 55 ml	2 x 10 ml	1000			
KIT for AU®400-640-2700-5400	IGLCR-B00/AU	1 x 55 ml	1 x 10 ml	400	36 months	CSF	3,375 - 249,78 (mg/l)
	IGLCR-H00/AU	2 x 55 ml	2 x 10 ml	800			
KIT for AU®5800	IGLCR-B00/AU	1 x 55 ml	1 x 10 ml	667	36 months	CSF	3,375 - 249,78 (mg/l)
	IGLCR-H00/AU	2 x 55 ml	2 x 10 ml	1333			
KIT for DxC 700 AU®	IGLCR-B00/AU	1 x 55 ml	1 x 10 ml	500	36 months	CSF	3,375 - 249,78 (mg/l)
	IGLCR-H00/AU	2 x 55 ml	2 x 10 ml	1000			

CALIBRATORS

KIT 5 LEV	IGREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	IGCOX-002	1x2 ml	24 months
M	IGCON-002	1x2 ml	24 months
L	IGCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS M

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	IMTUR-B00/AU	1x50 ml	1x10 ml	370	36 months	Serum	0.038 - 7 (g/l)
	IMTUR-H00/AU	2x50 ml	2x10 ml	741			
KIT for AU®400-640-2700-5400	IMTUR-B00/AU	1x50 ml	1x10 ml	299	36 months	Serum	0.038 - 7 (g/l)
	IMTUR-H00/AU	2x50 ml	2x10 ml	599			
KIT for AU®5800	IMTUR-B00/AU	1x50 ml	1x10 ml	500	36 months	Serum	0.038 - 7 (g/l)
	IMTUR-H00/AU	2x50 ml	2x10 ml	1000			
KIT for DxC 700 AU®	IMTUR-B00/AU	1x50 ml	1x10 ml	370	36 months	Serum	0.038 - 7 (g/l)
	IMTUR-H00/AU	2x50 ml	2x10 ml	741			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

LIPOPROTEIN (A)							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	LPTUR2-H00/AU	1x50 ml	1x6 ml	400	12 months	Serum & plasma	0,0235 - 0,946 (g/l)
KIT for AU®400-640-2700-5400	LPTUR2-H00/AU	1x50 ml	1x6 ml	240	12 months	Serum & plasma	0,0235 - 0,946 (g/l)
KIT for AU®5800	LPTUR2-H00/AU	1x50 ml	1x6 ml	462	12 months	Serum & plasma	0,0235 - 0,946 (g/l)
KIT for DxC 700 AU®	LPTUR2-H00/AU	1x50 ml	1x6 ml	400	12 months	Serum & plasma	0,0235 - 0,946 (g/l)

CALIBRATORS

KIT 4 LEV	LPREK-000	4x0,5 ml	18 months
-----------	-----------	----------	-----------

CONTROLS

M	LPCON-002	1x2 ml	24 months
---	-----------	--------	-----------

MICROALBUMIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	MATUR-B00/AU	1x50 ml	1x10 ml	417	36 months	Urine	1,86 - 400 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	833			
	MATUR-L00/AU	12x50 ml	12x10 ml	5000			
KIT for AU®400-640-2700-5400	MATUR-B00/AU	1x50 ml	1x10 ml	333	36 months	Urine	1,86 - 400 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	667			
	MATUR-L00/AU	12x50 ml	12x10 ml	4000			
KIT for AU®5800	MATUR-B00/AU	1x50 ml	1x10 ml	556	36 months	Urine	1,86 - 400 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	1111			
	MATUR-L00/AU	12x50 ml	12x10 ml	6667			
KIT for DxC 700 AU®	MATUR-B00/AU	1x50 ml	1x10 ml	417	36 months	Urine	1,86 - 400 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	833			
	MATUR-L00/AU	12x50 ml	12x10 ml	5000			

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MACOX-002	1x2 ml	24 months
M	MACON-002	1x2 ml	24 months
L	MACOS-002	1x2 ml	24 months

PREALBUMIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	PATUR-B00/AU	1x50 ml	1x10 ml	370	36 months	Serum & plasma	6.53 - 1000 (mg/l)
	PATUR-H00/AU	2x50 ml	2x10 ml	741			
KIT for AU®400-640-2700-5400	PATUR-B00/AU	1x50 ml	1x10 ml	299	36 months	Serum & plasma	6.53 - 1000 (mg/l)
	PATUR-H00/AU	2x50 ml	2x10 ml	599			
KIT for AU®5800	PATUR-B00/AU	1x50 ml	1x10 ml	500	36 months	Serum & plasma	6.53 - 1000 (mg/l)
	PATUR-H00/AU	2x50 ml	2x10 ml	1000			
KIT for DxC 700 AU®	PATUR-B00/AU	1x50 ml	1x10 ml	370	36 months	Serum & plasma	6.53 - 1000 (mg/l)
	PATUR-H00/AU	2x50 ml	2x10 ml	741			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

RETINOL BINDING PROTEIN (RBP)

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	RBTUR-B00/AU	1 x 30,5 ml	1 x 10 ml	263	36 months	Serum & plasma	7,81 - 141 (mg/l)
	RBTUR-H00/AU	2 x 30,5ml	2 x 10 ml	526			
	RBTUR-L00/AU	12 x 30,5 ml	12 x 10 ml	3158			
KIT for AU®400-640-2700-5400	RBTUR-B00/AU	1 x 30,5 ml	1 x 10 ml	263	36 months	Serum & plasma	7,81 - 141 (mg/l)
	RBTUR-H00/AU	2 x 30,5ml	2 x 10 ml	526			
	RBTUR-L00/AU	12 x 30,5 ml	12 x 10 ml	3158			
KIT for AU®5800	RBTUR-B00/AU	1 x 30,5 ml	1 x 10 ml	400	36 months	Serum & plasma	7,81 - 141 (mg/l)
	RBTUR-H00/AU	2 x 30,5ml	2 x 10 ml	800			
	RBTUR-L00/AU	12 x 30,5 ml	12 x 10 ml	4800			
KIT for DxC 700 AU®	RBTUR-B00/AU	1 x 30,5 ml	1 x 10 ml	333	36 months	Serum & plasma	7,81 - 141 (mg/l)
	RBTUR-H00/AU	2 x 30,5ml	2 x 10 ml	667			
	RBTUR-L00/AU	12 x 30,5 ml	12 x 10 ml	4000			

CALIBRATORS

KIT 5 LEV	RBREK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

H	RBCOX-002	1x2 ml	18 months
M	RBCON-002	1x2 ml	18 months
L	RBCOS-002	1x2 ml	18 months

RHEUMATOID FACTOR (RF) - Colloidal Gold

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	RFCOL-B00/AU	1x50 ml	1x20 ml	230	18 months	Serum & plasma	2.228 - 120 (kUI/l)
	RFCOL-H00/AU	2x50 ml	2x20 ml	460			
KIT for AU®400-640-2700-5400	RFCOL-B00/AU	1x50 ml	1x20 ml	200	18 months	Serum & plasma	2.228 - 120 (kUI/l)
	RFCOL-H00/AU	2x50 ml	2x20 ml	400			
KIT for AU®5800	RFCOL-B00/AU	1x50 ml	1x20 ml	333	18 months	Serum & plasma	2.228 - 120 (kUI/l)
	RFCOL-H00/AU	2x50 ml	2x20 ml	667			
KIT for DxC 700 AU®	RFCOL-B00/AU	1x50 ml	1x20 ml	333	18 months	Serum & plasma	2.228 - 120 (kUI/l)
	RFCOL-H00/AU	2x50 ml	2x20 ml	667			

CALIBRATORS

KIT 5 LEV	RFREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	RFCOX-002	1x2 ml	24 months
M	RFCON-002	1x2 ml	24 months
L	RFCOS-002	1x2 ml	24 months

TRANSFERRIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	TFTUR-H00/AU	1x60 ml	1x16,5 ml	236	36 months	Serum & plasma	0.076 - 8 (g/l)
	TFTUR-L00/AU	7x60 ml	7x16,5 ml	1650			
KIT for AU®400-640-2700-5400	TFTUR-H00/AU	1x60 ml	1x16,5 ml	295	36 months	Serum & plasma	0.076 - 8 (g/l)
	TFTUR-L00/AU	7x60 ml	7x16,5 ml	2063			
KIT for AU®5800	TFTUR-H00/AU	1x60 ml	1x16,5 ml	393	36 months	Serum & plasma	0.076 - 8 (g/l)
	TFTUR-L00/AU	7x60 ml	7x16,5 ml	2750			
KIT for AU®DXC	TFTUR-H00/AU	1x60 ml	1x16,5 ml	295	36 months	Serum & plasma	0.076 - 8 (g/l)
	TFTUR-L00/AU	7x60 ml	7x16,5 ml	2063			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

Please contact us for more information on the following tests:

CALPROTECTIN FECAL

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
----------	-----------	---------	--------	--------------	-----------	----------------------



-
- (*) Tests :** Amount of tests possible by kit. May vary from an instrument to another one.
- (**) Shelf Life :** Diagam provides to distributors a minimum of *12 months*. The time period written in this document shows the maximum shelf life time after manufacturing.
- (***) Sample :** Diagam may provide some applications for other samples, such as CSF, Urine, Veterinary samples etc.... Please contact your distributor or us directly at mail@diagam.com or dial +32 68 55 14 82.
- (****) Analytical range :** For general information purposes only. More details can be found in the package insert (instructions for use). Please contact your distributor or us directly at mail@diagam.com or dial +32 68 55 14 82 in order to get the updated data.



Medical device for in vitro diagnostics.

Contact information:

For orders, invoicing and technical queries, we kindly recommend you to use our general e-mail address that is mail@diagam.com.



DIAGAM SA, *Rue du Parc Industriel, 7822 Ghislenghien, BELGIUM*