

ALINITY LINE

SARS-CoV-2 SEROLOGICAL ASSAY **New!**

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
----------	-----------	---------	--------	-------------	----------	------------------

Please contact us for more information

SARS-CoV-2 SEROLOGICAL ASSAY based on PECIA Colloidal Gold



ALBUMIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	ALTUR-B00/ALI	1 x 20 ml R1 + 1 x 7,5 ml R2		100	36 months	serum & plasma	1,096 to 120,44 (g/L)
	ALTUR-H00/ALI	1 x 74 ml R1 + 1 x 19,5 ml R2		343			
	ALTUR-L00/ALI	10 x 74 ml R1 + 10 x 19,5 ml R2		3435			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	2-8 °C	24 months
------------------	-----------	--------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALBUMIN CSF							
REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	MATUR-B00/ALI	1 x 22,5 ml R1 + 1 x 7 ml R2		100	36 months	cerebrospinal fluid	10,86 - 494,33 mg/L
	MATUR-H00/ALI	1 x 70 ml R1 + 1 x 17,6 ml R2		342			
	MATUR-L00/ALI	5 x 70 ml R1 + 5 x 17,6 ml R2		1712			

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MACOZ-002	1x2 ml	24 months
H	MACOX-002	1x2 ml	24 months
M	MACON-002	1x2 ml	24 months
L	MACOS-002	1x2 ml	24 months

ALPHA-1 ACID-GLYCOPROTEIN (OROSOMUCOID)							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	GPTUR-B00/ALI	1 x 30,5 ml R1 + 1 x 8,6 ml R2		100	36 months	serum	0,143 - 2,42 g/L
	GPTUR-H00/ALI	1 x 70 ml R1 + 1 x 17,2 ml R2		243			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA-1 ANTITRYPSIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	ATTUR-B00/ALI	1 x 31 ml R1 + 1 x 8,6 ml R2	100	36 months	serum & plasma	0,007 to 3,5 (g/L)
	ATTUR-H00/ALI	1 x 70 ml R1 + 1 x 17,1 ml R2	242			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA-1 MICROGLOBULIN (URINE)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	A1DGX-B00/ALI	1 x 14,5 ml R1 + 1 x 14 ml R1	100	17 months	urine	1,29 to 166,6 (mg/L)
	A1DGX-H00/ALI	1 x 46 ml R1 + 1 x 45,5 ml R2	376			

CALIBRATORS

KIT 5 LEV	A1REK-000	5x1 ml	14 months
-----------	-----------	--------	-----------

CONTROLS

M	A1CON-002	1x2 ml	14 months
L	A1COS-002	1x2 ml	14 months

ALPHA-2 MACROGLOBULIN (SERUM/PLASMA and URINE)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
SERUM/PLASMA KIT	MGTUR-B00/ALI	1 x 30,5 ml R1 + 1 x 8,6 ml R2	100	36 months	serum & plasma	0,026 to 4,96 (g/L) serum
	MGTUR-H00/ALI	1 x 74 ml R1 + 1 x 18 ml R2	257			
	MGTUR-L00/ALI	5 x 74 ml R1 + 5 x 18 ml R2	1283			
URINARY KIT	MGTUR-B00/ALI	1 x 30,5 ml R1 + 1 x 8,6 ml R2	137	36 months	urine	3,36 - 199 (mg/L) urine
	MGTUR-H00/ALI	1 x 74 ml R1 + 1 x 18 ml R2	352			
	MGTUR-L00/ALI	5 x 74 ml R1 + 5 x 18 ml R2	1758			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

Serum H	MPCOX-002	1x2 ml	24 months
Serum M	MPCON-002	1x2 ml	24 months
Serum L	MPCOS-002	1x2 ml	24 months
Urinary LOW	MVCU1-002	1x2 ml	24 months
Urinary HIGH	MVCU2-002	1x2 ml	24 months

ANTI-STREPTOLYSIN O

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	ASLTX-B00/ALI	R1 1 x 29 ml R1 + 1 x 6,1 ml R1	100	20 months	serum & plasma	19,99 to 864,85 (IU/ml)
	ASLTX-H00/ALI	R2 1 x 70 ml R1 + 1 x 11,7 ml R2	259			

CALIBRATORS

KIT 5 LEV	ASREK-000	5x1 ml	18 months
------------------	-----------	--------	-----------

CONTROLS

M	ASCON-002	1x2 ml	12 months
----------	-----------	--------	-----------

APOLIPOPROTEIN A1

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	AATUR-B00/ALI	R1 1 x 31 ml R1 + 1 x 10,4 ml R1	101	36 months	serum & plasma	0,029 - 4 g/L
	AATUR-H00/ALI	R2 1 x 70 ml R1 + 1 x 21,4 ml R2	243			

CALIBRATORS

KIT 5 LEV	ABREK-000	5 x 0,5 ml	18 months
------------------	-----------	------------	-----------

CONTROLS

H	ABCOX-002	1x2 ml	18 months
M	ABCON-002	1x2 ml	18 months

APOLIPOPROTEIN B

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	ABTUR-B00/ALI	R1 1 x 21,6 ml R1 + 1 x 6,7 ml R1	98	36 months	serum & plasma	0,019 - 4 g/L
	ABTUR-H00/ALI	R2 1 x 70 ml R1 + 1 x 17,6 ml R2	360			

CALIBRATORS

KIT 5 LEV	ABREK-000	5 x 0,5 ml	18 months
------------------	-----------	------------	-----------

CONTROLS

H	ABCOX-002	1x2 ml	18 months
M	ABCON-002	1x2 ml	18 months

BETA2-MICROGLOBULIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1 R2				
KIT	B2DGX-B00/ALI	1 x 23 ml R1 + 1 x 11,8 ml R2	100	30 months	serum & plasma	0,132 - 95 mg/L
	B2DGX-H00/ALI	1 x 70 ml R1 + 1 x 33,5 ml R2	334			
	B2DGX-L00/ALI	5 x 70 ml R1 + 5 x 33,5 ml R2	1672			

CALIBRATORS

KIT 5 LEV	B2REK-000	5x1 ml	30 months
-----------	-----------	--------	-----------

CONTROLS

H	B2COX-001	1x1 ml	30 months
M	B2CON-003	1x3 ml	30 months
L	B2COS-003	1x3 ml	30 months

BETA2-MICROGLOBULIN URINARY

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1 R2				
KIT	BUDGX-B00/ALI	1 x 18,2 ml R1 + 1 x 9,7 ml R2	100	30 months	urine	0,055 - 9,3 mg/L
	BUDGX-H00/ALI	1 x 74 ml R1 + 1 x 35,7 ml R2	468			

CALIBRATORS

KIT 5 LEV	BUREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

M	B2CON-003	1x3 ml	30 months
L	B2COS-003	1x3 ml	30 months

CALPROTECTIN FECAL - Colloidal Gold

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1 R2				
KIT	CACOL-B00/ALI	1 x 25 ml R1 + 1 x 12 ml R2	205	18 months	faeces	10,93 to 1297,4 (mg/kg)
	CACOL-H00/ALI	2 x 25 ml R1 + 2 x 12 ml R2	409			
	CACOL-L00/ALI	3 x 25 ml R1 + 3 x 12 ml R2	614			

CALIBRATORS

KIT 6 LEV	CAREK-000	6x1 ml	12 months
-----------	-----------	--------	-----------

CONTROLS

H	CACOX-002	1x2 ml	12 months
M	CACON-002	1x2 ml	12 months
L	CACOS-002	1x2 ml	12 months

SAMPLE DILUTION BUFFER

KIT	SDBUF-B00	2x70 ml	15 months
	SDBUF-H00	8x70 ml	15 months
	SDBUF-L00	16x70 ml	15 months

CALPRO EASY EXTRACT™

KIT	CAL0510	50 devices
-----	---------	------------

CERULOPLASMIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	CETUR-B00/ALI	1 x 25,1 ml R1 + 1 x 7,5 ml R2		100	36 months	serum & plasma	9,48 - 1051,2 mg/L
	CETUR-H00/ALI	1 x 70 ml R1 + 1x 17,4 ml R2		301			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

COMPLEMENT C3							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	C3TUR-B00/ALI	1 x 24 ml R1 + 1 x 7,2 ml R2		100	36 months	serum & plasma	0,015 - 4,31 g/L
	C3TUR-H00/ALI	1 x 70 ml R1 + 1x 17,4 ml R2		322			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

COMPLEMENT C4							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	C4TUR-B00/ALI	1 x 25,1 ml R1 + 1 x 7,5 ml R2		100	36 months	serum & plasma	0,002 - 0,825 g/L
	C4TUR-H00/ALI	1 x 70 ml R1 + 1x 17,4 ml R2		301			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CRP GOLD LATEX

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i> <i>R2</i>				
KIT	CPDGX-B00/ALI	1 x 14,5 ml R1 + 1 x 14,1 ml R2	100	24 months	serum & plasma	0,031 - 300 mg/L
	CPDGX-H00/ALI	1 x 46 ml R1 + 1 x 45,5 ml R2	376			
	CPDGX-L00/ALI	5 x 46 ml R1 + 5 x 45,5 ml R2	1882			

CALIBRATORS

KIT 5 LEV	CPRLK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	24 months
M	CPCON-002	1x2 ml	24 months
L	CXCON-002	1x2 ml	24 months

CRP WIDE RANGE

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i> <i>R2</i>				
KIT	CPTXL-B00/ALI	1 x 31 ml R1 + 1 x 7,1 ml R2	100	36 months	serum & plasma	0,607 - 425 mg/L
	CPTXL-H00/ALI	1 x 74 ml R1 + 1 x 14,1 ml R2	256			
	CPTXL-L00/ALI	5 x 74 ml R1 + 5 x 14,1 ml R2	1281			

CALIBRATORS

KIT 5 LEV	CPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	24 months
M	CPCON-002	1x2 ml	24 months
L	CXCON-002	1x2 ml	24 months

CYSTATIN C - Colloidal Gold

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i> <i>R2</i>				
KIT	CYCOL-B00/ALI	1 x 23 ml R1 + 1 x 11,3 ml R2	100	24 months	serum & plasma	0,15 - 9,355 mg/L
	CYCOL-H00/ALI	1 x 70 ml R1 + 1 x 31,7 ml R2	334			
	CYCOL-L00/ALI	5 x 70 ml R1 + 5 x 31,7 ml R2	1672			

CALIBRATORS

KIT 5 LEV	CYREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	CYCOX-002	1x2 ml	24 months
M	CYCON-002	1x2 ml	24 months
L	CYCOX-002	1x2 ml	24 months

FERRITIN - Colloïdal Gold

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	FNCOL-B00/ALI	1 x 25,1 ml R1 + 1 x 7,5 ml R2	100	12 months	serum	6,344 - 500 µg/L
	FNCOL-H00/ALI	1 x 70 ml R1 + 1 x 17,4 ml R2	301			

CALIBRATORS

KIT 5 LEV	FNRGK-000	5x1 ml	18 months
------------------	-----------	--------	-----------

CONTROLS

M	FNCON-002	1x2 ml	18 months
----------	-----------	--------	-----------

HAPTOGLOBIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	HATUR-B00/ALI	1 x 30,5 ml R1 + 1 x 8,6 ml R2	100	36 months	serum & plasma	0,002 - 3,86 g/L
	HATUR-I00/ALI	5 x 70 ml R1 + 5 x 17,2 ml R2	1214			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS A

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	IATUR-B00/ALI	1 x 21,5 ml R1 + 1 x 7 ml R2	100	36 months	serum	0,02 - 8 g/L
	IATUR-H00/ALI	1 x 70 ml R1 + 1 x 18,5 ml R2	363			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS G

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
		R1 R2				
KIT	IGTUR-B00/ALI	1 x 21,6 ml R1 + 1 x 10,1 ml R2	100	36 months	serum	0,125 to 35 (g/L)
	IGTUR-H00/ALI	1 x 70 ml R1 + 1 x 29,5 ml R2	358			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS G (CSF)

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
		R1 R2				
KIT	IGLCR-B00/ALI	1 x 18,5 ml R1 + 1 x 5,9 ml R2	100	36 months	cerebrospinal fluid	3,375 to 249,78 (mg/L)
	IGLCR-H00/ALI	1 x 74 ml R1 + 1 x 18 ml R2	467			

CALIBRATORS

KIT 5 LEV	IGREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	IGCOX-002	1x2 ml	24 months
M	IGCON-002	1x2 ml	24 months
L	IGCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS M

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
		R1 R2				
KIT	IMTUR-B00/ALI	1 x 21,6 ml R1 + 1 x 7 ml R2	100	36 months	serum	0,038 - 7 g/L
	IMTUR-H00/ALI	1 x 70 ml R1 + 1 x 18,5 ml R2	363			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

LIPOPROTEIN A							
REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	LPTUR2-B00/ALI	1 x 25,5 ml R1	+ 1 x 5,9 ml R2	100	12 months	serum & plasma	0,0235 - 0,946 g/L
	LPTUR2-H00/ALI	1 x 70 ml R1	+ 1 x 12,5 ml R2	300			

CALIBRATORS

KIT 4 LEV	LPREK-000	4 x 0,5 ml	18 months
-----------	-----------	------------	-----------

CONTROLS

M	LPCON-002	1 x 2 ml	24 months
---	-----------	----------	-----------

MICROALBUMIN							
REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	MATUR-B00/ALI	1 x 22,5 ml R1	+ 1 x 7 ml R2	100	36 months	urine	1,860 - 400 mg/L
	MATUR-H00/ALI	1 x 70 ml R1	+ 1 x 17,6 ml R2	342			
	MATUR-L00/ALI	5 x 70 ml R1	+ 5 x 17,6 ml R2	1712			

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MACOX-002	1x2 ml	24 months
M	MACON-002	1x2 ml	24 months
L	MACOS-002	1x2 ml	24 months

PREALBUMIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	PATUR-B00/ALI	1 x 30,5 ml R1	+ 1 x 8,6 ml R2	100	36 months	serum & plasma	6,53 to 1000 (mg/L)
	PATUR-H00/ALI	1 x 74 ml R1	+ 1 x 18 ml R2	257			
	PATUR-L00/ALI	5 x 74 ml R1	+ 5 x 18 ml R2	1283			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

RETINOL BINDING PROTEIN (RBP)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	RBTUR-B00/ALI	1 x 20 ml R1 + 1 x 8,6 ml R2	100	36 months	serum & plasma	7,81 to 141 (mg/L)
	RBTUR-H00/ALI	1 x 70 ml R1 + 1 x 26,5 ml R2	398			
	RBTUR-L00/ALI	5 x 70 ml R1 + 5 x 26,5 ml R2	1992			

CALIBRATORS

KIT 5 LEV	RBREK-000	5x1 ml	18 months
------------------	-----------	--------	-----------

CONTROLS

H	RBCOX-002	1x2 ml	18 months
M	RBCON-002	1x2 ml	18 months
L	RBCOS-002	1x2 ml	18 months

RHEUMATOID FACTOR - Colloidal Gold

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	RFCOL-B00/ALI	1 x 25,5 ml R1 + 1 x 17,5 ml R2	100	18 months	serum & plasma	2,228 - 120 kIU/L
	RFCOL-H00/ALI	1 x 68 ml R1 + 1 x 46 ml R2	290			

CALIBRATORS

KIT 5 LEV	RFREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	RFCOX-002	1x2 ml	24 months
M	RFCON-002	1x2 ml	24 months
L	RFCOS-002	1x2 ml	24 months

TRANSFERRIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	TFTUR-B00/ALI	1 x 31 ml R1 + 1 x 10,8 ml R2	100	36 months	serum & plasma	0,076 - 8 g/L
	TFTUR-H00/ALI	1 x 70 ml R1 + 1 x 22,4 ml R2	243			
	TFTUR-L00/ALI	5 x 70 ml R1 + 5 x 22,4 ml R2	1213			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

Please contact us for more information on the following tests:

HDL CHOLESTEROL						
<i>REAGENTS</i>	<i>Reference</i>	<i>Content</i>	<i>Tests*</i>	<i>Shelf Life*</i>	<i>Sample**</i>	<i>Analytical Range</i>

IMMUNOGLOBULIN E						
<i>REAGENTS</i>	<i>Reference</i>	<i>Content</i>	<i>Tests*</i>	<i>Shelf Life*</i>	<i>Sample**</i>	<i>Analytical Range</i>

Your expert in specific proteins



-
- (*) Tests :** Amount of real tests. May vary from a analyzer to another one.
- (**) Shelf Life :** Diagam provides to distributors a minimum of **12 months** . The time period written in this document shows the maximum shelf life time after manufacturing.
- (***) Sample :** Diagam may provide some applications for other samples, such as CSF, Urine, Veterany samples etc.... Please contact your distributor or us directly at mail@diagam.com or dial **+32 68 55 14 82**.
- (****) Analytical range :** For general information purposes only. More details can be found in the package insert (instructions for use). Please contact your distributor or us directly at mail@diagam.com or dial **+32 68 55 14 82** in order to get the updated data.



Medical device for in vitro diagnosis.



DIAGAM SA, *Rue du Parc Industriel 40, 7822 Ghislenghien, BELGIUM*



The products fulfill the requirements of European directive 98/79 EC concerning the medical devices of in vitro diagnosis.

© Diagam. 2023. All rights reserved. All product names, registered trademarks, company names in this document remain the property of their respective owners. For general information purposes only. More details can be found in the package insert (instructions for use) and labelling. Not for the US market except for research use only.