

ARCHITECT LINE

New

SARS-CoV-2 SEROLOGICAL ASSAY

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
----------	-----------	---------	--------	--------------	-----------	----------------------

Please contact us for more information

SARS-CoV-2 SEROLOGICAL ASSAY based on PECIA Colloidal Gold



CATALOGUE ARCHITECT LINE

ALBUMIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	ALTUR-B00/ARC	1 x 40 ml R1 + 1 x 10 ml R2	169	36 months	serum or plasma	0.40 to 121.71 (g/l)
	ALTUR-H00/ARC	1 x 80 ml R1 + 1 x 20 ml R2	372			
	ALTUR-L00/ARC	7 x 80 ml R1 + 7 x 20 ml R2	2604			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALBUMIN CSF

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	MATUR-B00/ARC	1 x 22,5 ml R1 + 1 x 6,1 ml R2	100	36 months	cerebrospinal fluid	10.86 to 494.33 (mg/l)
	MATUR-H00/ARC	1 x 74 ml R1 + 1 x 17,5 ml R2	361			
	MATUR-L00/ARC	7 x 75 ml R1 + 7 x 15 ml R2	2126			

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MACOZ-002	1 x 2 ml	24 months
H	MACOX-002	1x2 ml	24 months
M	MACON-002	1x2 ml	24 months
L	MACOS-002	1x2 ml	24 months

ALPHA-1 GLYCOPROTEIN ACIDE/OROSOMUCOID

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	GPTUR-B00/ARC	1 x 37,5 ml R1 + 1 x 7,5 ml R2	97	36 months	serum	0.143 to 2.42 (g/l)
	GPTUR-H00/ARC	2 x 75 ml R1 + 2 x 15 ml R2	443			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA-1 ANTITRYPSIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	ATTUR-B00/ARC	R1 1 x 45 ml R1 + 1 x 9 ml R2	122	24 months	serum or plasma	0.007 to 3.5 (g/l)
	ATTUR-H00/ARC	R2 2 x 75 ml R1 + 2 x 15 ml R2	443			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA-1 MICROGLOBULIN (URINE)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	A1DGX-B00/ARC	R1 1 x 15 ml R1 + 1 x 15 ml R2	150	17 months	urine	1,29 to 166,6 (mg/l)
	A1DGX-H00/ARC	R2 1 x 55 ml R1 + 1 x 55 ml R2	456			

CALIBRATORS

KIT 5 LEV	A1REK-000	5x1 ml	14 months
------------------	-----------	--------	-----------

CONTROLS

M	A1CON-002	1x2 ml	14 months
L	A1COS-002	1x2 ml	14 months

ALPHA-2 MACROGLOBULIN - ALPHA-2 MACROGLOBULIN URINARY

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	MGTUR-H00/ARC	R1 1 x 55 ml R1 + 1 x 11 ml R2	155	36 months	serum/plasma or urine	0.026 to 4.96 (g/l) serum 3.36 to 199 (mg/l) urine
	MGTUR-L00/ARC	R2 10 x 55 ml R1 + 10 x 11 ml R2	1550			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months
URINARY LOW	MGCU1-002	1x2 ml	24 months
URINARY HIGH	MGCU2-002	1x2 ml	24 months

ANTI STREPTOLYSIN O (ASO)							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	ASLTX-B00/ARC	1 x 70 ml R1 + 1 x 8 ml R2		296	20 months	serum or plasma	19.99 to 864.85 (IU/ml)
	ASLTX-H00/ARC	2 x 87 ml R1 + 2 x 10 ml R2		741			
CALIBRATORS							
H	ASREK-000	5x1 ml			18 months		
CONTROLS							
M	ASCON-002	1x2 ml			12 months		

APOLIPOPROTEIN A1							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	AATUR-B00/ARC	1 x 38 ml R1 + 1 x 10 ml R2		107	36 months	serum or plasma	0.029 to 4.0 (g/l)
	AATUR-H00/ARC	2 x 75 ml R1 + 2 x 20 ml R2		474			
CALIBRATORS							
KIT 5 LEV	ABREK-000	5x0.5 ml			18 months		
CONTROLS							
H	ABCOX-002	1x2 ml			18 months		
M	ABCON-002	1x2 ml			18 months		

APOLIPOPROTEIN B							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	ABTUR-B00/ARC	1 x 41 ml R1 + 1 x 8 ml R2		151	36 months	serum or plasma	0.019 to 4.0 (g/l)
	ABTUR-H00/ARC	2 x 67 ml R1 + 2 x 13 ml R2		543			
CALIBRATORS							
KIT 5 LEV	ABREK-000	5x0.5 ml			18 months		
CONTROLS							
H	ABCOX-002	1x2 ml			18 months		
M	ABCON-002	1x2 ml			18 months		

BETA 2 MICROGLOBULIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	B2DGX-B00/ARC	1 x 27 ml R1 + 1 x 12 ml R2		111	30 months	serum or plasma	0,132 to 95 (mg/l)
	B2DGX-H00/ARC	2 x 45 ml R1 + 2 x 20 ml R2		396			
	B2DGX-L00/ARC	10 x 45 ml R1 + 10 x 20 ml R2		1981			
CALIBRATORS							
KIT 5 LEV	B2REK-000	5x1 ml			30 months		
CONTROLS							
H	B2COX-001	1 x 1 ml			30 months		
M	B2CON-003	1x3 ml			30 months		
L	B2COS-003	1x3 ml			30 months		
BETA 2 MICROGLOBULIN (URINE)							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	BUDGX-B00/ARC	1 x 23 ml R1 + 1 x 10 ml R2		117	30 months	Urine	0.055 to 9.3 (mg/l)
	BUDGX-H00/ARC	2 x 23 ml R1 + 2 x 10 ml R2		234			
CALIBRATORS							
KIT 5 LEV	BUREK-000	5x1 ml			24 months		
CONTROLS							
M	B2CON-003	1x3 ml			30 months		
L	B2COS-003	1x3 ml			30 months		
CERULOPLASMIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	CETUR-B00/ARC	1 x 40 ml R1 + 1 x 8 ml R2		128	36 months	serum or plasma	9.48 to 1051.20 (mg/l)
	CETUR-H00/ARC	2 x 60 ml R1 + 2 x 12 ml R2		419			
CALIBRATORS							
KIT 5 LEV	MPREK-000	5x1 ml			24 months		
CONTROLS							
H	MPCOX-002	1x2 ml			24 months		
M	MPCON-002	1x2 ml			24 months		
L	MPCOS-002	1x2 ml			24 months		

COMPLEMENT C3

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	C3TUR-B00/ARC	1 x 41 ml R1 + 1 x 8 ml R2		137	36 months	serum or plasma	0.015 to 4.31 (g/l)
	C3TUR-H00/ARC	2 x 56 ml R1 + 2 x 11 ml R2		405			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

COMPLEMENT C4

REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	C4TUR-B00/ARC	1 x 40 ml R1 + 1 x 8 ml R2		128	36 months	serum or plasma	0.002 to 0.825 (g/l)
	C4TUR-H00/ARC	2 x 60 ml R1 + 2 x 12 ml R2		419			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CRP GOLD LATEX

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	CPDGX-B00/ARC	1 x 15 ml R1 + 1 x 15 ml R2		105	24 months	serum or plasma	0.031 to 300 (mg/l)
	CPDGX-H00/ARC	1 x 55 ml R1 + 1 x 55 ml R2		456			
	CPDGX-L00/ARC	10 x 55 ml R1 + 10 x 55 ml R2		4561			

CALIBRATORS

KIT 5 LEV	CPRLK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	24 months
M	CPCON-002	1x2 ml	24 months
L	CXCON-002	1x2 ml	24 months

CRP XL WIDE RANGE

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	CPTXL-B00/ARC	1 x 55,5 ml R1 + 1 x 8 ml R2		140	36 months	serum or plasma	0.607 to 425 (mg/l)
	CPTXL-H00/ARC	2 x 85 ml R1 + 2 x 12 ml R2		459			
	CPTXL-L00/ARC	7 x 85 ml R1 + 7 x 12 ml R2		1607			

CALIBRATORS

KIT 5 LEV	CPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	24 months
M	CPCON-002	1x2 ml	24 months
L	CXCON-002	1x2 ml	24 months

CYSTATIN C - Colloidal Gold

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	CYCOL-B00/ARC	1 x 27 ml R1 + 1 x 11 ml R2		107	24 months	serum or plasma	0.15 to 9.355 (mg/l)
	CYCOL-H00/ARC	2 x 48 ml R1 + 2 x 20 ml R2		421			
	CYCOL-L00/ARC	10 x 48 ml R1 + 10 x 20 ml R2		2103			

CALIBRATORS

KIT 5 LEV	CYREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	CYCOX-002	1x2 ml	24 months
M	CYCON-002	1x2 ml	24 months
L	CYCOS-002	1x2 ml	24 months

FERRITIN - Colloidal Gold

REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	FNCOL-B00/ARC	1 x 40 ml R1 + 1 x 8 ml R2		128	12 months	serum	6.344 to 500 (µg/l)
	FNCOL-H00/ARC	2 x 55 ml R1 + 2 x 11 ml R2		378			

CALIBRATORS

KIT 5 LEV	FNRGK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

M	FNCON-002	1x2 ml	18 months
---	-----------	--------	-----------

HAPTOGLOBIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	HATUR-B00/ARC	1 x 37,5 ml R1 + 1 x 7,5 ml R2		97	36 months	serum or plasma	0.002 to 3.86 (g/l)
	HATUR-H00/ARC	2 x 75 ml R1 + 2 x 15 ml R2		443			
CALIBRATORS							
KIT 5 LEV	MPREK-000	5x1 ml			24 months		
CONTROLS							
H	MPCOX-002	1x2 ml			24 months		
M	MPCON-002	1x2 ml			24 months		
L	MPCOS-002	1x2 ml			24 months		

IMMUNOGLOBULINS A							
REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	IATUR-B00/ARC	1 x 38 ml R1 + 1 x 8 ml R2		144	36 months	serum	0.02 to 8.0 (g/l)
	IATUR-H00/ARC	2 x 54 ml R1 + 2 x 11 ml R2		425			
CALIBRATORS							
KIT 5 LEV	MPREK-000	5x1 ml			24 months		
CONTROLS							
H	MPCOX-002	1x2 ml			24 months		
M	MPCON-002	1x2 ml			24 months		
L	MPCOS-002	1x2 ml			24 months		

IMMUNOGLOBULINS G							
REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	IGTUR-B00/ARC	1 x 25 ml R1 + 1 x 9,5 ml R2		104	36 months	serum	0.125 to 35 (g/l)
	IGTUR-H00/ARC	2 x 50 ml R1 + 2 x 19 ml R2		461			
CALIBRATORS							
KIT 5 LEV	MPREK-000	5x1 ml			24 months		
CONTROLS							
H	MPCOX-002	1x2 ml			24 months		
M	MPCON-002	1x2 ml			24 months		
L	MPCOS-002	1x2 ml			24 months		

IMMUNOGLOBULINS G (CSF)

REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	IGLCR-B00/ARC	1 x 18,5 ml R1 + 1 x 5 ml R2		100	36 months	cerebrospinal fluid	3,375 to 249,78 (mg/l)
	IGLCR-H00/ARC	1 x 74 ml R1 + 1 x 17 ml R2		464			

CALIBRATORS

KIT 5 LEV	IGREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	IGCOX-002	1x2 ml	24 months
M	IGCON-002	1x2 ml	24 months
L	IGCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS M

REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	IMTUR-B00/ARC	1 x 38 ml R1 + 1 x 8 ml R2		144	36 months	serum	0.038 to 7.0 (g/l)
	IMTUR-H00/ARC	2 x 49,5 ml R1 + 2 x 10,5 ml R2		402			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

LIPOPROTEIN (A)

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	LPTUR2-B00/ARC	1 x 64 ml R1 + 1 x 8 ml R2		191	12 months	serum or plasma	0.0235 to 0.946 (g/l)
	LPTUR2-H00/ARC	2 x 80 ml R1 + 2 x 10 ml R2		503			

CALIBRATORS

KIT 4 LEV	LPREK-000	4x0,5 ml	18 months
------------------	-----------	----------	-----------

CONTROLS

M	LPCON-002	1x2 ml	24 months
----------	-----------	--------	-----------

MICROALBUMIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	MATUR-B00/ARC	1 x 22,5 ml R1 + 1 x 6,1 ml R2	100	36 months	urine	1.86 to 400 (mg/l)
	MATUR-H00/ARC	1 x 74 ml R1 + 1 x 17,5 ml R2	361			
	MATUR-L00/ARC	7 x 75 ml R1 + 7 x 15 ml R2	2126			

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MACOX-002	1x2 ml	24 months
M	MACON-002	1x2 ml	24 months
L	MACOS-002	1x2 ml	24 months

PREALBUMIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	PATUR-H00/ARC	1 x 55 ml R1 + 1 x 11 ml R2	155	36 months	serum/plasma	6.53 to 1000 (mg/l)

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

RETINOL BINDING PROTEIN (RBP)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	RBTUR-B00/ARC	1 x 31 ml R1 + 1 x 10 ml R2	138	36 months	serum or plasma	7,81 to 141 (mg/l)
	RBTUR-H00/ARC	2 x 31 ml R1 + 2 x 10 ml R2	277			
	RBTUR-L00/ARC	10 x 31 ml R1 + 10 x 10 ml R2	1383			

CALIBRATORS

KIT 5 LEV	RBREK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

H	RBCOX-002	1x2 ml	18 months
M	RBCON-002	1x2 ml	18 months
L	RBCOS-002	1x2 ml	18 months

RHEUMATOID FACTOR (RF) - Colloidal Gold

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	RFCOL-B00/ARC	1 x 30 ml R1 + 1 x 20 ml R2		122	18 months	serum and plasma	2.228 to 120 (kIU/l)
	RFCOL-H00/ARC	2 x 60 ml R1 + 2 x 40 ml R2		512			

CALIBRATORS

KIT 5 LEV	RFREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	RFCOX-002	1x2 ml	24 months
M	RFCON-002	1x2 ml	24 months
L	RFCOS-002	1x2 ml	24 months

TRANSFERRIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	TFTUR-B00/ARC	1 x 37,5 ml R1 + 1 x 10,5 ml R2		108	36 months	serum and plasma	0.076 to 8.0 (g/l)
	TFTUR-H00/ARC	2 x 72 ml R1 + 2 x 20 ml R2		449			
	TFTUR-L00/ARC	7 x 72 ml R1 + 7 x 20 ml R2		1570			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

Please contact us for more information on the following tests:

CALPROTECTIN FECAL

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
----------	-----------	---------	--------	--------------	-----------	----------------------

IMMUNOGLOBULIN E

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
----------	-----------	---------	--------	--------------	-----------	----------------------

Your expert in specific proteins



(*) Tests : Amount of real tests. May vary from a analyzer to another one.

() Shelf Life :** Diagam provides to distributors a minimum of **12 months**. The time period written in this document shows the maximum shelf life time after manufacturing.

(*) Sample :** Diagam may provide some applications for other samples, such as CSF, Urin, Veterinary samples etc.... Please contact your distributor or us directly at mail@diagam.com or dial **+32 68 55 14 82**.

(**)**
Analytical range : For general information purposes only. More details can be found in the package insert (instructions for use). Please contact your distributor or us directly at mail@diagam.com or dial **+32 68 55 14 82** in order to get the updated data.



Medical device for in vitro diagnosis.



DIAGAM SA, *Rue du Parc Industriel 40, 7822 Ghislenghien, BELGIUM*



The products fulfill the requirements of European directive 98/79 EC concerning the medical devices of in vitro diagnosis.