

Bulk OEM**New****SARS-CoV-2 SEROLOGICAL ASSAY**

REAGENTS	<i>Reference</i>	<i>Content</i>	<i>T°C</i>	<i>Shelf Life*</i>	<i>Sample**</i>	<i>Analytical Range***</i>	<i>BUF/Ab Ratio</i>
-----------------	------------------	----------------	------------	--------------------	-----------------	----------------------------	---------------------

Please contact us for more information

SARS-CoV-2 SEROLOGICAL ASSAY based on PECIA Colloidal Gold

ALBUMIN							
REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	serum & plasma	1.096 to 120.44 (g/l)	150/40 (µl)
Ab	ALTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

ALBUMIN (CSF)							
REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio

BUF	I7BUF-999	Bulk	2-25 °C	36 months	CSF	10.86 to 494.33 (mg/l)	250/50 (µl)
Ab	MATUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MAREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MACOZ-999	Bulk	2-8 °C	24 months
H	MACOX-999	Bulk	2-8 °C	24 months
M	MACON-999	Bulk	2-8 °C	24 months
L	MACOS-999	Bulk	2-8 °C	24 months

ALPHA-1-ANTITRYPSIN							
REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio

BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum & plasma	0.007 to 3.5 (g/l)	250/50 (µl)
Ab	ATTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

ALPHA1-GLYCOPROTEIN ACID / OROSOMUCOID

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	GYBUF-999	Bulk	2-25 °C	36 months	serum	0.143 to 2.42 (g/l)	250/50 (µl)
Ab	GPTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

ALPHA-1 MICROGLOBULIN (URINE)

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	A1BUX-999	Bulk	2-8 °C	17 months	urine	1.29 to 166.6 (mg/l)	100/100 (µl)
Ab	A1DGX-H99	Bulk	2-8 °C	17 months			

CALIBRATORS

KIT 5 LEV	A1REK-999	Bulk	2-8 °C	14 months
-----------	-----------	------	--------	-----------

CONTROLS

M	A1CON-999	Bulk	2-8 °C	14 months
L	A1COS-999	Bulk	2-8 °C	14 months

ALPHA2-MACROGLOBULIN - ALPHA2-MACROGLOBULIN URINARY

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum & plasma	0.026 to 4.96 (g/l)	250/50 (µl)
Ab	MGTUR-H99	Bulk	2-8 °C	36 months	urine	3.36 to 199 (mg/l)	

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months
URINARY LOW	MGCU1-999	Bulk	2-8 °C	24 months
URINARY HIGH	MGCU2-999	Bulk	2-8 °C	24 months

ANTI STREPTOLYSIN O (ASO)

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	ASBUF-999	Bulk	2-8 °C	20 months	serum & plasma	19.992 to 864.85 (IU/ml)	235/30 (µl)
Ab	ASLTX-H99	Bulk	2-8 °C	20 months			

CALIBRATORS

KIT 5 LEV	ASREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

M	ASCQN-999	Bulk	2-8 °C	18 months
---	-----------	------	--------	-----------

ANTITHROMBIN III

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TIBUF-999	Bulk	2-25 °C	36 months	plasma	28.15 to 806.3 (mg/l)	250/50 (µl)
Ab	A3TUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

APOLIPOPROTEIN A1

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum & plasma	0.029 to 4 (g/l)	250/66 (µl)
Ab	AATUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	ABREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

H	ABCOX-999	Bulk	2-8 °C	18 months
M	ABCON-999	Bulk	2-8 °C	18 months

APOLIPOPROTEIN B

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum & plasma	0.019 to 4.0 (g/l)	250/50 (µl)
Ab	ABTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	ABREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

H	ABCOX-999	Bulk	2-8 °C	18 months
M	ABCON-999	Bulk	2-8 °C	18 months

BETA2-MICROGLOBULIN

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	B2BUX-999	Bulk	2-8 °C	30 months	serum & plasma	0,132 to 95 (mg/l)	180/80 (µl)
Ab	B2DGX- H99	Bulk	2-8 °C	30 months			

CALIBRATORS

KIT 5 LEV	B2REK-999	Bulk	2-8 °C	30 months
-----------	-----------	------	--------	-----------

CONTROLS

H	B2COX-001	lyophilized	2-8 °C	30 months
M	B2CON-003	lyophilized	2-8 °C	30 months
L	B2COS-003	lyophilized	2-8 °C	30 months

BETA2-MICROGLOBULIN (URINE)

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	BUBUX-999	Bulk	2-8 °C	30 months	urine	0.055 to 9.3 (mg/l)	135/60 (µl)
Ab	BUDGX- H99	Bulk	2-8 °C	30 months			

CALIBRATORS

KIT 5 LEV	BUREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

M	B2CON-003	lyophilized	2-8 °C	30 months
L	B2COS-003	lyophilized	2-8 °C	30 months

CALPROTECTIN FECAL - Colloidal Gold

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	faeces	10.93 to 1297.4 (mg/kg)	180/75 (µl)
Ab	CACOL- H99	Bulk	2-8 °C	18 months			

CALIBRATORS

KIT 6 LEV	CAREK-999	Bulk	2-8 °C	12 months
-----------	-----------	------	--------	-----------

CONTROLS

H	CACOX-999	Bulk	2-8 °C	12 months
M	CACON-999	Bulk	2-8 °C	12 months
L	CACOS-999	Bulk	2-8 °C	12 months

SAMPLE DILUTION BUFFER

KIT	SDBUF-999	Bulk	2-8 °C	15 months
-----	-----------	------	--------	-----------

CALPRO EASY EXTRACT™

KIT	CAL0510	50 devices
-----	---------	------------

CERULOPLASMIN

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	MEBUF-999	Bulk	2-8 °C	36 months	serum & plasma	9.48 to 1051.2 (mg/l)	250/50 (µl)
Ab	CETUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

COMPLEMENT C3

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	serum & plasma	0.015 to 4.31 (g/l)	250/50 (µl)
Ab	C3TUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

COMPLEMENT C4

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	serum & plasma	0.002 to 0.825 (g/l)	250/50 (µl)
Ab	C4TUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

CRP GOLD LATEX

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	CPBUF-999	Bulk	2-8 °C	24 months	serum & plasma	0.031 to 300 (mg/l)	100/100 (µl)
Ab	CPDGX-H99	Bulk	2-8 °C	24 months			

CALIBRATORS

KIT 5 LEV	CPRLK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	CPCOX-999	Bulk	2-8 °C	24 months
M	CPCON-999	Bulk	2-8 °C	24 months
L	CXCOS-999	Bulk	2-8 °C	24 months

CRP XL WIDE RANGE

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	PITRX-999	Bulk	2-25 °C	36 months	serum & plasma	0.607 to 425 (mg/l)	250/36 (µl)
Ab	CPTXL-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	CPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	CPCOX-999	Bulk	2-8 °C	24 months
M	CPCON-999	Bulk	2-8 °C	24 months
L	CXCON-999	Bulk	2-8 °C	24 months

CYSTATIN C - Colloidal Gold

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum & plasma	0.15 to 9.355 (mg/l)	180/75 (µl)
Ab	CYCOL-H99	Bulk	2-8 °C	24 months			

CALIBRATORS

KIT 5 LEV	CYREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	CYCOX-999	Bulk	2-8 °C	24 months
M	CYCON-999	Bulk	2-8 °C	24 months
L	CYCOS-999	Bulk	2-8 °C	24 months

FERRITIN - Colloidal Gold

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	FDBUF-999	Bulk	2-8 °C	18 months	serum	6.344 to 500 (µg/l)	250/50 (µl)
Ab	FNCOL-H99	Bulk	2-8 °C	12 months			

CALIBRATORS

KIT 5 LEV	FNRGK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

M	FNCON-999	Bulk	2-8 °C	18 months
---	-----------	------	--------	-----------

FIBRINOGEN							
REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	plasma	0.05 to 10 (g/l)	250/50 (µl)
Ab	FITUR-H99	Bulk	2-8 °C	36 months			
CALIBRATORS							
M	FIREM-001	lyophilized	2-8 °C	28 months			
CONTROLS							
M	FICON-001	lyophilized	2-8 °C	28 months			
HAPTOGLOBIN							
REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	I7TRX-999	Bulk	2-25 °C	36 months	serum & plasma	0.002 to 3.86 (g/l)	250/50 (µl)
Ab	HATUR-H99	Bulk	2-8 °C	36 months			
CALIBRATORS							
KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months			
CONTROLS							
H	MPCOX-999	Bulk	2-8 °C	24 months			
M	MPCON-999	Bulk	2-8 °C	24 months			
L	MPCOS-999	Bulk	2-8 °C	24 months			
IMMUNOGLOBULINS A							
REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TIBUF-999	Bulk	2-25 °C	36 months	serum	0.02 to 8.0 (g/l)	250/50 (µl)
Ab	IATUR-H99	Bulk	2-8 °C	36 months			
CALIBRATORS							
KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months			
CONTROLS							
H	MPCOX-999	Bulk	2-8 °C	24 months			
M	MPCON-999	Bulk	2-8 °C	24 months			
L	MPCOS-999	Bulk	2-8 °C	24 months			

IMMUNOGLOBULINS E

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	IEBUX-999	Bulk	2-8 °C	20 months	serum & plasma	21.0 to 1436 (UI/ml)	140/70 (µl)
Ab	IEDGX-H99	Bulk	2-8 °C	20 months			

CALIBRATORS

KIT 5 LEV	IEREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

M	IECON-003	lyophilized	2-8 °C	18 months
L	IECOS-003	lyophilized	2-8 °C	18 months

IMMUNOGLOBULINS G

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TIBUF-999	Bulk	2-25 °C	36 months	serum	0.125 to 35 (g/l)	250/95 (µl)
Ab	IGTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

IMMUNOGLOBULINS G (CSF)

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TIBUF-999	Bulk	2-25 °C	36 months	CSF	3.375 to 249.78 (mg/l)	135/25 (µl)
Ab	IGLCR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	IGREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	IGCOX-999	Bulk	2-8 °C	24 months
M	IGCON-999	Bulk	2-8 °C	24 months
L	IGCOS-999	Bulk	2-8 °C	24 months

IMMUNOGLOBULINS M

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TIBUF-999	Bulk	2-25 °C	36 months	serum	0.038 to 7.0 (g/l)	250/50 (µl)
Ab	IMTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

LIPOPROTEIN (A)

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum & plasma	0.0235 to 0.946 (g/l)	200/25 (µl)
Ab	LPTUR2-H99	Bulk	2-8 °C	12 months			

CALIBRATORS

KIT 4 LEV	LPREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

M	LPCON-999	Bulk	2-8 °C	24 months
---	-----------	------	--------	-----------

MICROALBUMIN

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	urine	1.86 to 400 (mg/l)	250/50 (µl)
Ab	MATUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MAREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MACOX-999	Bulk	2-8 °C	24 months
M	MACON-999	Bulk	2-8 °C	24 months
L	MACOS-999	Bulk	2-8 °C	24 months

PREALBUMIN							
REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	serum & plasma	6.53 to 1000 (mg/l)	250/50 (µl)
Ab	PATUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

RETINOL BINDING PROTEIN							
REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio

BUF	PITRX-999	Bulk	2-25 °C	36 months	serum & plasma	7.81 to 141 (mg/l)	75/25 (µl)
Ab	RBTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	RBREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

H	RBCOX-999	Bulk	2-8 °C	18 months
M	RBCON-999	Bulk	2-8 °C	18 months
L	RBCOS-999	Bulk	2-8 °C	18 months

RHEUMATOID FACTOR (RF) - Colloidal Gold							
REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio

BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum & plasma	2.228 to 120 (kUI/l)	150/100 (µl)
Ab	RFCOL-H99	Bulk	2-8 °C	18 months			

CALIBRATORS

KIT 5 LEV	RFREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	RFCOX-999	Bulk	2-8 °C	24 months
M	RFCON-999	Bulk	2-8 °C	24 months
L	RFCOS-999	Bulk	2-8 °C	24 months

TRANSFERRIN							
REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	17BUF-999	Bulk	2-25 °C	36 months	serum & plasma	0.076 to 8.0 (g/l)	250/70 (µl)
Ab	TFTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
------------------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

Please contact us for more information on the following tests:

HDL CHOLESTEROL							
REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio



(*) Shelf Life : Diagam provides to distributors a minimum of **12 months** . The time period written in this document shows the maximum shelf life time after manufacturing.

() Sample :** Diagam may provide some applications for other samples, such as CSF, Urin, Veterinary samples etc.... Please contact your distributor or us directly at mail@diagam.com or dial **+32 68 55 14 82**.

(*) Analytical range :** For general information purposes only. More details can be found in the package insert (instructions for use). Please contact your distributor or us directly at mail@diagam.com or dial **+32 68 55 14 82** in order to get the updated data.



Medical device for in vitro diagnosis.



DIAGAM SA, *Rue du Parc Industriel, 7822 Ghislenghien BELGIUM*



The products fulfill the requirements of European directive 98/79 EC concerning the medical devices of in vitro diagnosis.