



BECKMAN COULTER SYSTEMS LINE



CATALOGUE 2025

ALBUMIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	ALTUR-B00/AU	1x50 ml	1x15 ml	417	36 months	Serum & plasma	1,096 - 120.44 (g/l)
	ALTUR-H00/AU	2x50 ml	2x15 ml	833			
	ALTUR-L00/AU	12x50 ml	12x15 ml	5000			
KIT for AU®400-640-2700-5400	ALTUR-B00/AU	1x50 ml	1x15 ml	222	36 months	Serum & plasma	1,096 - 120.44 (g/l)
	ALTUR-H00/AU	2x50 ml	2x15 ml	444			
	ALTUR-L00/AU	12x50 ml	12x15 ml	2667			
KIT for AU®5800	ALTUR-B00/AU	1x50 ml	1x15 ml	510	36 months	Serum & plasma	1,096 - 120.44 (g/l)
	ALTUR-H00/AU	2x50 ml	2x15 ml	1020			
	ALTUR-L00/AU	12x50 ml	12x15 ml	6122			
KIT for DxC 700 AU®	ALTUR-B00/AU	1x50 ml	1x15 ml	510	36 months	Serum & plasma	1,096 - 120.44 (g/l)
	ALTUR-H00/AU	2x50 ml	2x15 ml	1020			
	ALTUR-L00/AU	12x50 ml	12x15 ml	6122			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALBUMIN CSF

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	MATUR-B00/AU	1x50 ml	1x10 ml	417	36 months	CSF	10,86 - 494,33 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	833			
	MATUR-L00/AU	12x50 ml	12x10 ml	5000			
KIT for AU®400-640-2700-5400	MATUR-B00/AU	1x50 ml	1x10 ml	333	36 months	CSF	10,86 - 494,33 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	667			
	MATUR-L00/AU	12x50 ml	12x10 ml	4000			
KIT for AU®5800	MATUR-B00/AU	1x50 ml	1x10 ml	556	36 months	CSF	10,86 - 494,33 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	1111			
	MATUR-L00/AU	12x50 ml	12x10 ml	6667			
KIT for DxC 700 AU®	MATUR-B00/AU	1x50 ml	1x10 ml	417	36 months	CSF	10,86 - 494,33 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	833			
	MATUR-L00/AU	12x50 ml	12x10 ml	5000			

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MACOZ-002	1x2 ml	24 months
H	MACOX-002	1x2 ml	24 months
M	MACON-002	1x2 ml	24 months
L	MACOS-002	1x2 ml	24 months

CATALOGUE 2025

ALPHA-1 GLYCOPROTEIN ACIDE/OROSOMUCOID

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	GPTUR-B00/AU	1x60 ml	1x12 ml	300	36 months	Serum	0,143 - 2,42 (g/l)
	GPTUR-H00/AU	2x60 ml	2x12 ml	600			
KIT for AU®400-640-2700-5400	GPTUR-B00/AU	1x60 ml	1x12 ml	240	36 months	Serum	0,143 - 2,42 (g/l)
	GPTUR-H00/AU	2x60 ml	2x12 ml	480			
KIT for AU®5800	GPTUR-B00/AU	1x60 ml	1x12 ml	400	36 months	Serum	0,143 - 2,42 (g/l)
	GPTUR-H00/AU	2x60 ml	2x12 ml	800			
KIT for DxC 700 AU®	GPTUR-B00/AU	1x60 ml	1x12 ml	400	36 months	Serum	0,143 - 2,42 (g/l)
	GPTUR-H00/AU	2x60 ml	2x12 ml	800			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA-1 ANTITRYPSIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	ATTUR-H00/AU	1x60 ml	1x12 ml	300	36 months	Serum & plasma	0.007 - 3.50 (g/l)
KIT for AU®400-640-2700-5400	ATTUR-H00/AU	1x60 ml	1x12 ml	240	36 months	Serum & plasma	0.007 - 3.50 (g/l)
KIT for AU®5800	ATTUR-H00/AU	1x60 ml	1x12 ml	400	36 months	Serum & plasma	0.007 - 3.50 (g/l)
KIT for DxC 700 AU®	ATTUR-H00/AU	1x60 ml	1x12 ml	400	36 months	Serum & plasma	0.007 - 3.50 (g/l)

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CATALOGUE 2025

ALPHA-1 MICROGLOBULIN (URINE)

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®400-480-640-680-2700-5400-5800	A1DGX-B00/AU	1 x 15 ml	1 x 15 ml	250	17 months	Serum & plasma	1,29 - 166,6 (mg/l)
	A1DGX-H00/AU	2 x 20 ml	2 x 20 ml	667			
	A1DGX-L00/AU	12 x 20 ml	12 x 20 ml	4000			
KIT for DxC 700 AU®	A1DGX-B00/AU	1 x 15 ml	1 x 15 ml	200	17 months	Serum & plasma	1,29 - 166,6 (mg/l)
	A1DGX-H00/AU	2 x 20 ml	2 x 20 ml	400			

CALIBRATORS

KIT 5 LEV	A1REK-000	5x1 ml	14 months
-----------	-----------	--------	-----------

CONTROLS

M	A1CON-002	1x2 ml	14 months
L	A1COS-002	1x2 ml	14 months

ALPHA-2 MACROGLOBULIN - ALPHA-2 MACROGLOBULIN URINARY

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	MGTUR-H00/AU	1x60 ml	1x12 ml	300	36 months	Serum & plasma or urine	0,026 - 4,96 (g/l) S/P 3,36 - 199 (mg/l) urine
	MGTUR-L00/AU	12x60 ml	12x12 ml	3600			
KIT for AU®400-640-2700-5400	MGTUR-H00/AU	1x60 ml	1x12 ml	240	36 months	Serum & plasma or urine	0,026 - 4,96 (g/l) S/P 3,36 - 199 (mg/l) urine
	MGTUR-L00/AU	12x60 ml	12x12 ml	2880			
KIT for AU®5800	MGTUR-H00/AU	1x60 ml	1x12 ml	400	36 months	Serum & plasma or urine	0,026 - 4,96 (g/l) S/P 3,36 - 199 (mg/l) urine
	MGTUR-L00/AU	12x60 ml	12x12 ml	4800			
KIT for DxC 700 AU®	MGTUR-H00/AU	1x60 ml	1x12 ml	400	36 months	Serum & plasma or urine	0,026 - 4,96 (g/l) S/P 3,36 - 199 (mg/l) urine
	MGTUR-L00/AU	12x60 ml	12x12 ml	4800			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months
URINARY LOW	MGCU1-002	1x2 ml	24 months
URINARY HIGH	MGCU2-002	1x2 ml	24 months

CATALOGUE 2025

ANTI STREPTOLYSIN O (ASO)

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	ASLTX-B00/AU	1x60 ml	1x6.5 ml	295	20 months	Serum & plasma	19.99 to 864.85 (IU/ml)
	ASLTX-H00/AU	2x60 ml	2x6.5 ml	591			
KIT for AU®400-640-2700-5400	ASLTX-B00/AU	1x60 ml	1x6.5 ml	241	20 months	Serum & plasma	19.99 to 864.85 (IU/ml)
	ASLTX-H00/AU	2x60 ml	2x6.5 ml	481			
KIT for AU®5800	ASLTX-B00/AU	1x60 ml	1x6.5 ml	406	20 months	Serum & plasma	19.99 to 864.85 (IU/ml)
	ASLTX-H00/AU	2x60 ml	2x6.5 ml	813			
KIT for DxC 700 AU®	ASLTX-B00/AU	1x60 ml	1x6.5 ml	406	20 months	Serum & plasma	19.99 to 864.85 (IU/ml)
	ASLTX-H00/AU	2x60 ml	2x6.5 ml	813			

CALIBRATORS

H	ASREK-000	1x1 or (1x5) ml	18 months
---	-----------	-----------------	-----------

CONTROLS

M	ASCON-002	1x2 ml	20 months
---	-----------	--------	-----------

APOLIPOPROTEIN A1

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	AATUR-B00/AU	1x50 ml	1x13 ml	245	36 months	Serum & plasma	0,029 - 4 (g/l)
	AATUR-H00/AU	2x50 ml	2x13 ml	491			
KIT for AU®400-640-2700-5400	AATUR-B00/AU	1x50 ml	1x13 ml	197	36 months	Serum & plasma	0,029 - 4 (g/l)
	AATUR-H00/AU	2x50 ml	2x13 ml	394			
KIT for AU®5800	AATUR-B00/AU	1x50 ml	1x13 ml	325	36 months	Serum & plasma	0,029 - 4 (g/l)
	AATUR-H00/AU	2x50 ml	2x13 ml	650			
KIT for DxC 700 AU®	AATUR-B00/AU	1x50 ml	1x13 ml	325	36 months	Serum & plasma	0,029 - 4 (g/l)
	AATUR-H00/AU	2x50 ml	2x13 ml	650			

CALIBRATORS

KIT 5 LEV	ABREK-000	5x0,5 ml	18 months
-----------	-----------	----------	-----------

CONTROLS

H	ABCOX-002	1x2 ml	18 months
M	ABCON-002	1x2 ml	18 months

CATALOGUE 2025

APOLIPOPROTEIN B

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	ABTUR-B00/AU	1x50 ml	1x10 ml	286	36 months	Serum & plasma	0.019 - 4 (g/l)
	ABTUR-H00/AU	2x50 ml	2x10 ml	571			
KIT for AU®400-640-2700-5400	ABTUR-B00/AU	1x50 ml	1x10 ml	200	36 months	Serum & plasma	0.019 - 4 (g/l)
	ABTUR-H00/AU	2x50 ml	2x10 ml	400			
KIT for AU®5800	ABTUR-B00/AU	1x50 ml	1x10 ml	333	36 months	Serum & plasma	0.019 - 4 (g/l)
	ABTUR-H00/AU	2x50 ml	2x10 ml	667			
KIT for DxC 700 AU®	ABTUR-B00/AU	1x50 ml	1x10 ml	333	36 months	Serum & plasma	0.019 - 4 (g/l)
	ABTUR-H00/AU	2x50 ml	2x10 ml	667			

CALIBRATORS

KIT 5 LEV	ABREK-000	5x0,5 ml	18 months
-----------	-----------	----------	-----------

CONTROLS

H	ABCOX-002	1x2 ml	18 months
M	ABCON-002	1x2 ml	18 months

BETA-2 MICROGLOBULIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	B2DGX-B00/AU	1x18 ml	1x8 ml	125	30 months	serum & plasma	0,132 to 95 (mg/l)
	B2DGX-H00/AU	2x44 ml	2x19 ml	594			
	B2DGX-L00/AU	8x44 ml	8x19 ml	2375			
KIT for AU®400-640-2700-5400	B2DGX-B00/AU	1x18 ml	1x8 ml	100	30 months	serum & plasma	0,132 to 95 (mg/l)
	B2DGX-H00/AU	2x44 ml	2x19 ml	475			
	B2DGX-L00/AU	8x44 ml	8x19 ml	1900			
KIT for AU®5800	B2DGX-B00/AU	1x18 ml	1x8 ml	167	30 months	serum & plasma	0,132 to 95 (mg/l)
	B2DGX-H00/AU	2x44 ml	2x19 ml	792			
	B2DGX-L00/AU	8x44 ml	8x19 ml	3167			
KIT for DxC 700 AU®	B2DGX-B00/AU	1x18 ml	1x8 ml	167	30 months	serum & plasma	0,132 to 95 (mg/l)
	B2DGX-H00/AU	2x44 ml	2x19 ml	792			
	B2DGX-L00/AU	8x44 ml	8x19 ml	3167			

CALIBRATORS

KIT 5 LEV	B2REK-000	5x1 ml	30 months
-----------	-----------	--------	-----------

CONTROLS

H	B2COX-001	1 x 1 ml	30 months
M	B2CON-003	1x3 ml	30 months
L	B2COS-003	1x3 ml	30 months

CATALOGUE 2025

BETA-2 MICROGLOBULIN (URINE)

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	BUDGX-B00/AU	1 x 28 ml	1 x 12 ml	125	30 months	Urine	0,055 to 9.3 (mg/l)
	BUDGX-H00/AU	2 x 28 ml	2 x 12 ml	594			
KIT for AU®400-640-2700-5400	BUDGX-B00/AU	1 x 28 ml	1 x 12 ml	100	30 months	Urine	0,055 to 9.3 (mg/l)
	BUDGX-H00/AU	2 x 28 ml	2 x 12 ml	475			
KIT for AU®5800	BUDGX-B00/AU	1 x 28 ml	1 x 12 ml	167	30 months	Urine	0,055 to 9.3 (mg/l)
	BUDGX-H00/AU	2 x 28 ml	2 x 12 ml	792			
KIT for DxC 700 AU®	BUDGX-B00/AU	1 x 28 ml	1 x 12 ml	300	30 months	Urine	0,055 to 9.3 (mg/l)
	BUDGX-H00/AU	2 x 28 ml	2 x 12 ml	600			

CALIBRATORS

KIT 5 LEV	BUREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

M	B2CON-003	1x3 ml	30 months
L	B2COS-003	1x3 ml	30 months

CALPROTECTIN FECAL - Colloidal Gold

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	CACOL-B00/AU	1x20 ml	1x8,2 ml	222	18 months	faeces	10,9 mg/kg - 1297,4 (mg/kg)
	CACOL-H00/AU	2x20 ml	2x8,2 ml	443			
	CACOL-H00/AU	4x20 ml	4x8,2 ml	886			
KIT for AU®400-640-2700-5400	CACOL-B00/AU	1x20 ml	1x8,2 ml	222	18 months	faeces	10,9 mg/kg - 1297,4 (mg/kg)
	CACOL-H00/AU	2x20 ml	2x8,2 ml	443			
	CACOL-H00/AU	4x20 ml	4x8,2 ml	886			
KIT for AU®5800	CACOL-B00/AU	1x20 ml	1x8,2 ml	222	18 months	faeces	10,9 mg/kg - 1297,4 (mg/kg)
	CACOL-H00/AU	2x20 ml	2x8,2 ml	443			
	CACOL-H00/AU	4x20 ml	4x8,2 ml	886			
KIT for DxC 700 AU®	CACOL-B00/AU	1x20 ml	1x8,2 ml	222	18 months	faeces	10,9 mg/kg - 1297,4 (mg/kg)
	CACOL-H00/AU	2x20 ml	2x8,2 ml	443			
	CACOL-H00/AU	4x20 ml	4x8,2 ml	886			

CALIBRATORS

KIT 5 LEV	CAREK-000	5x1 ml	12 months
-----------	-----------	--------	-----------

CONTROLS

H	CACOX-002	1x2 ml	12 months
M	CACON-002	1x2 ml	12 months
L	CACOS-002	1x2 ml	12 months

CATALOGUE 2025

CERULOPLASMIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	CETUR-B00/AU	1x30 ml	1x6 ml	250	36 months	Serum & plasma	9,48 - 1051,2 (mg/l)
	CETUR-H00/AU	1x60 ml	1x12 ml	500			
KIT for AU®400-640-2700-5400	CETUR-B00/AU	1x30 ml	1x6 ml	200	36 months	Serum & plasma	9,48 - 1051,2 (mg/l)
	CETUR-H00/AU	1x60 ml	1x12 ml	400			
KIT for AU®5800	CETUR-B00/AU	1x30 ml	1x6 ml	333	36 months	Serum & plasma	9,48 - 1051,2 (mg/l)
	CETUR-H00/AU	1x60 ml	1x12 ml	667			
KIT for DxC 700 AU®	CETUR-B00/AU	1x30 ml	1x6 ml	250	36 months	Serum & plasma	9,48 - 1051,2 (mg/l)
	CETUR-H00/AU	1x60 ml	1x12 ml	500			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

COMPLEMENT C3

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	C3TUR-B00/AU	1x50 ml	1x10 ml	400	36 months	Serum & plasma	0.015 - 4.31 (g/l)
	C3TUR-H00/AU	2x50 ml	2x10 ml	800			
KIT for AU®400-640-2700-5400	C3TUR-B00/AU	1x50 ml	1x10 ml	263	36 months	Serum & plasma	0.015 - 4.31 (g/l)
	C3TUR-H00/AU	2x50 ml	2x10 ml	526			
KIT for AU®5800	C3TUR-B00/AU	1x50 ml	1x10 ml	500	36 months	Serum & plasma	0.015 - 4.31 (g/l)
	C3TUR-H00/AU	2x50 ml	2x10 ml	1000			
KIT for DxC 700 AU®	C3TUR-B00/AU	1x50 ml	1x10 ml	400	36 months	Serum & plasma	0.015 - 4.31 (g/l)
	C3TUR-H00/AU	2x50 ml	2x10 ml	800			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CATALOGUE 2025

COMPLEMENT C4

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	C4TUR-B00/AU	1x50 ml	1x10 ml	400	36 months	Serum & plasma	0,002 - 0,825 (g/l)
	C4TUR-H00/AU	2x50 ml	2x10 ml	800			
KIT for AU®400-640-2700-5400	C4TUR-B00/AU	1x50 ml	1x10 ml	278	36 months	Serum & plasma	0,002 - 0,825 (g/l)
	C4TUR-H00/AU	2x50 ml	2x10 ml	556			
KIT for AU®5800	C4TUR-B00/AU	1x50 ml	1x10 ml	500	36 months	Serum & plasma	0,002 - 0,825 (g/l)
	C4TUR-H00/AU	2x50 ml	2x10 ml	1000			
KIT for DxC 700 AU®	C4TUR-B00/AU	1x50 ml	1x10 ml	400	36 months	Serum & plasma	0,002 - 0,825 (g/l)
	C4TUR-H00/AU	2x50 ml	2x10 ml	800			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CRP XL WIDE RANGE

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	CPTXL-B00/AU	1x60 ml	1x12 ml	120	36 months	Serum & plasma	0,607 - 425 (mg/l)
	CPTXL-H00/AU	2x60 ml	2x12 ml	923			
	CPTXL-L00/AU	12x60 ml	12x12 ml	5938			
KIT for AU®400-640-2700-5400	CPTXL-B00/AU	1x60 ml	1x12 ml	333	36 months	Serum & plasma	0,607 - 425 (mg/l)
	CPTXL-H00/AU	2x60 ml	2x12 ml	667			
	CPTXL-L00/AU	12x60 ml	12x12 ml	4000			
KIT for AU®5800	CPTXL-B00/AU	1x60 ml	1x12 ml	556	36 months	Serum & plasma	0,607 - 425 (mg/l)
	CPTXL-H00/AU	2x60 ml	2x12 ml	1111			
	CPTXL-L00/AU	12x60 ml	12x12 ml	6667			
KIT for DxC 700 AU®	CPTXL-B00/AU	1x60 ml	1x12 ml	462	36 months	Serum & plasma	0,607 - 425 (mg/l)
	CPTXL-H00/AU	2x60 ml	2x12 ml	923			
	CPTXL-L00/AU	12x60 ml	12x12 ml	5538			

CALIBRATORS

KIT 5 LEV	CPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	24 months
M	CPCON-002	1x2 ml	24 months
L	CXCON-002	1x2 ml	24 months

CATALOGUE 2025

CRP GOLD LATEX

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	CPDGX-B00/AU	1x20 ml	1x20 ml	200	24 months	Serum & plasma	0,031 to 300 (mg/l)
	CPDGX-H00/AU	2x50 ml	2x50 ml	1000			
	CPDGX-L00/AU	8x50 ml	8x50 ml	4000			
KIT for AU®400-640-2700-5400	CPDGX-B00/AU	1x20 ml	1x20 ml	200	24 months	Serum & plasma	0,031 to 300 (mg/l)
	CPDGX-H00/AU	2x50 ml	2x50 ml	1000			
	CPDGX-L00/AU	8x50 ml	8x50 ml	4000			
KIT for AU®5800	CPDGX-B00/AU	1x20 ml	1x20 ml	333	24 months	Serum & plasma	0,031 to 300 (mg/l)
	CPDGX-H00/AU	2x50 ml	2x50 ml	1667			
	CPDGX-L00/AU	8x50 ml	8x50 ml	6667			
KIT for DxC 700 AU®	CPDGX-B00/AU	1x20 ml	1x20 ml	333	24 months	Serum & plasma	0,031 to 300 (mg/l)
	CPDGX-H00/AU	2x50 ml	2x50 ml	1667			
	CPDGX-L00/AU	8x50 ml	8x50 ml	6667			

CALIBRATORS

KIT 5 LEV	CPRLK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	24 months
M	CPCON-002	1x2 ml	24 months
L	CXCXN-002	1x2 ml	24 months

CYSTATIN C - Colloidal Gold

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	CYCOL-B00/AU	1x20 ml	1x7.5 ml	125	24 months	Serum & plasma	0.15 to 9.355 (mg/l)
	CYCOL-H00/AU	2x50 ml	2x15 ml	500			
	CYCOL-L00/AU	12x50 ml	12x15 ml	3000			
KIT for AU®400-640-2700-5400	CYCOL-B00/AU	1x20 ml	1x7.5 ml	100	24 months	Serum & plasma	0.15 to 9.355 (mg/l)
	CYCOL-H00/AU	2x50 ml	2x15 ml	400			
	CYCOL-L00/AU	12x50 ml	12x15 ml	2400			
KIT for AU®5800	CYCOL-B00/AU	1x20 ml	1x7.5 ml	167	24 months	Serum & plasma	0.15 to 9.355 (mg/l)
	CYCOL-H00/AU	2x50 ml	2x15 ml	667			
	CYCOL-L00/AU	12x50 ml	12x15 ml	4000			
KIT for DxC 700 AU®	CYCOL-B00/AU	1x20 ml	1x7.5 ml	167	24 months	Serum & plasma	0.15 to 9.355 (mg/l)
	CYCOL-H00/AU	2x50 ml	2x15 ml	667			
	CYCOL-L00/AU	12x50 ml	12x15 ml	4000			

CALIBRATORS

KIT 6 LEV	CYREK-000	6x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	CYCOX-002	1x2 ml	24 months
M	CYCON-002	1x2 ml	24 months
L	CYCOS-002	1x2 ml	24 months

CATALOGUE 2025

FERRITIN - Colloidal Gold

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	FNCOL-B00/AU	1x50 ml	1x9 ml	360	12 months	Serum	6,344 - 500 (µg/l)
	FNCOL-H00/AU	2x50 ml	2x9 ml	720			
KIT for AU®400-640-2700-5400	FNCOL-B00/AU	1x50 ml	1x9 ml	300	12 months	Serum	6,344 - 500 (µg/l)
	FNCOL-H00/AU	2x50 ml	2x9 ml	600			
KIT for AU®5800	FNCOL-B00/AU	1x50 ml	1x9 ml	450	12 months	Serum	6,344 - 500 (µg/l)
	FNCOL-H00/AU	2x50 ml	2x9 ml	900			
KIT for DxC 700 AU®	FNCOL-B00/AU	1x50 ml	1x9 ml	450	12 months	Serum	6,344 - 500 (µg/l)
	FNCOL-H00/AU	2x50 ml	2x9 ml	900			

CALIBRATORS

KIT 5 LEV	FNRGK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

M	FNCON-002	1x2 ml	18 months
---	-----------	--------	-----------

HAPTOGLOBIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	HATUR-B00/AU	1x30 ml	1x6 ml	150	36 months	Serum & plasma	0.002 - 3.86 (g/l)
	HATUR-H00/AU	1x60 ml	1x12 ml	300			
KIT for AU®400-640-2700-5400	HATUR-B00/AU	1x30 ml	1x6 ml	120	36 months	Serum & plasma	0.002 - 3.86 (g/l)
	HATUR-H00/AU	1x60 ml	1x12 ml	240			
KIT for AU®5800	HATUR-B00/AU	1x30 ml	1x6 ml	200	36 months	Serum & plasma	0.002 - 3.86 (g/l)
	HATUR-H00/AU	1x60 ml	1x12 ml	400			
KIT for DxC 700 AU®	HATUR-B00/AU	1x30 ml	1x6 ml	200	36 months	Serum & plasma	0.002 - 3.86 (g/l)
	HATUR-H00/AU	1x60 ml	1x12 ml	400			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CATALOGUE 2025

IMMUNOGLOBULINS A

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	IATUR-B00/AU	1x50 ml	1x10 ml	370	36 months	Serum	0.02 - 8 (g/l)
	IATUR-H00/AU	2x50 ml	2x10 ml	741			
KIT for AU®400-640-2700-5400	IATUR-B00/AU	1x50 ml	1x10 ml	299	36 months	Serum	0.02 - 8 (g/l)
	IATUR-H00/AU	2x50 ml	2x10 ml	599			
KIT for AU®5800	IATUR-B00/AU	1x50 ml	1x10 ml	500	36 months	Serum	0.02 - 8 (g/l)
	IATUR-H00/AU	2x50 ml	2x10 ml	1000			
KIT for DxC 700 AU®	IATUR-B00/AU	1x50 ml	1x10 ml	370	36 months	Serum	0.02 - 8 (g/l)
	IATUR-H00/AU	2x50 ml	2x10 ml	741			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS E

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	IEDGX-B00/AU	1x20 ml	1x10 ml	217	20 months	Serum & plasma	21,0 - 1436 UI/ml
	IEDGX-H00/AU	2x30 ml	2x15 ml	652			
KIT for AU®400-640-2700-5400	IEDGX-B00/AU	1x20 ml	1x10 ml	167	20 months	Serum & plasma	21,0 - 1436 UI/ml
	IEDGX-H00/AU	2x30 ml	2x15 ml	500			
KIT for AU®5800	IEDGX-B00/AU	1x20 ml	1x10 ml	250	20 months	Serum & plasma	21,0 - 1436 UI/ml
	IEDGX-H00/AU	2x30 ml	2x15 ml	750			
KIT for DxC 700 AU®	IEDGX-B00/AU	1x20 ml	1x10 ml	250	20 months	Serum & plasma	21,0 - 1436 UI/ml
	IEDGX-H00/AU	2x30 ml	2x15 ml	750			

CALIBRATORS

KIT 5 LEV	IEREK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

M	IECON-003	1x3 ml	18 months
L	IECOS-003	1x3 ml	18 months

CATALOGUE 2025

IMMUNOGLOBULINS G

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	IGTUR-B00/AU	1x50 ml	1x18 ml	353	36 months	Serum	0.125 - 35 (g/l)
	IGTUR-H00/AU	2x50 ml	2x18 ml	706			
KIT for AU®400-640-2700-5400	IGTUR-B00/AU	1x50 ml	1x18 ml	286	36 months	Serum	0.125 - 35 (g/l)
	IGTUR-H00/AU	2x50 ml	2x18 ml	571			
KIT for AU®5800	IGTUR-B00/AU	1x50 ml	1x18 ml	474	36 months	Serum	0.125 - 35 (g/l)
	IGTUR-H00/AU	2x50 ml	2x18 ml	947			
KIT for DxC 700 AU®	IGTUR-B00/AU	1x50 ml	1x18 ml	353	36 months	Serum	0.125 - 35 (g/l)
	IGTUR-H00/AU	2x50 ml	2x18 ml	706			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS G (CSF)

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	IGLCR-B00/AU	1 x 55 ml	1 x 10 ml	500	36 months	CSF	3,375 - 249,78 (mg/l)
	IGLCR-H00/AU	2 x 55 ml	2 x 10 ml	1000			
KIT for AU®400-640-2700-5400	IGLCR-B00/AU	1 x 55 ml	1 x 10 ml	400	36 months	CSF	3,375 - 249,78 (mg/l)
	IGLCR-H00/AU	2 x 55 ml	2 x 10 ml	800			
KIT for AU®5800	IGLCR-B00/AU	1 x 55 ml	1 x 10 ml	667	36 months	CSF	3,375 - 249,78 (mg/l)
	IGLCR-H00/AU	2 x 55 ml	2 x 10 ml	1333			
KIT for DxC 700 AU®	IGLCR-B00/AU	1 x 55 ml	1 x 10 ml	500	36 months	CSF	3,375 - 249,78 (mg/l)
	IGLCR-H00/AU	2 x 55 ml	2 x 10 ml	1000			

CALIBRATORS

KIT 5 LEV	IGREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	IGCOX-002	1x2 ml	24 months
M	IGCON-002	1x2 ml	24 months
L	IGCOS-002	1x2 ml	24 months

CATALOGUE 2025

IMMUNOGLOBULINS M

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	IMTUR-B00/AU	1x50 ml	1x10 ml	370	36 months	Serum	0.038 - 7 (g/l)
	IMTUR-H00/AU	2x50 ml	2x10 ml	741			
KIT for AU®400-640-2700-5400	IMTUR-B00/AU	1x50 ml	1x10 ml	299	36 months	Serum	0.038 - 7 (g/l)
	IMTUR-H00/AU	2x50 ml	2x10 ml	599			
KIT for AU®5800	IMTUR-B00/AU	1x50 ml	1x10 ml	500	36 months	Serum	0.038 - 7 (g/l)
	IMTUR-H00/AU	2x50 ml	2x10 ml	1000			
KIT for DxC 700 AU®	IMTUR-B00/AU	1x50 ml	1x10 ml	370	36 months	Serum	0.038 - 7 (g/l)
	IMTUR-H00/AU	2x50 ml	2x10 ml	741			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

LIPOPROTEIN (A)

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	LPTUR2-H00/AU	1x50 ml	1x6 ml	400	12 months	Serum & plasma	0,0235 - 0,946 (g/l)
KIT for AU®400-640-2700-5400	LPTUR2-H00/AU	1x50 ml	1x6 ml	240	12 months	Serum & plasma	0,0235 - 0,946 (g/l)
KIT for AU®5800	LPTUR2-H00/AU	1x50 ml	1x6 ml	462	12 months	Serum & plasma	0,0235 - 0,946 (g/l)
KIT for DxC 700 AU®	LPTUR2-H00/AU	1x50 ml	1x6 ml	400	12 months	Serum & plasma	0,0235 - 0,946 (g/l)

CALIBRATORS

KIT 4 LEV	LPREK-000	4x0,5 ml	18 months
-----------	-----------	----------	-----------

CONTROLS

M	LPCON-002	1x2 ml	24 months
---	-----------	--------	-----------

CATALOGUE 2025

MICROALBUMIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	MATUR-B00/AU	1x50 ml	1x10 ml	417	36 months	Urine	1,86 - 400 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	833			
	MATUR-L00/AU	12x50 ml	12x10 ml	5000			
KIT for AU®400-640-2700-5400	MATUR-B00/AU	1x50 ml	1x10 ml	333	36 months	Urine	1,86 - 400 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	667			
	MATUR-L00/AU	12x50 ml	12x10 ml	4000			
KIT for AU®5800	MATUR-B00/AU	1x50 ml	1x10 ml	556	36 months	Urine	1,86 - 400 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	1111			
	MATUR-L00/AU	12x50 ml	12x10 ml	6667			
KIT for DxC 700 AU®	MATUR-B00/AU	1x50 ml	1x10 ml	417	36 months	Urine	1,86 - 400 (mg/l)
	MATUR-H00/AU	2x50 ml	2x10 ml	833			
	MATUR-L00/AU	12x50 ml	12x10 ml	5000			

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MACOX-002	1x2 ml	24 months
M	MACON-002	1x2 ml	24 months
L	MACOS-002	1x2 ml	24 months

PREALBUMIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	PATUR-B00/AU	1x50 ml	1x10 ml	370	36 months	Serum & plasma	6.53 - 1000 (mg/l)
	PATUR-H00/AU	2x50 ml	2x10 ml	741			
KIT for AU®400-640-2700-5400	PATUR-B00/AU	1x50 ml	1x10 ml	299	36 months	Serum & plasma	6.53 - 1000 (mg/l)
	PATUR-H00/AU	2x50 ml	2x10 ml	599			
KIT for AU®5800	PATUR-B00/AU	1x50 ml	1x10 ml	500	36 months	Serum & plasma	6.53 - 1000 (mg/l)
	PATUR-H00/AU	2x50 ml	2x10 ml	1000			
KIT for DxC 700 AU®	PATUR-B00/AU	1x50 ml	1x10 ml	370	36 months	Serum & plasma	6.53 - 1000 (mg/l)
	PATUR-H00/AU	2x50 ml	2x10 ml	741			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CATALOGUE 2025

RETINOL BINDING PROTEIN (RBP)

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	RBTUR-B00/AU	1 x 30,5 ml	1 x 10 ml	263	36 months	Serum & plasma	7,81 - 141 (mg/l)
	RBTUR-H00/AU	2 x 30,5ml	2 x 10 ml	526			
	RBTUR-L00/AU	12 x 30,5 ml	12 x 10 ml	3158			
KIT for AU®400-640-2700-5400	RBTUR-B00/AU	1 x 30,5 ml	1 x 10 ml	263	36 months	Serum & plasma	7,81 - 141 (mg/l)
	RBTUR-H00/AU	2 x 30,5ml	2 x 10 ml	526			
	RBTUR-L00/AU	12 x 30,5 ml	12 x 10 ml	3158			
KIT for AU®5800	RBTUR-B00/AU	1 x 30,5 ml	1 x 10 ml	400	36 months	Serum & plasma	7,81 - 141 (mg/l)
	RBTUR-H00/AU	2 x 30,5ml	2 x 10 ml	800			
	RBTUR-L00/AU	12 x 30,5 ml	12 x 10 ml	4800			
KIT for DxC 700 AU®	RBTUR-B00/AU	1 x 30,5 ml	1 x 10 ml	333	36 months	Serum & plasma	7,81 - 141 (mg/l)
	RBTUR-H00/AU	2 x 30,5ml	2 x 10 ml	667			
	RBTUR-L00/AU	12 x 30,5 ml	12 x 10 ml	4000			

CALIBRATORS

KIT 5 LEV	RBREK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

H	RBCOX-002	1x2 ml	18 months
M	RBCON-002	1x2 ml	18 months
L	RBCOS-002	1x2 ml	18 months

RHEUMATOID FACTOR (RF) - Colloidal Gold

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	RFCOL-B00/AU	1x50 ml	1x20 ml	230	18 months	Serum & plasma	2.228 - 120 (kUI/l)
	RFCOL-H00/AU	2x50 ml	2x20 ml	460			
KIT for AU®400-640-2700-5400	RFCOL-B00/AU	1x50 ml	1x20 ml	200	18 months	Serum & plasma	2.228 - 120 (kUI/l)
	RFCOL-H00/AU	2x50 ml	2x20 ml	400			
KIT for AU®5800	RFCOL-B00/AU	1x50 ml	1x20 ml	333	18 months	Serum & plasma	2.228 - 120 (kUI/l)
	RFCOL-H00/AU	2x50 ml	2x20 ml	667			
KIT for DxC 700 AU®	RFCOL-B00/AU	1x50 ml	1x20 ml	333	18 months	Serum & plasma	2.228 - 120 (kUI/l)
	RFCOL-H00/AU	2x50 ml	2x20 ml	667			

CALIBRATORS

KIT 5 LEV	RFREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	RFCON-002	1x2 ml	24 months
M	RFCON-002	1x2 ml	24 months
L	RFCON-002	1x2 ml	24 months

CATALOGUE 2025

TRANSFERRIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT for AU®480-680	TFTUR-H00/AU	1x60 ml	1x16,5 ml	236	36 months	Serum & plasma	0.076 - 8 (g/l)
	TFTUR-L00/AU	7x60 ml	7x16,5 ml	1650			
KIT for AU®400-640-2700-5400	TFTUR-H00/AU	1x60 ml	1x16,5 ml	295	36 months	Serum & plasma	0.076 - 8 (g/l)
	TFTUR-L00/AU	7x60 ml	7x16,5 ml	2063			
KIT for AU®5800	TFTUR-H00/AU	1x60 ml	1x16,5 ml	393	36 months	Serum & plasma	0.076 - 8 (g/l)
	TFTUR-L00/AU	7x60 ml	7x16,5 ml	2750			
KIT for AU®DXC	TFTUR-H00/AU	1x60 ml	1x16,5 ml	295	36 months	Serum & plasma	0.076 - 8 (g/l)
	TFTUR-L00/AU	7x60 ml	7x16,5 ml	2063			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

Your expert in specific proteins



-
- (*) Tests :** Amount of tests possible by kit. May vary from an instrument to another one.
- (**) Shelf Life :** Diagam provides to distributors a minimum of *12 months*. The time period written in this document shows the maximum shelf life time after manufacturing.
- (***) Sample :** Diagam may provide some applications for other samples, such as CSF, Urine, Veterany samples etc.... Please contact your distributor or us directly at mail@diagam.com or dial +32 68 55 14 82.
- (****) Analytical range :** For general information purposes only. More details can be found in the package insert (instructions for use). Please contact your distributor or us directly at mail@diagam.com or dial +32 68 55 14 82 in order to get the updated data.



Medical device for in vitro diagnostics.

**Contact
information:**

For orders, invoicing and technical queries, we kindly recommend you to use our general e-mail address that is mail@diagam.com.



DIAGAM SA, Rue du Parc Industriel, 7822 Ghislenghien, BELGIUM

© Diagam. 2025. All rights reserved. All product names, registered trademarks, company names in this document remain the property of their respective owners. For general information purposes only. More details can be found in the package insert (instructions for use) and labelling. Not for the US market except for research use only.