



ARCHITECT LINE



CATALOGUE 2025

ALBUMIN						
REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>			
	ALTUR-B00/ARC	1 x 40 ml R1 + 1 x 10 ml R2	169	36 months	serum or plasma	0.40 to 121.71 (g/l)
	ALTUR-H00/ARC	1 x 80 ml R1 + 1 x 20 ml R2	372			
	ALTUR-L00/ARC	7 x 80 ml R1 + 7 x 20 ml R2	2604			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALBUMIN CSF						
REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>			
	MATUR-B00/ARC	1 x 22,5 ml R1 + 1 x 6,1 ml R2	100	36 months	cerebrospinal fluid	10.86 to 494.33 (mg/l)
	MATUR-H00/ARC	1 x 74 ml R1 + 1 x 17,5 ml R2	361			
	MATUR-L00/ARC	7 x 74 ml R1 + 7 x 17,5 ml R2	2525			

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MACOZ-002	1 x 2 ml	24 months
H	MACOX-002	1x2 ml	24 months
M	MACON-002	1x2 ml	24 months
L	MACOS-002	1x2 ml	24 months

ALPHA-1 GLYCOPROTEIN ACIDE/OROSOMUCOID						
REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>			
	GPTUR-B00/ARC	1 x 37,5 ml R1 + 1 x 7,5 ml R2	97	36 months	serum	0.143 to 2.42 (g/l)
	GPTUR-H00/ARC	2 x 75 ml R1 + 2 x 15 ml R2	443			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CATALOGUE 2025

ALPHA-1 ANTITRYPSIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	ATTUR-B00/ARC	1 x 45 ml R1 + 1 x 9 ml R2	122	24 months	serum or plasma	0.007 to 3.5 (g/l)
	ATTUR-H00/ARC	2 x 75 ml R1 + 2 x 15 ml R2	443			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA-1 MICROGLOBULIN (URINE)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	A1DGX-B00/ARC	1 x 15 ml R1 + 1 x 15 ml R2	150	17 months	urine	1,29 to 166,6 (mg/l)
	A1DGX-H00/ARC	1 x 55 ml R1 + 1 x 55 ml R2	456			

CALIBRATORS

KIT 5 LEV	A1REK-000	5x1 ml	14 months
-----------	-----------	--------	-----------

CONTROLS

M	A1CON-002	1x2 ml	14 months
L	A1COS-002	1x2 ml	14 months

ALPHA-2 MACROGLOBULIN - ALPHA-2 MACROGLOBULIN URINARY

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2			
KIT	MGTUR-H00/ARC	1 x 55 ml R1 + 1 x 11 ml R2	155	36 months	serum/plasma or urine	0.026 to 4.96 (g/l) serum 3.36 to 199 (mg/l) urine
	MGTUR-L00/ARC	10 x 55 ml R1 + 10 x 11 ml R2	1550			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months
URINARY LOW	MGCU1-002	1x2 ml	24 months
URINARY HIGH	MGCU2-002	1x2 ml	24 months

CATALOGUE 2025

ANTI STREPTOLYSIN O (ASO)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	ASLTX-B00/ARC	R1 1 x 70 ml R1 + 1 x 8 ml R2	296	20 months	serum or plasma	19.99 to 864.85 (IU/ml)
	ASLTX-H00/ARC	R2 2 x 87 ml R1 + 2 x 10 ml R2	741			

CALIBRATORS

H	ASREK-000	5x1 ml	18 months
---	-----------	--------	-----------

CONTROLS

M	ASCON-002	1x2 ml	12 months
---	-----------	--------	-----------

APOLIPOPROTEIN A1

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	AATUR-B00/ARC	R1 1 x 38 ml R1 + 1 x 10 ml R2	107	36 months	serum or plasma	0.029 to 4.0 (g/l)
	AATUR-H00/ARC	R2 2 x 75 ml R1 + 2 x 20 ml R2	474			

CALIBRATORS

KIT 5 LEV	ABREK-000	5x0.5 ml	18 months
-----------	-----------	----------	-----------

CONTROLS

H	ABCOX-002	1x2 ml	18 months
M	ABCON-002	1x2 ml	18 months

APOLIPOPROTEIN B

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	ABTUR-B00/ARC	R1 1 x 41 ml R1 + 1 x 8 ml R2	151	36 months	serum or plasma	0.019 to 4.0 (g/l)
	ABTUR-H00/ARC	R2 2 x 67 ml R1 + 2 x 13 ml R2	543			

CALIBRATORS

KIT 5 LEV	ABREK-000	5x0.5 ml	18 months
-----------	-----------	----------	-----------

CONTROLS

H	ABCOX-002	1x2 ml	18 months
M	ABCON-002	1x2 ml	18 months

CATALOGUE 2025

BETA 2 MICROGLOBULIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>			
	B2DGX-B00/ARC	1 x 27 ml R1 + 1 x 12 ml R2	111	30 months	serum or plasma	0,132 to 95 (mg/l)
	B2DGX-H00/ARC	2 x 45 ml R1 + 2 x 20 ml R2	396			
	B2DGX-L00/ARC	10 x 45 ml R1 + 10 x 20 ml R2	1981			

CALIBRATORS

KIT 5 LEV	B2REK-000	5x1 ml	30 months
-----------	-----------	--------	-----------

CONTROLS

H	B2COX-001	1 x 1 ml	30 months
M	B2CON-003	1x3 ml	30 months
L	B2COS-003	1x3 ml	30 months

BETA 2 MICROGLOBULIN (URINE)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>			
	BUDGX-B00/ARC	1 x 23 ml R1 + 1 x 10 ml R2	117	30 months	Urine	0.06 to 9.3 (mg/l)
	BUDGX-H00/ARC	2 x 23 ml R1 + 2 x 10 ml R2	234			

CALIBRATORS

KIT 5 LEV	BUREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

M	B2CON-003	1x3 ml	30 months
L	B2COS-003	1x3 ml	30 months

CALPROTECTIN FECAL - Colloidal Gold

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>			
	CACOL-B00/ARC	1 x 18 ml R1 + 1 x 7,5 ml R2	108	18 months	faeces	10.9 mg/kg – 1297.4 (mg/kg)
	CACOL-H00/ARC	3 x 18 ml R1 + 3 x 7,5 ml R2	325			
	CACOL-L00/ARC	7 x 18 ml R1 + 7 x 7,5 ml R2	759			

CALIBRATORS

KIT 5 LEV	CAREK-000	5x1 ml	12 months
-----------	-----------	--------	-----------

CONTROLS

H	CACOX-002	1x2 ml	12 months
M	CACON-002	1x2 ml	12 months
L	CACOS-002	1x2 ml	12 months

CATALOGUE 2025

CERULOPLASMIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	CETUR-B00/ARC	R1 1 x 40 ml R1 + 1 x 8 ml R2	128	36 months	serum or plasma	9.48 to 1051.20 (mg/l)
	CETUR-H00/ARC	R2 2 x 60 ml R1 + 2 x 12 ml R2	419			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

COMPLEMENT C3

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	C3TUR-B00/ARC	R1 1 x 41 ml R1 + 1 x 8 ml R2	137	36 months	serum or plasma	0.015 to 4.31 (g/l)
	C3TUR-H00/ARC	R2 2 x 56 ml R1 + 2 x 11 ml R2	405			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

COMPLEMENT C4

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
KIT	C4TUR-B00/ARC	R1 1 x 40 ml R1 + 1 x 8 ml R2	128	36 months	serum or plasma	0.002 to 0.825 (g/l)
	C4TUR-H00/ARC	R2 2 x 60 ml R1 + 2 x 12 ml R2	419			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CATALOGUE 2025

CRP GOLD LATEX

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>			
	CPDGX-B00/ARC	1 x 15 ml R1 + 1 x 15 ml R2	105	24 months	serum or plasma	0.031 to 300 (mg/l)
	CPDGX-H00/ARC	1 x 55 ml R1 + 1 x 55 ml R2	456			
	CPDGX-L00/ARC	10 x 55 ml R1 + 10 x 55 ml R2	4561			

CALIBRATORS

KIT 5 LEV	CPRLK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	24 months
M	CPCON-002	1x2 ml	24 months
L	CXCON-002	1x2 ml	24 months

CRP XL WIDE RANGE

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>			
	CPTXL-B00/ARC	1 x 55,5 ml R1 + 1 x 8 ml R2	140	36 months	serum or plasma	0.607 to 425 (mg/l)
	CPTXL-H00/ARC	2 x 85 ml R1 + 2 x 12 ml R2	459			
	CPTXL-L00/ARC	7 x 85 ml R1 + 7 x 12 ml R2	1607			

CALIBRATORS

KIT 5 LEV	CPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	24 months
M	CPCON-002	1x2 ml	24 months
L	CXCON-002	1x2 ml	24 months

CYSTATIN C - Colloidal Gold

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>			
	CYCOL-B00/ARC	1 x 27 ml R1 + 1 x 11 ml R2	107	24 months	serum or plasma	0.15 to 9.355 (mg/l)
	CYCOL-H00/ARC	2 x 48 ml R1 + 2 x 20 ml R2	421			
	CYCOL-L00/ARC	10 x 48 ml R1 + 10 x 20 ml R2	2103			

CALIBRATORS

KIT 5 LEV	CYREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	CYCOX-002	1X2 ml	24 months
M	CYCON-002	1x2 ml	24 months
L	CYCOS-002	1x2 ml	24 months

CATALOGUE 2025

FERRITIN - Colloidal Gold

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
KIT	FNCOL-B00/ARC FNCOL-H00/ARC	R1 1 x 40 ml R1 + 1 x 8 ml R2	R2 128	12 months	serum	6.344 to 500 (µg/l)
		2 x 55 ml R1 + 2 x 11 ml R2	378			

CALIBRATORS

KIT 5 LEV	FNRGK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

M	FNCON-002	1x2 ml	18 months
---	-----------	--------	-----------

HAPTOGLOBIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	HATUR-B00/ARC HATUR-H00/ARC	R1 1 x 37,5 ml R1 + 1 x 7,5 ml R2	R2 97	36 months	serum or plasma	0.002 to 3.86 (g/l)
		2 x 75 ml R1 + 2 x 15 ml R2	443			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS A

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
KIT	IATUR-B00/ARC IATUR-H00/ARC	R1 1 x 38 ml R1 + 1 x 8 ml R2	R2 144	36 months	serum	0.02 to 8.0 (g/l)
		2 x 54 ml R1 + 2 x 11 ml R2	425			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS G

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
KIT	IGTUR-B00/ARC IGTUR-H00/ARC	R1 1 x 25 ml R1 + 1 x 9,5 ml R2	R2 104	36 months	serum	0.125 to 35 (g/l)
		2 x 50 ml R1 + 2 x 19 ml R2	461			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CATALOGUE 2025

IMMUNOGLOBULINS G (CSF)

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
KIT	IGLCR-B00/ARC	R1 1 x 18,5 ml R1 + 1 x 5 ml R2	100	36 months	cerebrospinal fluid	3,375 to 249,78 (mg/l)
	IGLCR-H00/ARC	R2 1 x 74 ml R1 + 1 x 17 ml R2	464			

CALIBRATORS

KIT 5 LEV	IGREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	IGCOX-002	1x2 ml	24 months
M	IGCON-002	1x2 ml	24 months
L	IGCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS M

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
KIT	IMTUR-B00/ARC	R1 1 x 38 ml R1 + 1 x 8 ml R2	144	36 months	serum	0.038 to 7.0 (g/l)
	IMTUR-H00/ARC	R2 2 x 49,5 ml R1 + 2 x 10,5 ml R2	402			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

LIPOPROTEIN (A)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	LPTUR2-B00/ARC	R1 1 x 64 ml R1 + 1 x 8 ml R2	191	12 months	serum or plasma	0.0235 to 0.946 (g/l)
	LPTUR2-H00/ARC	R2 2 x 80 ml R1 + 2 x 10 ml R2	503			

CALIBRATORS

KIT 4 LEV	LPREK-000	4x0,5 ml	18 months
-----------	-----------	----------	-----------

CONTROLS

M	LPCON-002	1x2 ml	24 months
---	-----------	--------	-----------

CATALOGUE 2025

MICROALBUMIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>	36 months	urine	1.86 to 400 (mg/l)
	MATUR-B00/ARC	1 x 22,5 ml R1 + 1 x 6,1 ml R2	100			
	MATUR-H00/ARC	1 x 74 ml R1 + 1 x 17,5 ml R2	361			
	MATUR-L00/ARC	7 x 74 ml R1 + 7 x 17,5 ml R2	2525			

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MACOX-002	1x2 ml	24 months
M	MACON-002	1x2 ml	24 months
L	MACOS-002	1x2 ml	24 months

PREALBUMIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>	36 months	serum/plasma	6.53 to 1000 (mg/l)
	PATUR-H00/ARC	1 x 55 ml R1 + 1 x 11 ml R2	155			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

RETINOL BINDING PROTEIN (RBP)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>	36 months	serum or plasma	7,81 to 141 (mg/l)
	RBTUR-B00/ARC	1 x 31 ml R1 + 1 x 10 ml R2	138			
	RBTUR-H00/ARC	2 x 31 ml R1 + 2 x 10 ml R2	277			
	RBTUR-L00/ARC	10 x 31 ml R1 + 10 x 10 ml R2	1383			

CALIBRATORS

KIT 5 LEV	RBREK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

H	RBCOX-002	1x2 ml	18 months
M	RBCON-002	1x2 ml	18 months
L	RBCOS-002	1x2 ml	18 months

CATALOGUE 2025

RHEUMATOID FACTOR (RF) - Colloidal Gold

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>			
	RFCOL-B00/ARC	1 x 30 ml R1 + 1 x 20 ml R2	122	18 months	serum and plasma	2.228 to 120 (kIU/l)
	RFCOL-H00/ARC	2 x 60 ml R1 + 2 x 40 ml R2	512			

CALIBRATORS

KIT 5 LEV	RFREK-000	5x1 ml		24 months
-----------	-----------	--------	--	-----------

CONTROLS

H	RFCOX-002	1x2 ml		24 months
M	RFCON-002	1x2 ml		24 months
L	RFCOS-002	1x2 ml		24 months

TRANSFERRIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>			
	TFTUR-B00/ARC	1 x 37,5 ml R1 + 1 x 10,5 ml R2	108	36 months	serum and plasma	0.076 to 8.0 (g/l)
	TFTUR-H00/ARC	2 x 72 ml R1 + 2 x 20 ml R2	449			
	TFTUR-L00/ARC	7 x 72 ml R1 + 7 x 20 ml R2	1570			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml		24 months
-----------	-----------	--------	--	-----------

CONTROLS

H	MPCOX-002	1x2 ml		24 months
M	MPCON-002	1x2 ml		24 months
L	MPCOS-002	1x2 ml		24 months

Please contact us for more information on the following tests:

IMMUNOGLOBULIN E

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
----------	-----------	---------	--------	--------------	-----------	----------------------

Your expert in specific proteins



(*) Tests : Amount of real tests. May vary from a analyzer to another one.

() Shelf Life :** Diagam provides to distributors a minimum of **12 months**. The time period written in this document shows the maximum shelf life time after manufacturing.

(*) Sample :** Diagam may provide some applications for other samples, such as CSF, Urin, Veterany samples etc.... Please contact your distributor or us directly at mail@diagam.com or dial **+32 68 55 14 82**.

(**) Analytical range :** For general information purposes only. More details can be found in the package insert (instructions for use). Please contact your distributor or us directly at mail@diagam.com or dial **+32 68 55 14 82** in order to get the updated data.



Medical device for in vitro diagnosis.



DIAGAM SA, *Rue du Parc Industriel 40, 7822 Ghislenghien, BELGIUM*



The products fulfill the requirements of European directive 98/79 EC concerning the medical devices of in vitro diagnosis.