



ATELLICA® LINE



CATALOGUE 2025

CATALOGUE

ALBUMIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>	36 months	serum or plasma	0.815 to 124.94 (g/l)
	ALTUR-B00/ATE	1 x 1 x 13,2 ml R1 + 1 x 1 x 5,2 ml R2	101			
	ALTUR-H00/ATE	1 x 2 x 24,5 ml R1 + 1 x 2 x 8,7 ml R2	395			
	ALTUR-L00/ATE	5 x 2 x 24,5 ml R1 + 5 x 2 x 8,7 ml R2	1973			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALBUMIN CSF

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>	36 months	cerebrospinal fluid	10.86 to 494.33 (mg/l)
	MATUR-B00/ATE	1 x 1 x 13,1 ml R1 + 1 x 1 x 4,6 ml R2	99			
	MATUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 7,5 ml R2	385			

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

H	MACOZ-002	1x2 ml	24 months
H	MACOX-002	1x2 ml	24 months
M	MACON-002	1x2 ml	24 months
L	MACOS-002	1x2 ml	24 months

ALPHA-1 ANTITRYPSIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>	36 months	serum or plasma	0.007 to 3.5 (g/l)
	ATTUR-B00/ATE	1 x 1 x 13,2 ml R1 + 1 x 1 x 4,6 ml R2	99			
	ATTUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 7,5 ml R2	385			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CATALOGUE 2025

ALPHA-1 MICROGLOBULIN URINARY

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>	17 months	urine	1.29 to 166.6 (mg/l)
	A1DGX-B00/ATE	1 x 1 x 11 ml R1 + 1 x 1 x 11 ml R1	100			
	A1DGX-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 24 ml R2	475			

CALIBRATORS

KIT 5 LEV	A1REK-000	5x1 ml	14 months
-----------	-----------	--------	-----------

CONTROLS

M	A1CON-002	1x2 ml	14 months
L	A1COS-002	1x2 ml	14 months

ALPHA-2 MACROGLOBULIN SERUM and URINARY

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>	36 months	serum or plasma or urine	0.026 to 4.96 (g/l) serum 3.36 to 199 (mg/l) urine
	MGTUR-B00/ATE	1 x 1 x 13,1 ml R1 + 1 x 1 x 4,6 ml R2	99			
	MGTUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 7,5 ml R2	384			
	MGTUR-L00/ATE	5 x 2 x 24 ml R1 + 5 x 2 x 7,5 ml R2	1920			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months
URINARY LOW	MGCU1-002	1x2 ml	24 months
URINARY HIGH	MGCU2-002	1x2 ml	24 months

APOLIPOPROTEIN A1

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>	36 months	serum or plasma	0.029 to 4.0 (g/l)
	AATUR-B00/ATE	1 x 1 x 13,2 ml R1 + 1 x 1 x 5,2 ml R1	99			
	AATUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 8,7 ml R2	383			

CALIBRATORS

KIT 5 LEV	ABREK-000	5x0,5 ml	18 months
-----------	-----------	----------	-----------

CONTROLS

H	ABCOX-002	1x2 ml	18 months
M	ABCON-002	1x2 ml	18 months

CATALOGUE 2025

APOLIPOPROTEIN B

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>	36 months	serum or plasma	0.019 to 4.0 (g/l)
	ABTUR-B00/ATE	1 x 1 x 13,2 ml R1 + 1 x 1 x 4,6 ml R1	99			
	ABTUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 7,5 ml R2	385			

CALIBRATORS

KIT 5 LEV	ABREK-000	5 x 0,5 ml	18 months
-----------	-----------	------------	-----------

CONTROLS

H	ABCOX-002	1x2 ml	18 months
M	ABCON-002	1x2 ml	18 months

BETA2-MICROGLOBULIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>	30 months	serum or plasma	0.132 to 95 (mg/l)
	B2DGX-B00/ATE	1 x 1 x 13 ml R1 + 1 x 1 x 6,7 ml R2	99			
	B2DGX-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 12,7 ml R2	427			
	B2DGX-L00/ATE	5 x 2 x 24 ml R1 + 5 x 2 x 12,7 ml R2	2137			

CALIBRATORS

KIT 5 LEV	B2REK-000	5x1 ml	30 months
-----------	-----------	--------	-----------

CONTROLS

H	B2COX-001	1x1 ml	30 months
M	B2CON-003	1x3 ml	30 months
L	B2COS-003	1x3 ml	30 months

BETA2-MICROGLOBULIN URINARY

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT		<i>R1</i>	<i>R2</i>	30 months	Urine	0.055 to 9.3 (mg/l)
	BUDGX-B00/ATE	1 x 1 x 17 ml R1 + 1 x 1 x 8,9 ml R2	165			
	BUDGX-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 12,3 ml R2	480			

CALIBRATORS

KIT 5 LEV	BUREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

M	B2CON-003	1x3 ml	30 months
L	B2COS-003	1x3 ml	30 months

CATALOGUE 2025

CALPROTECTIN FECAL - Colloidal Gold

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i>	<i>R2</i>			
KIT	CACOL-B00/ATE	1 x 1 x 19 ml R1 + 1 x 1 x 10 ml R2	200	18 months	Faeces	10.93 to 1297.4 (mg/kg)
	CACOL-H00/ATE	1 x 2 x 19 ml R1 + 1 x 2 x 10 ml R2	400			
	CACOL-L00/ATE	5 x 2 x 19 ml R1 + 5 x 2 x 10 ml R2	2000			

CALIBRATORS

KIT 6 LEV	CAREK-000	6x1 ml	12 months
-----------	-----------	--------	-----------

CONTROLS

H	CACOX-002	1x2 ml	12 months
M	CACON-002	1x2 ml	12 months
L	CACOS-002	1x2 ml	12 months

CALPRO EASY EXTRACT™

KIT	CAL0510	50 devices
-----	---------	------------

CERULOPLASMIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i>	<i>R2</i>			
KIT	CETUR-B00/ATE	1 x 1 x 13,2 ml R1 + 1 x 1 x 4,6 ml R2	99	36 months	serum or plasma	9.48 to 1051.20 (mg/l)
	CETUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 7,5 ml R2	385			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	2-8 °C	24 months
M	MPCON-002	1x2 ml	2-8 °C	24 months
L	MPCOS-002	1x2 ml	2-8 °C	24 months

COMPLEMENT C3

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i>	<i>R2</i>			
KIT	C3TUR-B00/ATE	1 x 1 x 13,2 ml R1 + 1 x 1 x 4,6 ml R2	99	36 months	serum or plasma	0.015 to 4.31 (g/l)
	C3TUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 7,5 ml R2	385			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CATALOGUE 2025

COMPLEMENT C4

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i>	<i>R2</i>			
KIT	C4TUR-B00/ATE	1 x 1 x 13,2 ml R1 + 1 x 1 x 4,6 ml R2	99	36 months	serum or plasma	0.002 to 0.825 (g/l)
	C4TUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 7,5 ml R2	385			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CRP GOLD LATEX

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i>	<i>R2</i>			
KIT	CPDGX-B00/ATE	1 x 1 x 7 ml R1 + 1 x 1 x 7 ml R2	105	24 months	serum or plasma	0.031 to 300 (mg/l)
	CPDGX-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 24 ml R2	859			
	CPDGX-L00/ATE	5 x 2 x 24 ml R1 + 5 x 2 x 24 ml R2	4294			

CALIBRATORS

KIT 5 LEV	CPRLK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	24 months
M	CPCON-002	1x2 ml	24 months
L	CXCON-002	1x2 ml	24 months

CYSTATIN C

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i>	<i>R2</i>			
KIT	CYCOL-B00/ATE	1 x 1 x 10,1 ml R1 + 1 x 1 x 5,7 ml R2	100	24 months	serum or plasma	0.15 to 9.36 (mg/l)
	CYCOL-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 12,3 ml R2	517			

CALIBRATORS

KIT 5 LEV	CYREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	CYCOX-002	1x2 ml	24 months
M	CYCON-002	1x2 ml	24 months
L	CYCOX-002	1x2 ml	24 months

CATALOGUE 2025

HAPTOGLOBIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i>	<i>R2</i>			
KIT	HATUR-B00/ATE	1 x 1 x 13,1 ml R1 + 1 x 1 x 4,6 ml R2	99	36 months	serum or plasma	0.002 to 3.86 (g/l)
	HATUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 7,5 ml R2	385			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS A

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i>	<i>R2</i>			
KIT	IATUR-B00/ATE	1 x 1 x 13,1 ml R1 + 1 x 1 x 4,6 ml R2	99	36 months	serum	0.02 to 8.0 (g/l)
	IATUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 7,5 ml R2	385			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS G

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
		<i>R1</i>	<i>R2</i>			
KIT	IGTUR-B00/ATE	1 x 1 x 13,1 ml R1 + 1 x 1 x 6,7 ml R2	100	36 months	serum	0.125 to 35 (g/l)
	IGTUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 7,6 ml R2	388			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULINS M

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
		<i>R1</i>	<i>R2</i>			
KIT	IMTUR-B00/ATE	1 x 1 x 24 ml R1 + 1 x 1 x 7,5 ml R2	192	36 months	serum	0.038 to 7.250 (g/l)
	IMTUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 7,5 ml R2	384			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CATALOGUE 2025

LIPOPROTEIN A

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
		<i>R1</i>	<i>R2</i>			
KIT	LPTUR2-B00/ATE	1 x 1 x 13,1 ml R1 + 1 x 1 x 3,8 ml R2	97	12 months	serum or plasma	0.0235 to 0.946 (g/l)
	LPTUR2-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 6,1 ml R2	387			

CALIBRATORS

KIT 5 LEV	LPREK-000	5x1 ml		18 months
-----------	-----------	--------	--	-----------

CONTROLS

M	LPCON-002	1x2 ml		24 months
---	-----------	--------	--	-----------

MICROALBUMIN

REAGENTS	Reference	Content	Tests*	Shelf Life*	Sample**	Analytical Range
		<i>R1</i>	<i>R2</i>			
KIT	MATUR-B00/ATE	1 x 1 x 13,1 ml R1 + 1 x 1 x 4,6 ml R2	99	36 months	urine	1.86 to 400 (mg/l)
	MATUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 7,5 ml R2	385			

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml		24 months
-----------	-----------	--------	--	-----------

CONTROLS

H	MACOX-002	1x2 ml		24 months
M	MACON-002	1x2 ml		24 months
L	MACOS-002	1x2 ml		24 months

PREALBUMIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i>	<i>R2</i>			
KIT	PATUR-B00/ATE	1 x 1 x 13,1 ml R1 + 1 x 1 x 4,7 ml R2	100	36 months	serum or plasma	6.53 to 1000 (mg/l)
	PATUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 7 ml R2	263			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml		24 months
-----------	-----------	--------	--	-----------

CONTROLS

H	MPCOX-002	1x2 ml		24 months
M	MPCON-002	1x2 ml		24 months
L	MPCOS-002	1x2 ml		24 months

CATALOGUE 2025

RETINOL BINDING PROTEIN (RBP)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i>	<i>R2</i>			
KIT	RBTUR-B00/ATE	1 x 1 x 11 ml R1 + 1 x 1 x 5,3 ml R2	98	36 months	serum or plasma	7,81 to 141 (mg/l)
	RBTUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 10,6 ml R2	470			

CALIBRATORS

KIT 5 LEV	RBREK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

H	RBCOX-002	1x2 ml	18 months
M	RBCON-002	1x2 ml	18 months
L	RBCOS-002	1x2 ml	18 months

RHEUMATOID FACTOR

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i>	<i>R2</i>			
KIT	RFCOL-B00/ATE	1 x 1 x 10,5 ml R1 + 1 x 1 x 7,8 ml R2	100	18 months	serum or plasma	2.228 to 120 (kIU/l)
	RFCOL-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 17,3 ml R2	502			

CALIBRATORS

KIT 5 LEV	RFREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	RFCOX-002	1x2 ml	24 months
M	RFCON-002	1x2 ml	24 months
L	RFCOS-002	1x2 ml	24 months

TRANSFERRIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
		<i>R1</i>	<i>R2</i>			
KIT	TFTUR-B00/ATE	1 x 1 x 13,5 ml R1 + 1 x 1 x 5,7 ml R2	103	36 months	serum or plasma	0.027 to 8.19 (g/l)
	TFTUR-H00/ATE	1 x 2 x 24 ml R1 + 1 x 2 x 9,2 ml R2	388			
	TFTUR-L00/ATE	5 x 2 x 24 ml R1 + 5 x 2 x 9,2 ml R2	1940			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CATALOGUE 2025

Please contact us for more information on the following tests:

ANTISTREPTOLYSIN O						
REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
CRP XL Wide Range						
REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
FERRITIN						
REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
ALPHA-1-GLYCOPROTEIN ACID - OROSOMUCOID						
REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
HDL CHOLESTEROL						
REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
IMMUNOGLOBULIN E						
REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****

CATALOGUE 2025

Your expert in specific proteins



(*) Tests : Amount of real tests. May vary from a analyzer to another one.

() Shelf Life :** Diagam provides to distributors a minimum of *12 months*. The time period written in this document shows the maximum shelf life time after manufacturing.

(*) Sample :** Diagam may provide some applications for other samples, such as CSF, Urine, Veterinary samples etc.... Please contact your distributor or us directly at mail@diagam.com or dial +32 68 55 14 82.

(**) Analytical range :** For general information purposes only. More details can be found in the package insert (instructions for use). Please contact your distributor or us directly at mail@diagam.com or dial +32 68 55 14 82 in order to get the updated data.



Medical device for in vitro diagnosis.



DIAGAM SA, *Rue du Parc Industriel 40, 7822 Ghislenghien, BELGIUM*



The products fulfill the requirements of European directive 98/79 EC concerning the medical devices of in vitro diagnosis.

© Diagam. 2025. All rights reserved. All product names, registered trademarks, company names in this document remain the property of their respective owners. For general information purposes only. More details can be found in the package insert (instructions for use) and labelling. Not for the US market except for research use only.