



Bulk OEM

CATALOGUE 2025

ALBUMIN

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	serum & plasma	1.096 to 120.44 (g/l)	150/40 (µl)
Ab	ALTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

ALBUMIN (CSF)

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	CSF	10.86 to 494.33 (mg/l)	250/50 (µl)
Ab	MATUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MAREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MACOZ-999	Bulk	2-8 °C	24 months
H	MACOX-999	Bulk	2-8 °C	24 months
M	MACON-999	Bulk	2-8 °C	24 months
L	MACOS-999	Bulk	2-8 °C	24 months

ALPHA-1-ANTITRYPSIN

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum & plasma	0.007 to 3.5 (g/l)	250/50 (µl)
Ab	ATTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

CATALOGUE 2025

ALPHA1-GLYCOPROTEIN ACID / OROSOMUCOID

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	GYBUF-999	Bulk	2-25 °C	36 months	serum	0.143 to 2.42 (g/l)	250/50 (µl)
Ab	GPTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

ALPHA-1 MICROGLOBULIN (URINE)

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	A1BUX-999	Bulk	2-8 °C	17 months	urine	1.29 to 166.6 (mg/l)	100/100 (µl)
Ab	A1DGX-H99	Bulk	2-8 °C	17 months			

CALIBRATORS

KIT 5 LEV	A1REK-999	Bulk	2-8 °C	14 months
-----------	-----------	------	--------	-----------

CONTROLS

M	A1CON-999	Bulk	2-8 °C	14 months
L	A1COS-999	Bulk	2-8 °C	14 months

ALPHA2-MACROGLOBULIN - ALPHA2-MACROGLOBULIN URINARY

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum & plasma	0.026 to 4.96 (g/l)	250/50 (µl)
Ab	MGTUR-H99	Bulk	2-8 °C	36 months	urine	3.36 to 199 (mg/l)	

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months
URINARY LOW	MGCU1-999	Bulk	2-8 °C	24 months
URINARY HIGH	MGCU2-999	Bulk	2-8 °C	24 months

CATALOGUE 2025

ANTI STREPTOLYSIN O (ASO)

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	ASBUF-999	Bulk	2-8 °C	20 months	serum & plasma	19.992 to 864.85 (IU/ml)	235/30 (µl)
Ab	ASLTX-H99	Bulk	2-8 °C	20 months			

CALIBRATORS

KIT 5 LEV	ASREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

M	ASCON-999	Bulk	2-8 °C	18 months
---	-----------	------	--------	-----------

APOLIPOPROTEIN A1

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum & plasma	0.029 to 4 (g/l)	250/66 (µl)
Ab	AATUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	ABREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

H	ABCOX-999	Bulk	2-8 °C	18 months
M	ABCON-999	Bulk	2-8 °C	18 months

APOLIPOPROTEIN B

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum & plasma	0.019 to 4.0 (g/l)	250/50 (µl)
Ab	ABTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	ABREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

H	ABCOX-999	Bulk	2-8 °C	18 months
M	ABCON-999	Bulk	2-8 °C	18 months

CATALOGUE 2025

BETA2-MICROGLOBULIN

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	B2BUX-999	Bulk	2-8 °C	30 months	serum & plasma	0,132 to 95 (mg/l)	180/80 (µl)
Ab	B2DGX- H99	Bulk	2-8 °C	30 months			

CALIBRATORS

KIT 5 LEV	B2REK-999	Bulk	2-8 °C	30 months
-----------	-----------	------	--------	-----------

CONTROLS

H	B2COX-001	lyophilized	2-8 °C	30 months
M	B2CON-003	lyophilized	2-8 °C	30 months
L	B2COS-003	lyophilized	2-8 °C	30 months

BETA2-MICROGLOBULIN (URINE)

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	BUBUX-999	Bulk	2-8 °C	30 months	urine	0.055 to 9.3 (mg/l)	135/60 (µl)
Ab	BUDGX- H99	Bulk	2-8 °C	30 months			

CALIBRATORS

KIT 5 LEV	BUREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

M	B2CON-003	lyophilized	2-8 °C	30 months
L	B2COS-003	lyophilized	2-8 °C	30 months

CATALOGUE 2025

CALPROTECTIN FECAL - Colloidal Gold

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	faeces	10.93 to 1297.4 (mg/kg)	180/75 (µl)
Ab	CACOL- H99	Bulk	2-8 °C	18 months			

CALIBRATORS

KIT 6 LEV	CAREK-999	Bulk	2-8 °C	12 months
-----------	-----------	------	--------	-----------

CONTROLS

H	CACOX-999	Bulk	2-8 °C	12 months
M	CACON-999	Bulk	2-8 °C	12 months
L	CACOS-999	Bulk	2-8 °C	12 months

SAMPLE DILUTION BUFFER

KIT	SDBUF-999	Bulk	2-8 °C	15 months
-----	-----------	------	--------	-----------

CALPRO EASY EXTRACT™

KIT	CAL0510	50 devices
-----	---------	------------

CERULOPLASMIN

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	MEBUF-999	Bulk	2-8 °C	36 months	serum & plasma	9.48 to 1051.2 (mg/l)	250/50 (µl)
Ab	CETUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

CATALOGUE 2025

COMPLEMENT C3

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	serum & plasma	0.015 to 4.31 (g/l)	250/50 (µl)
Ab	C3TUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

COMPLEMENT C4

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	serum & plasma	0.002 to 0.825 (g/l)	250/50 (µl)
Ab	C4TUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

CRP GOLD LATEX

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	CPBUF-999	Bulk	2-8 °C	24 months	serum & plasma	0.031 to 300 (mg/l)	100/100 (µl)
Ab	CPDGX-H99	Bulk	2-8 °C	24 months			

CALIBRATORS

KIT 5 LEV	CPRLK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	CPCOX-999	Bulk	2-8 °C	24 months
M	CPCON-999	Bulk	2-8 °C	24 months
L	CXCON-999	Bulk	2-8 °C	24 months

CATALOGUE 2025

CRP XL WIDE RANGE

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	PITRX-999	Bulk	2-25 °C	36 months	serum & plasma	0.607 to 425 (mg/l)	250/36 (µl)
Ab	CPTXL-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	CPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	CPCOX-999	Bulk	2-8 °C	24 months
M	CPCON-999	Bulk	2-8 °C	24 months
L	CXCON-999	Bulk	2-8 °C	24 months

CYSTATIN C - Colloidal Gold

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum & plasma	0.15 to 9.355 (mg/l)	180/75 (µl)
Ab	CYCOL-H99	Bulk	2-8 °C	24 months			

CALIBRATORS

KIT 5 LEV	CYREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	CYCOX-999	Bulk	2-8 °C	24 months
M	CYCON-999	Bulk	2-8 °C	24 months
L	CYCOS-999	Bulk	2-8 °C	24 months

FERRITIN - Colloidal Gold

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	FDBUF-999	Bulk	2-8 °C	18 months	serum	6.344 to 500 (µg/l)	250/50 (µl)
Ab	FNCOL-H99	Bulk	2-8 °C	12 months			

CALIBRATORS

KIT 5 LEV	FNRGK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

M	FNCON-999	Bulk	2-8 °C	18 months
---	-----------	------	--------	-----------

CATALOGUE 2025

FIBRINOGEN

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	plasma	0.05 to 10 (g/l)	250/50 (µl)
Ab	FITUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

M	FIREM-001	lyophilized	2-8 °C	28 months
---	-----------	-------------	--------	-----------

CONTROLS

M	FICON-001	lyophilized	2-8 °C	28 months
---	-----------	-------------	--------	-----------

HAPTOGLOBIN

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	I7TRX-999	Bulk	2-25 °C	36 months	serum & plasma	0.002 to 3.86 (g/l)	250/50 (µl)
Ab	HATUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

IMMUNOGLOBULINS A

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TIBUF-999	Bulk	2-25 °C	36 months	serum	0.02 to 8.0 (g/l)	250/50 (µl)
Ab	IATUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

CATALOGUE 2025

IMMUNOGLOBULINS E

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	IEBUX-999	Bulk	2-8 °C	20 months	serum & plasma	21.0 to 1436 (UI/ml)	140/70 (µl)
Ab	IEDGX-H99	Bulk	2-8 °C	20 months			

CALIBRATORS

KIT 5 LEV	IEREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

M	IECON-003	lyophilized	2-8 °C	18 months
L	IECOS-003	lyophilized	2-8 °C	18 months

IMMUNOGLOBULINS G

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TIBUF-999	Bulk	2-25 °C	36 months	serum	0.125 to 35 (g/l)	250/95 (µl)
Ab	IGTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

IMMUNOGLOBULINS G (CSF)

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TIBUF-999	Bulk	2-25 °C	36 months	CSF	3.375 to 249.78 (mg/l)	135/25 (µl)
Ab	IGLCR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	IGREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	IGCOX-999	Bulk	2-8 °C	24 months
M	IGCON-999	Bulk	2-8 °C	24 months
L	IGCOS-999	Bulk	2-8 °C	24 months

CATALOGUE 2025

IMMUNOGLOBULINS M

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TIBUF-999	Bulk	2-25 °C	36 months	serum	0.038 to 7.0 (g/l)	250/50 (µl)
Ab	IMTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

LIPOPROTEIN (A)

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum & plasma	0.0235 to 0.946 (g/l)	200/25 (µl)
Ab	LPTUR2-H99	Bulk	2-8 °C	12 months			

CALIBRATORS

KIT 4 LEV	LPREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

M	LPCON-999	Bulk	2-8 °C	24 months
---	-----------	------	--------	-----------

MICROALBUMIN

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	urine	1.86 to 400 (mg/l)	250/50 (µl)
Ab	MATUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MAREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MACOX-999	Bulk	2-8 °C	24 months
M	MACON-999	Bulk	2-8 °C	24 months
L	MACOS-999	Bulk	2-8 °C	24 months

CATALOGUE 2025

PREALBUMIN

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	serum & plasma	6.53 to 1000 (mg/l)	250/50 (µl)
Ab	PATUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

RETINOL BINDING PROTEIN

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	PITRX-999	Bulk	2-25 °C	36 months	serum & plasma	7.81 to 141 (mg/l)	75/25 (µl)
Ab	RBTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	RBREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

H	RBCOX-999	Bulk	2-8 °C	18 months
M	RBCON-999	Bulk	2-8 °C	18 months
L	RBCOS-999	Bulk	2-8 °C	18 months

RHEUMATOID FACTOR (RF) - Colloidal Gold

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum & plasma	2.228 to 120 (kUI/l)	150/100 (µl)
Ab	RFCOL-H99	Bulk	2-8 °C	18 months			

CALIBRATORS

KIT 5 LEV	RFREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	RFCOX-999	Bulk	2-8 °C	24 months
M	RFCON-999	Bulk	2-8 °C	24 months
L	RFCOS-999	Bulk	2-8 °C	24 months

CATALOGUE 2025

TRANSFERRIN

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	serum & plasma	0.076 to 8.0 (g/l)	250/70 (µl)
Ab	TFTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

Please contact us for more information on the following tests:

HDL CHOLESTEROL

REAGENTS	Reference	Content	T°C	Shelf Life*	Sample**	Analytical Range***	BUF / Ab Ratio
----------	-----------	---------	-----	-------------	----------	---------------------	----------------

Your expert in specific proteins



(*) Shelf Life : Diagam provides to distributors a minimum of *12 months* . The time period written in this document shows the maximum shelf life time after manufacturing.

() Sample :** Diagam may provide some applications for other samples, such as CSF, Urin, Veterany samples etc.... Please contact your distributor or us directly at mail@diagam.com or dial +32 68 55 14 82.

(*) Analytical range :** For general information purposes only. More details can be found in the package insert (instructions for use). Please contact your distributor or us directly at mail@diagam.com or dial +32 68 55 14 82 in order to get the updated data.



Medical device for in vitro diagnosis.



DIAGAM SA, *Rue du Parc Industriel, 7822 Ghislenghien BELGIUM*



The products fulfill the requirements of European directive 98/79 EC concerning the medical devices of in vitro diagnosis.

© Diagam. 2025. All rights reserved. All product names, registered trademarks, company names in this document remain the property of their respective owners. For general information purposes only. More details can be found in the package insert (instructions for use) and labelling. Not for the US market except for research use only.