



Your expert in specific proteins

CATALOGUE

Reagents, Calibrators & Controls

Immunoturbidimetry

Vitros® 5,1 FS
Vitros® 4600
Vitros® 5600




ALBUMIN
ALPHA 1 - ACID GLYCOPROTEIN
ALPHA 1 - ANTITRYPSIN
ALPHA 1 - MICROGLOBULIN
ALPHA 2 - MACROGLOBULIN
ANTITHROMBIN III
APOLIPOPROTEIN A1
APOLIPOPROTEIN B
ASO

BETA 2 - MICROGLOBULIN
CALPROTECTIN **(new)**
CERULOPLASMIN
C3
C4
CRP
CYSTATIN C
FERRITIN
FIBRINOGEN
HAPTOGLOBIN


Ig A
Ig E **(new)**
Ig G
Ig M
LIPOPROTEIN A
MICROALBUMIN
PREALBUMIN
RETINOL BINDING PROTEIN
RHEUMATOID FACTOR
TRANSFERRIN





New

Calprotectin


REAGENTS *Reference* *Content* *Shelf Life**  *Sample*** *Analytical Range* *Tests*

Please contact us for more information

CALPROTECTIN TEST

New

SARS-CoV-2 SEROLOGICAL ASSAY

REAGENTS *Reference* *Content* *Shelf Life**  *Sample*** *Analytical Range* *Tests*

Please contact us for more information

SARS-CoV-2 SEROLOGICAL ASSAY based on PETIA Colloidal Gold

ALBUMIN							
REAGENTS	Reference	Content		Shelf Life*	Sample**	Analytical Range	Tests
		R1	R2				
MPA KIT	ALTUR-B00/UDA	2x16 ml	2x3 ml	36 months	Serum and plasma	0.187 to 119 (g/l)	150
	ALTUR-H00/UDA	8x16 ml	8x3 ml				600

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALBUMIN CSF

REAGENTS	Reference	Content		Shelf Life*	Sample**	Analytical Range	Tests
		R1	R2				
UDA KIT	MATUR-B00/UDA	2x14 ml	2x2,8 ml	36 months	CSF	4.6 to 400 (mg/l)	150
	MATUR-H00/UDA	8x14 ml	8x2,8 ml				600
	MATUR-L00/UDA	20x14 ml	20x2,8 ml				1500

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

H	MACOZ-002	1x2 ml	18 months
H	MACOX-002	1x2 ml	18 months
M	MACON-002	1x2 ml	18 months
L	MACOS-002	1x2 ml	18 months

ALPHA1-GLYCOPROTEIN ACID/OROSOMUCOID

REAGENTS	Reference	Content		Shelf Life*	Sample**	Analytical Range	Tests
		R1	R2				
MPA KIT	GPTUR-B00/UDA	2x18 ml	2x3,6 ml	36 months	Serum	0.043 to 2.0 (g/l)	150
	GPTUR-H00/UDA	8x18 ml	8x3,6 ml				600
	GPTUR-L00/UDA	20x18 ml	20x3,6 ml				1500

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA1-ANTITRYPSIN

REAGENTS	Reference	Content		Shelf Life*	Sample**	Analytical Range	Tests
		R1	R2				
UDA KIT	ATTUR-B00/UDA	2x18 ml	2x4,4 ml	24 months	Serum and plasma	0.007 to 3.5 (g/l)	150
	ATTUR-H00/UDA	8x18 ml	8x4,4 ml				600

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA1-MICROGLOBULIN (URINE)

Contact us

ALPHA2-MACROGLOBULIN - ALPHA2-MACROGLOBULIN URINARY

REAGENTS	Reference	Content		Shelf Life*	Sample**	Analytical Range	Tests
		R1	R2				
MPA KIT	MGTUR-B00/UDA	2x17 ml	2x3,4 ml	36 months	Serum and plasma or urine	0.005 to 5.43 (g/l)	150
	MGTUR-H00/UDA	8x17 ml	8x3,4 ml				600

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months
URINARY LOW	MGCU1-002	1x2 ml	12 months
URINARY HIGH	MGCU2-002	1x2 ml	12 months

ANTITHROMBIN III

REAGENTS	Reference	Content		Shelf Life*	Sample**	Analytical Range	Tests
		R1	R2				
MPA KIT	A3TUR-B00/UDA	2x18 ml	2x3,6 ml	36 months	Plasma	42.8 to 900 (mg/l)	150
	A3TUR-H00/UDA	8x18 ml	8x3,6 ml				600

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

For VITROS® 5,1 FS Chemistry, VITROS® 4600 Chemistry and VITROS® 5600 Integrated Systems

BETA2-MICROGLOBULIN

REAGENTS	Reference	Content		Shelf Life*	Sample**	Analytical Range	Tests
		R1	R2				
MPA KIT	B2DGX-B00/UDA	2x16,5 ml	2x8 ml	15 months	Serum and plasma	0.48 to 82 (mg/l)	160
	B2DGX-H00/UDA	8x16,5 ml	8x8 ml				640
	B2DGX-L00/UDA	20x16,5 ml	20x8 ml				1600

CALIBRATORS

KIT 5 LEV	B2REK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

H	B2COX-001	1 x 1 ml	18 months
M	B2CON-003	1x3 ml	18 months
L	B2COS-003	1x3 ml	18 months

BETA2-MICROGLOBULIN (URINE)

REAGENTS	Reference	Content		Shelf Life*	Sample**	Analytical Range	Tests
		R1	R2				
MPA KIT	BUDGX-B00/UDA	Contact us			Urine	0.055 to 9.3 (mg/l)	160
	BUDGX-H00/UDA						640
	BUDGX-L00/UDA						1600

CALIBRATORS

KIT 5 LEV	BUREK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

M	B2CON-003	1x3 ml	18 months
L	B2COS-003	1x3 ml	18 months

CERULOPLASMIN

REAGENTS	Reference	Content		Shelf Life*	Sample**	Analytical Range	Tests
		R1	R2				
MPA KIT	CETUR-B00/UDA	2x14 ml	2x2,8 ml	36 months	Serum / plasma	3.1 to 910 (mg/l)	150
	CETUR-H00/UDA	8x14 ml	8x2,8 ml				600

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

For VITROS® 5,1 FS Chemistry, VITROS® 4600 Chemistry and VITROS® 5600 Integrated Systems

CRP GOLD LATEX

REAGENTS	Reference	Content		Shelf Life*	Sample**	Analytical Range	Tests
		R1	R2				
UDA KIT	CPDGX-B00/UDA	2x11 ml	2x9 ml	24 months	Serum and plasma	0.24 to 330 (mg/l)	160
	CPDGX-H00/UDA	8x11 ml	8x9 ml				640
	CPDGX-L00/UDA	20x11 ml	20x9 ml				1600

CALIBRATORS

KIT 5 LEV	CPRLK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	18 months
M	CPCON-002	1x2 ml	18 months
L	CXCEN-002	1x2 ml	18 months

CRP XL WIDE RANGE

REAGENTS	Reference	Content		Shelf Life*	Sample**	Analytical Range	Tests
		R1	R2				
UDA KIT	CPTXL-B00/UDA	2x17 ml	2x2,5 ml	36 months	Serum and plasma	1.09 to 300 (mg/l)	150
	CPTXL-H00/UDA	8x17 ml	8x2,5 ml				600
	CPTXL-L00/UDA	20x17 ml	20x2,5 ml				1500

CALIBRATORS

KIT 5 LEV	CPREK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	18 months
M	CPCON-002	1x2 ml	18 months
L	CXCEN-002	1x2 ml	18 months

CYSTATIN C - Colloidal Gold

REAGENTS	Reference	Content		Shelf Life*	Sample**	Analytical Range	Tests
		R1	R2				
UDA KIT	CYCOL-B00/UDA	2x15 ml	2x6 ml	6 months	Serum and plasma	0.22 to 7.58 (mg/l)	150
	CYCOL-H00/UDA	8x15 ml	8x6 ml				600

CALIBRATORS

KIT 5 LEV	CYREK-000	5x1 ml	12 months
-----------	-----------	--------	-----------

CONTROLS

H	CYCOX-002	1x2 ml	12 months
M	CYCON-002	1x2 ml	12 months
L	CYCOS-002	1x2 ml	12 months

For VITROS® 5,1 FS Chemistry, VITROS® 4600 Chemistry and VITROS® 5600 Integrated Systems

FIBRINOGEN

REAGENTS	Reference	Content		Shelf Life*	Sample**	Analytical Range	Tests
		R1	R2				
MPA KIT	FITUR-B00/UDA	2x17 ml	2x4,4 ml	36 months	Plasma	0.072 to 10 (g/l)	150
	FITUR-H00/UDA	8x17 ml	8x4,4 ml				600

CALIBRATORS

M	FIREM-001	1x1 ml	12 months
---	-----------	--------	-----------

CONTROLS

M	FICON-001	1x1 ml	12 months
---	-----------	--------	-----------

HAPTOGLOBIN

REAGENTS	Reference	Content		Shelf Life*	Sample**	Analytical Range	Tests
		R1	R2				
UDA KIT	HATUR-B00/UDA	2x17 ml	2x4,4 ml	36 months	Serum / plasma	0.011 to 4.0 (g/l)	150
	HATUR-H00/UDA	8x17 ml	8x4,4 ml				600

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

IMMUNOGLOBULIN G (CSF)

Contact us

LIPOPROTEIN (A)

REAGENTS	Reference	Content		Shelf Life*	Sample**	Analytical Range	Tests
		R1	R2				
UDA KIT	LPTUR2-B00/UDA	2x15 ml	2x1,8 ml	12 months	Serum and plasma	Contact us	150
	LPTUR2-H00/UDA	8x15 ml	8x1,8 ml				600


CALIBRATORS

KIT 4 LEV	LPREK-000	4x0,5 ml	10 months
-----------	-----------	----------	-----------

CONTROLS

M	LPCON-002	1x2 ml	12 months
---	-----------	--------	-----------

MICROALBUMIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	Tests
		R1	R2				
UDA KIT	MATUR-B00/UDA	2x14 ml	2x2,8 ml	36 months	Urine	4.6 to 400 (mg/l)	150
	MATUR-H00/UDA	8x14 ml	8x2,8 ml				600
	MATUR-L00/UDA	20x14 ml	20x2,8 ml				1500

CALBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	2-8 °C	18 months
-----------	-----------	--------	--------	-----------

CONTROLS

H	MACOX-002	1x2 ml	2-8 °C	18 months
M	MACON-002	1x2 ml		18 months
L	MACOS-002	1x2 ml		18 months


RETINOL BINDING PROTEIN

Contact us



Your expert in specific proteins

(*) Shelf Life :

 DiAgam provides to distributors a minimum of 12 months. The time period written in this document shows the maximum shelf life time after manufacturing.

() Sample :**

DiAgam can provide some applications for other samples, such as CSF, Urin, Veterinary samples etc.... Please contact your distributor or us directly at mail@diagam.com or dial +32 68 55 14 82.

(*) Analytical range :**

For general information purposes only. More details can be found in the package insert (instructions for use). Please contact your distributor or us directly at mail@diagam.com or dial +32 68 55 14 82 in order to get the updated data.



Medical device for in vitro diagnosis.



DIAGAM SA, Rue du Parc Industriel 40, 7822 Ghislenghien, BELGIUM



The products fulfill the requirements of European directive 98/79 EC concerning the medical devices of in vitro diagnosis.

© DiAgam. 2021. All rights reserved. All product names, registered trademarks, company names in this document remain the property of their respective owners. For general information purposes only. More details can be found in the package insert (instructions for use) and labelling. Not for the US market except for research use only.