



Your expert in specific proteins

CATALOGUE

Reagents, Calibrators & Controls

Immunoturbidimetry

BULK OEM




ALBUMIN
ALPHA 1 - ACID GLYCOPROTEIN
ALPHA 1 - ANTITRYPSIN
ALPHA 1 - MICROGLOBULIN
ALPHA 2 - MACROGLOBULIN
ANTITHROMBIN III
APOLIPOPROTEIN A1
APOLIPOPROTEIN B
ASO

BETA 2 - MICROGLOBULIN
CALPROTECTIN (new)
CERULOPLASMIN
C3
C4
CRP
CYSTATIN C
FERRITIN
FIBRINOGEN
HAPTOGLOBIN

Ig A
Ig E (new)
Ig G
Ig M
LIPOPROTEIN A
MICROALBUMIN
PREALBUMIN
RETINOL BINDING PROTEIN
RHEUMATOID FACTOR
TRANSFERRIN





New

Calprotectin

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
----------	-----------	---------	---	---	----------	------------------	----------------

Please contact us for more information

CALPROTECTIN TEST



New

SARS-CoV-2 SEROLOGICAL ASSAY

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range
----------	-----------	---------	--	---	----------	------------------

Please contact us for more information

SARS-CoV-2 SEROLOGICAL ASSAY based on PETIA Colloidal Gold

ALBUMIN							
REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	17BUF-999	Bulk	2-25 °C	36 months	serum / plasma	1.096 to 120.44 (g/l)	150/40 (µl)
Ab	ALTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

ALBUMIN (CSF)

[Contact us](#)

ALPHA1-ANTITRYPSIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum / plasma	0.007 to 3.5 (g/l)	250/50 (µl)
Ab	ATTUR-H99	Bulk	2-8 °C	24 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

ALPHA1-GLYCOPROTEIN ACID / OROSOMUCOID

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	GYBUF-999	Bulk	2-25 °C	36 months	serum	0.143 to 2.42 (g/l)	250/50 (µl)
Ab	GPTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

ALPHA-1 MICROGLOBULIN (URINE)

[Contact us](#)

ALPHA2-MACROGLOBULIN - ALPHA2-MACROGLOBULIN URINARY

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
----------	-----------	---------	---	---	----------	------------------	----------------

BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum / plasma or urine	0.026 to 4.96 (g/l)	250/50 (µl)
Ab	MGTUR-H99	Bulk	2-8 °C	36 months		serum 3.36 to 199 (mg/l) urine	

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months
URINARY LOW	MGCU1-999	Bulk	2-8 °C	12 months
URINARY HIGH	MGCU2-999	Bulk	2-8 °C	12 months

ANTI STREPTOLYSIN O (ASO)

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
----------	-----------	---------	---	---	----------	------------------	----------------

BUF	ASBUF-999	Bulk	2-8 °C	20 months	serum	19.992 to 864.85	235/30 (µl)
Ab	ASLTX-H99	Bulk	2-8 °C	20 months		(IU/ml)	

CALIBRATORS

H	ASREH-999	Bulk	2-8 °C	12 months
---	-----------	------	--------	-----------

CONTROLS

M	ASCON-999	Bulk	2-8 °C	12 months
---	-----------	------	--------	-----------

ANTITHROMBIN III

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
----------	-----------	---------	---	---	----------	------------------	----------------

BUF	TIBUF-999	Bulk	2-25 °C	36 months	plasma	28.15 to 806.3 (mg/l)	250/50 (µl)
Ab	A3TUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS



KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

CATALOGUE - Bulk OEM

APOLIPOPROTEIN A1

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum / plasma	0.029 to 4 (g/l)	250/66 (µl)
Ab	AATUR-H99	Bulk	2-8 °C	36 months			



CALIBRATORS

KIT 5 LEV	ABREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

H	ABCOX-999	Bulk	2-8 °C	18 months
M	ABCON-999	Bulk	2-8 °C	18 months

APOLIPOPROTEIN B

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum / plasma	0.019 to 4.0 (g/l)	250/50 (µl)
Ab	ABTUR-H99	Bulk	2-8 °C	36 months			



CALIBRATORS

KIT 5 LEV	ABREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

H	ABCOX-999	Bulk	2-8 °C	18 months
M	ABCON-999	Bulk	2-8 °C	18 months

BETA2-MICROGLOBULIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	B2BUX-999	Bulk	2-8 °C	15 months	serum / plasma	0,132 to 95 (mg/l)	180/80 (µl)
Ab	B2DGX- H99	Bulk	2-8 °C	15 months			



CALIBRATORS

KIT 5 LEV	B2REK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

H	B2COX-001	lyophilized	2-8 °C	18 months
M	B2CON-003	lyophilized	2-8 °C	18 months
L	B2COS-003	lyophilized	2-8 °C	18 months



BETA2-MICROGLOBULIN (URINE)

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	BUBUX-999				Urine	0,055 - 9,3 mg/L	135/60 (µl)
Ab	BUDGX- H99						

Contact us

CATALOGUE - Bulk OEM

CERULOPLASMIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	MEBUF-999	Bulk	2-8 °C	36 months	serum / plasma	9.48 to 1051.2 (mg/l)	250/50 (µl)
Ab	CETUR-H99	Bulk	2-8 °C	36 months			



CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

COMPLEMENT C3

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	serum / plasma	0.015 to 4.31 (g/l)	250/50 (µl)
Ab	C3TUR-H99	Bulk	2-8 °C	36 months			



CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

COMPLEMENT C4

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	serum / plasma	0.002 to 0.825 (g/l)	250/50 (µl)
Ab	C4TUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS



KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

CATALOGUE - Bulk OEM

CRP GOLD LATEX

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	CPBUF-999	Bulk	2-8 °C	24 months	serum / plasma	0.031 to 300 (mg/l)	100/100 (µl)
Ab	CPDGX-H99	Bulk	2-8 °C	24 months			



CALIBRATORS

KIT 5 LEV	CPRLK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

H	CPCOX-999	Bulk	2-8 °C	18 months
M	CPCON-999	Bulk	2-8 °C	18 months
L	CXCON-999	Bulk	2-8 °C	18 months

CRP XL WIDE RANGE

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	PITRX-999	Bulk	2-25 °C	36 months	serum / plasma	0.607 to 425 (mg/l)	250/36 (µl)
Ab	CPTXL-H99	Bulk	2-8 °C	36 months			



CALIBRATORS

KIT 5 LEV	CPREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

H	CPCOX-999	Bulk	2-8 °C	18 months
M	CPCON-999	Bulk	2-8 °C	18 months
L	CXCON-999	Bulk	2-8 °C	18 months

CRP XS ULTRA - Colloidal Gold

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	PITRX-999	Bulk	2-25 °C	36 months	serum / plasma	0.471 to 285.05 (mg/l)	250/20 (µl)
Ab	CPTXS-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	CPR(EK/LK)-999	Bulk	2-8 °C	18 months
-----------	----------------	------	--------	-----------

CONTROLS

H	CPCOX-999	Bulk	2-8 °C	18 months
M	CPCON-999	Bulk	2-8 °C	18 months
L	CXCON-999	Bulk	2-8 °C	18 months

CATALOGUE - Bulk OEM

CYSTATIN C - Colloidal Gold

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum / plasma	0.15 to 9.355 (mg/l)	180/75 (µl)
Ab	CYCOL-H99	Bulk	2-8 °C	6 months			

CALIBRATORS

KIT 5 LEV	CYREK-999	Bulk	2-8 °C	12 months
-----------	-----------	------	--------	-----------

CONTROLS

H	CYCOX-999	Bulk	2-8 °C	12 months
M	CYCON-999	Bulk	2-8 °C	12 months
L	CYCOS-999	Bulk	2-8 °C	12 months

FERRITIN - Colloidal Gold

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	FDBUF-999	Bulk	2-8 °C	18 months	serum	6.344 to 500 (µg/l)	250/50 (µl)
Ab	FNCOL-H99	Bulk	2-8 °C	12 months			



CALIBRATORS

KIT 5 LEV	FNRGK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

M	FNCON-999	Bulk	2-8 °C	18 months
---	-----------	------	--------	-----------

FIBRINOGEN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	plasma	0.05 to 10 (g/l)	250/80 (µl)
Ab	FITUR-H99	Bulk	2-8 °C	36 months			



CALIBRATORS

M	FIREM-001	lyophilized	2-8 °C	12 months
---	-----------	-------------	--------	-----------

CONTROLS

M	FICON-001	lyophilized	2-8 °C	12 months
---	-----------	-------------	--------	-----------

HAPTOGLOBIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	I7TRX-999	Bulk	2-25 °C	36 months	serum	0.002 to 3.86 (g/l)	250/50 (µl)
Ab	HATUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS



KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

CATALOGUE - Bulk OEM

IMMUNOGLOBULINS A

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	TIBUF-999	Bulk	2-25 °C	36 months	serum	0.02 to 8.0 (g/l)	250/50 (µl)
Ab	IATUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

IMMUNOGLOBULINS E

[Contact us](#)

IMMUNOGLOBULINS G

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	TIBUF-999	Bulk	2-25 °C	36 months	serum	0.125 to 35 (g/l)	250/95 (µl)
Ab	IGTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------



CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

IMMUNOGLOBULINS G (CSF)

[Contact us](#)

IMMUNOGLOBULINS M

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	TIBUF-999	Bulk	2-25 °C	36 months	serum	0.038 to 7.0 (g/l)	250/50 (µl)
Ab	IMTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS



KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

CATALOGUE - Bulk OEM

LIPOPROTEIN (A)

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum / plasma	0.0235 to 0.946 (g/l)	200/25 (µl)
Ab	LPTUR2-H99	Bulk	2-8 °C	12 months			



CALIBRATORS

KIT 5 LEV	LPREK-999	Bulk	2-8 °C	10 months
-----------	-----------	------	--------	-----------

CONTROLS

M	LPCON-999	Bulk	2-8 °C	12 months
---	-----------	------	--------	-----------

MICROALBUMIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	urine	1.86 to 400 (mg/l)	250/50 (µl)
Ab	MATUR-H99	Bulk	2-8 °C	36 months			



CALIBRATORS

KIT 5 LEV	MAREK-999	Bulk	2-8 °C	18 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MACOX-999	Bulk	2-8 °C	18 months
M	MACON-999	Bulk	2-8 °C	18 months
L	MACOS-999	Bulk	2-8 °C	18 months

PREALBUMIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	serum / plasma	6.53 to 1000 (mg/l)	250/50 (µl)
Ab	PATUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS



H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months

RETINOL BINDING PROTEIN

Contact us

CATALOGUE - Bulk OEM

RHEUMATOID FACTOR (RF) - Colloidal Gold

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	TRBUF-999	Bulk	2-8 °C	36 months	serum	2.092 to 120 (kUI/l)	150/100 (µl)
Ab	RFCOL-H99	Bulk	2-8 °C	18 months			



CALIBRATORS

KIT 5 LEV	RFREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

M	RFCON-999	Bulk	2-8 °C	24 months
---	-----------	------	--------	-----------

TRANSFERRIN

REAGENTS	Reference	Content		Shelf Life* 	Sample**	Analytical Range	BUF / Ab Ratio
BUF	I7BUF-999	Bulk	2-25 °C	36 months	serum	0.076 to 8.0 (g/l)	250/70 (µl)
Ab	TFTUR-H99	Bulk	2-8 °C	36 months			

CALIBRATORS

KIT 5 LEV	MPREK-999	Bulk	2-8 °C	24 months
-----------	-----------	------	--------	-----------

CONTROLS

H	MPCOX-999	Bulk	2-8 °C	24 months
M	MPCON-999	Bulk	2-8 °C	24 months
L	MPCOS-999	Bulk	2-8 °C	24 months



Your expert in specific proteins

(*) Shelf Life : DiAgam provides to distributors a minimum of *12 months (except for Cystatin-c, Ferritin and Lp(a): < 12 months)*. The time period written in this document shows the maximum shelf life time after manufacturing.

() Sample :** DiAgam may provide some applications for other samples, such as CSF, Urin, Veterinary samples etc.... Please contact your distributor or us directly at mail@diagam.com or dial *+32 68 55 14 82*.

(*) Analytical range :** For general information purposes only. More details can be found in the package insert (instructions for use). Please contact your distributor or us directly at mail@diagam.com or dial *+32 68 55 14 82* in order to get the updated data.

(**)** Bulk formats or any production require a minimum volume of order of 250 mL for calibrators and controls. If less than 250 mL of calibrators and controls are ordered, DiAgam may refuse to process the order.



Medical device for in vitro diagnosis.



DIAGAM SA, *Rue du Parc Industriel, 7822 Ghislenghien BELGIUM*



The products fulfill the requirements of European directive 98/79 EC concerning the medical devices of in vitro diagnosis.

© DiAgam. 2021. All rights reserved. All product names, registered trademarks, company names in this document remain the property of their respective owners. For general information purposes only. More details can be found in the package insert (instructions for use) and labelling. Not for the US market except for research use only.

