

ADVIA® LINE

**New**

SARS-CoV-2 SEROLOGICAL ASSAY

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
----------	-----------	---------	--------	--------------	-----------	----------------------

*Please contact us for more information*

SARS-CoV-2 SEROLOGICAL ASSAY based on PECIA Colloidal Gold



ALBUMIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	ALTUR-B00/ADV	1 x 30 ml R1 + 1 x 8 ml R2		300	36 months	serum & plasma	0.815 to 124.94 (g/l)
	ALTUR-L00/ADV	2 x 70 ml R1 + 2 x 18 ml R2		1400			
	ALTUR-L00/ADV	5 x 70 ml R1 + 5 x 18 ml R2		3500			

**CALIBRATORS**

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALBUMIN CSF							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	MATUR-B00/ADV	1 x 46 ml R1 + 1 x 8 ml R2		400	36 months	cerebrospinal fluid	10.86 to 494.33 (mg/l)
	MATUR-L00/ADV	2 x 50 ml R1 + 1 x 15 ml R2		750			

**CALIBRATORS**

KIT 5 LEV	MAREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	MACOZ-002	1x2 ml	24 months
H	MACOX-002	1x2 ml	24 months
M	MACON-002	1x2 ml	24 months
L	MACOS-002	1x2 ml	24 months

ALPHA-1 GLYCOPROTEIN ACIDE/OROSOMUCOID							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	GPTUR-B00/ADV	1 x 41 ml R1 + 1 x 8 ml R2		267	36 months	serum	0.143 to 2.42 (g/l)
	GPTUR-H00/ADV	2 x 70 ml R1 + 2 x 12 ml R2		800			
	GPTUR-C00/ADV	2 x 35 ml R1 + 2 x 8 ml R2		467			

**CALIBRATORS**

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

**ALPHA-1 ANTITRYPSIN**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	ATTUR-B00/ADV	1 x 46 ml R1 + 1 x 8 ml R2		400	36 months	serum & plasma	0.007 to 3.5 (g/l)
	ATTUR-H00/ADV	2 x 35 ml R1 + 2 x 9 ml R2		700			
	ATTUR-L00/ADV	1 x 70 ml R1 + 1 x 14 ml R2		700			

**CALIBRATORS**

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

**ALPHA-1 MICROGLOBULIN**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	A1DGX-B00/ADV	1 x 15 ml R1 + 1 x 15 ml R2		375	17 months	urine	1.29 to 166.6 (mg/l)
	A1DGX-H00/ADV	2 x 20 ml R1 + 2 x 20 ml R2		1000			

**CALIBRATORS**

KIT 5 LEV	A1REK-000	5x1 ml	14 months
-----------	-----------	--------	-----------

**CONTROLS**

M	A1CON-002	1x2 ml	14 months
L	A1COS-002	1x2 ml	14 months

**ALPHA-2 MACROGLOBULIN - ALPHA-2 MACROGLOBULIN URINARY**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	MGTUR-H00/ADV	1 x 50 ml R1 + 1 x 10 ml R2		500	36 months	serum & plasma urine	0.026 to 4.96 (g/l) 3.36 to 199 (mg/l)
	MGTUR-C00/ADV	2 x 35 ml R1 + 2 x 9 ml R2		700			
	MGTUR-L00/ADV	12 x 70 ml R1 + 12 x 12 ml R2		7200			

**CALIBRATORS**

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months
URINARY LOW	MGCU1-002	1x2 ml	24 months
URINARY HIGH	MGCU2-002	1x2 ml	24 months

**ANTI STREPTOLYSIN O (ASO)**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	ASLTX-B00/ADV	1 x 62 ml R1 + 1 x 8 ml R2		615	20 months	serum & plasma	19.99 to 864.85 (IU/ml)
	ASLTX-H00/ADV	2 x 70 ml R1 + 2 x 15 ml R2		1400			

**CALIBRATORS**

H	ASREK-000	5x1 ml	18 months
---	-----------	--------	-----------

**CONTROLS**

M	ASCON-002	1x2 ml	18 months
---	-----------	--------	-----------

**APOLIPOPROTEIN A1**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	AATUR-B00/ADV	1 x 46 ml R1 + 1 x 12 ml R2		444	36 months	serum & plasma	0.029 to 4.0 (g/l)
	AATUR-H00/ADV	2 x 70 ml R1 + 2 x 20 ml R2		1400			

**CALIBRATORS**

KIT 5 LEV	ABREK-000	5x0.5 ml	18 months
-----------	-----------	----------	-----------

**CONTROLS**

H	ABCOX-002	1x2 ml	18 months
M	ABCON-002	1x2 ml	18 months

**APOLIPOPROTEIN B**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	ABTUR-B00/ADV	1 x 46 ml R1 + 1 x 8 ml R2		400	36 months	serum & plasma	0.019 to 4.0 (g/l)
	ABTUR-H00/ADV	2 x 70 ml R1 + 2 x 18 ml R2		1400			

**CALIBRATORS**

KIT 5 LEV	ABREK-000	5x0.5 ml	18 months
-----------	-----------	----------	-----------

**CONTROLS**

H	ABCOX-002	1x2 ml	18 months
M	ABCON-002	1x2 ml	18 months

**BETA 2 MICROGLOBULIN**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	B2DGX-B00/ADV	1 x 20 ml R1 + 1 x 8 ml R2		200	30 months	serum & plasma	0.132 to 95 (mg/l)
	B2DGX-H00/ADV	1 x 44 ml R1 + 1 x 19 ml R2		475			
	B2DGX-L00/ADV	12 x 44 ml R1 + 12 x 19 ml R2		5700			

**CALIBRATORS**

KIT 5 LEV	B2REK-000	5x1 ml	30 months
-----------	-----------	--------	-----------

**CONTROLS**

H	B2COX-001	1 x 1 ml	30 months
M	B2CON-003	1x3 ml	30 months
L	B2COS-003	1x3 ml	30 months

**BETA 2 MICROGLOBULIN (Urine)**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	BUDGX-B00/ADV BUDGX-H00/ADV BUDGX-L00/ADV	Contact us			30 months	Urine	0.055 to 9.3 (mg/l)

**CALIBRATORS**

KIT 5 LEV	BUREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

M	B2CON-003	1x3 ml	30 months
L	B2COS-003	1x3 ml	30 months

**CALPROTECTIN FECAL**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	CACOL-B00/ADV CACOL-H00/ADV	Contact us		101 303	18 months	Faecal samples	10,93 to 1297,4 (mg/kg)

**CALIBRATORS**

KIT 6 LEV	CAREK-000	6x1 ml	12 months
-----------	-----------	--------	-----------

**CONTROLS**

H	CACOX-002	1x2 ml	12 months
M	CACON-002	1x2 ml	12 months
L	CACOS-002	1x2 ml	12 months

**SAMPLE DILUTION BUFFER**

KIT	SDBUF-B00	2x70 ml	15 months
	SDBUF-H00	8x70ml	15 months
	SDBUF-L00	16x70 ml	15 months

**CALPRO EASY EXTRACT™**

KIT	CAL0510	50 devices
-----	---------	------------

CERULOPLASMIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	CETUR-B00/ADV	1 x 46 ml R1 + 1 x 8 ml R2		400	36 months	serum & plasma	9.48 to 1051.20 (mg/l)
	CETUR-H00/ADV	2 x 35 ml R1 + 2 x 9 ml R2		700			
	CETUR-L00/ADV	2 x 70 ml R1 + 2 x 18 ml R2		1400			

**CALIBRATORS**

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

COMPLEMENT C3							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	C3TUR-B00/ADV	1 x 46 ml R1 + 1 x 8 ml R2		400	36 months	serum & plasma	0.015 to 4.31 (g/l)
	C3TUR-C00/ADV	2 x 35 ml R1 + 2 x 9 ml R2		700			
	C3TUR-H00/ADV	2 x 70 ml R1 + 2 x 18 ml R2		1400			

**CALIBRATORS**

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

COMPLEMENT C4							
REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	C4TUR-B00/ADV	1 x 46 ml R1 + 1 x 8 ml R2		400	36 months	serum & plasma	0.002 to 0.825 (g/l)
	C4TUR-C00/ADV	2 x 35 ml R1 + 2 x 9 ml R2		700			
	C4TUR-H00/ADV	2 x 70 ml R1 + 2 x 18 ml R2		1400			

**CALIBRATORS**

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

CRP GOLD LATEX							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	CPDGX-B00/ADV	1 x 10 ml R1 + 1 x 10 ml R2		250	24 months	serum & plasma	0.031 to 300 (mg/l)
	CPDGX-H00/ADV	2 x 20 ml R1 + 2 x 20 ml R2		1000			
	CPDGX-L00/ADV	12 x 20 ml R1 + 12 x 20 ml R2		6000			

**CALIBRATORS**

KIT 5 LEV	CPRLK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	CPCOX-002	1x2 ml	24 months
M	CPCON-002	1x2 ml	24 months
L	CXCON-002	1x2 ml	24 months

CRP XL WIDE RANGE							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	CPTXL-B00/ADV	1 x 57 ml R1 + 1 x 8 ml R2		615	36 months	serum & plasma	0.607 to 425 (mg/l)
	CPTXL-H00/ADV	2 x 70 ml R1 + 2 x 10 ml R2		1538			
	CPTXL-L00/ADV	12 x 70 ml R1 + 12 x 15 ml R2		9333			

**CALIBRATORS**

KIT 5 LEV	CPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	CPCOX-002	1x2 ml	24 months
M	CPCON-002	1x2 ml	24 months
L	CXCON-002	1x2 ml	24 months

CYSTATIN C - Colloidal Gold							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	CYCOL-B00/ADV	1 x 25,5 ml R1 + 1 x 9 ml R2		300	24 months	serum & plasma	0.15 to 9.36 (mg/l)
	CYCOL-H00/ADV	1 x 50 ml R1 + 1 x 20,5 ml R2		680			

**CALIBRATORS**

KIT 5 LEV	CYREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	CYCOX-002	1x2 ml	24 months
M	CYCON-002	1x2 ml	24 months
L	CYCOS-002	1x2 ml	24 months

**FERRITIN - Colloidal Gold**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	FNCOL-B00/ADV	1 x 40 ml R1 + 1 x 8 ml R2		533	12 months	serum	6.344 to 500 (µg/l)
	FNCOL-H00/ADV	2 x 70 ml R1 + 2 x 14 ml R2		1867			

**CALIBRATORS**

KIT 5 LEV	FNRGK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

**CONTROLS**

M	FNCON-002	1x2 ml	18 months
---	-----------	--------	-----------

**HAPTOGLOBIN**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	HATUR-B00/ADV	1 x 40 ml R1 + 1 x 8 ml R2		400	36 months	serum & plasma	0.002 to 3.86 (g/l)
	HATUR-C00/ADV	2 x 35 ml R1 + 2 x 9 ml R2		700			
	HATUR-H00/ADV	2 x 70 ml R1 + 2 x 18 ml R2		1400			

**CALIBRATORS**

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

**HDL CHOLESTEROL**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	HDLC-KIT/ADV	12 x 60 ml R1 + 12 x 20 ml R2		Contact us	serum	1 to 180 (mg/dl)	

**CALIBRATORS**

KIT 1 LEV	BIOCAL	5x3 ml	Contact us
-----------	--------	--------	------------

**CONTROLS**

Control +	BIOPATH	6x5 ml	Contact us
Control -	BIONORM	6x5 ml	Contact us



**IMMUNOGLOBULINS A**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	IATUR-B00/ADV	1 x 46 ml R1 + 1 x 8 ml R2		400	36 months	serum	0.02 to 8.0 (g/l)
	IATUR-H00/ADV	2 x 70 ml R1 + 2 x 18 ml R2		1400			

**CALIBRATORS**

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

**IMMUNOGLOBULINS E** **NEW**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	IEDGX-B00/ADV	Contact us		20 months	serum & plasma	21 to 1436 (U/ml)	
	IEDGX-H00/ADV						

**CALIBRATORS**

KIT 5 LEV	IEREK-000	5x1 ml	20 months
-----------	-----------	--------	-----------

**CONTROLS**

M	IECON-002	1x2 ml	20 months
L	IECOS-002	1x2 ml	20 months

**IMMUNOGLOBULINS G**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	IGTUR-B00/ADV	1 x 23 ml R1 + 1 x 8 ml R2		200	36 months	serum	0.125 to 35 (g/l)
	IGTUR-H00/ADV	2 x 50 ml R1 + 2 x 20 ml R2		1000			

**CALIBRATORS**

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

**IMMUNOGLOBULINS G (CSF)**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	IGLCR-B00/ADV	1 x 55 ml R1 + 1 x 10 ml R2		667	24 months	cerebrospinal fluid	3,375 to 249,78 (mg/l)
	IGLCR-H00/ADV	2 x 55 ml R1 + 2 x 10 ml R2		1333			

**CALIBRATORS**

KIT 5 LEV	IGREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	IGCOX-002	1x2 ml	24 months
M	IGCON-002	1x2 ml	24 months
L	IGCOS-002	1x2 ml	24 months

**IMMUNOGLOBULINS M**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	IMTUR-B00/ADV	1 x 46 ml R1 + 1 x 8 ml R2		400	36 months	serum	0.038 to 7.0 (g/l)
	IMTUR-H00/ADV	2 x 70 ml R1 + 2 x 18 ml R2		1400			

**CALIBRATORS**

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

**LIPOPROTEIN (A)**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	LPTUR2-B00/ADV	1 x 70 ml R1 + 1 x 9,1 ml R2		700	12 months	serum & plasma	0.0235 to 0.946 (g/l)
	LPTUR2-H00/ADV	2 x 70 ml R1 + 2 x 9,1 ml R2		1400			

**CALIBRATORS**

KIT 4 LEV	LPREK-000	4x0,5 ml	18 months
-----------	-----------	----------	-----------

**CONTROLS**

M	LPCON-002	1x2 ml	24 months
---	-----------	--------	-----------

MICROALBUMIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	MATUR-B00/ADV	1 x 46 ml R1 + 1 x 8 ml R2		400	36 months	urine	1.86 to 400 (mg/l)
	MATUR-L00/ADV	2 x 50 ml R1 + 1 x 15 ml R2		750			

#### CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

#### CONTROLS

H	MACOX-002	1x2 ml	24 months
M	MACON-002	1x2 ml	24 months
L	MACOS-002	1x2 ml	24 months

PREALBUMIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	PATUR-B00/ADV	1 x 41 ml R1 + 1 x 8 ml R2		267	36 months	serum & plasma	6.53 to 1000 (mg/l)
	PATUR-C00/ADV	2 x 35 ml R1 + 2 x 8 ml R2		467			
	PATUR-H00/ADV	2 x 70 ml R1 + 2 x 15 ml R2		933			

#### CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

#### CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

RETINOL BINDING PROTEIN (RBP)							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	RBTUR-B00/ADV	1 x 31 ml R1 + 1 x 10 ml R2		400	36 months	serum & plasma	7,81 to 141 (mg/l)
	RBTUR-H00/ADV	2 x 31 ml R1 + 2 x 10 ml R2		800			
	RBTUR-L00/ADV	12 x 31 ml R1 + 12 x 10 ml R2		4800			

#### CALIBRATORS

KIT 5 LEV	RBREK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

#### CONTROLS

H	RBCOX-002	1x2 ml	18 months
M	RBCON-002	1x2 ml	18 months
L	RBCOS-002	1x2 ml	18 months

**RHEUMATOID FACTOR (RF) - Colloidal Gold**

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	RFCOL-B00/ADV	1 x 18 ml R1 + 1 x 10 ml R2		200	18 months	serum & plasma	2.228 to 120 (kIU/l)
	RFCOL-C00/ADV	2 x 40 ml R1 + 2 x 20 ml R2		800			
	RFCOL-H00/ADV	2 x 34,5 ml R1 + 2 x 20 ml R2		800			

**CALIBRATORS**

KIT 5 LEV	RFREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	RFCOX-002	1x2 ml	24 months
M	RFCON-002	1x2 ml	24 months
L	RFCOS-002	1x2 ml	24 months

**TRANSFERRIN**

REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	TFTUR-C00/ADV	1 x 30 ml R1 + 1 x 8 ml R2		286	36 months	serum & plasma	0.027 to 8.19 (g/l)
	TFTUR-H00/ADV	2 x 70 ml R1 + 2 x 20 ml R2		1400			

**CALIBRATORS**

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

**CONTROLS**

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

# The Diagam logo features a stylized red and black 'D' icon followed by the word 'Diagam' in a bold, black, sans-serif font.

---

**(\*) Tests :** Amount of theoretical tests. May vary from a analyzer to another one.

**(\*\*) Shelf Life :** Diagam provides to distributors a minimum of *12 months* . The time period written in this document shows the maximum shelf life time after manufacturing.

**(\*\*\*) Sample :** Diagam may provide some applications for other samples, such as LCR, Urin, Veterinary samples etc.... Please contact your distributor or us directly at [mail@diagam.com](mailto:mail@diagam.com) or dial *+32 68 55 14 82*.

**(\*\*\*\*)**  
**Analytical range :** For general information purposes only. More details can be found in the package insert (instructions for use). Please contact your distributor or us directly at [mail@diagam.com](mailto:mail@diagam.com) or dial *+32 68 55 14 82* in order to get the updated data.

 **IVD** Medical device for in vitro diagnosis.



DIAGAM SA, *Rue du Parc Industriel 40, 7822 Ghislenghien, BELGIUM*



The products fulfill the requirements of European directive 98/79 EC concerning the medical devices of in vitro diagnosis.

© Diagam. 2023. All rights reserved. All product names, registered trademarks, company names in this document remain the property of their respective owners. For general information purposes only. More details can be found in the package insert (instructions for use) and labelling. Not for the US market except for research use only.