

ARCHITECT LINE

New

SARS-CoV-2 SEROLOGICAL ASSAY

<i>REAGENTS</i>	<i>Reference</i>	<i>Content</i>	<i>Tests*</i>	<i>Shelf Life**</i>	<i>Sample***</i>	<i>Analytical Range****</i>
-----------------	------------------	----------------	---------------	---------------------	------------------	-----------------------------

Please contact us for more information

SARS-CoV-2 SEROLOGICAL ASSAY based on PECIA Colloidal Gold



CATALOGUE ARCHITECT LINE

ALBUMIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	ALTUR-B00/ARC	R1 1 x 40 ml R1 + 1 x 10 ml R2	169	36 months	serum or plasma	0.40 to 121.71 (g/l)
	ALTUR-H00/ARC	R2 1 x 80 ml R1 + 1 x 20 ml R2	372			
	ALTUR-L00/ARC	7 x 80 ml R1 + 7 x 20 ml R2	2604			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALBUMIN CSF

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	MATUR-B00/ARC	R1 1 x 22,5 ml R1 + 1 x 6,1 ml R2	100	36 months	cerebrospinal fluid	10.86 to 494.33 (mg/l)
	MATUR-H00/ARC	R2 1 x 74 ml R1 + 1 x 17,5 ml R2	361			
	MATUR-L00/ARC	7 x 75 ml R1 + 7 x 15 ml R2	2126			

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MACOZ-002	1 x 2 ml	24 months
H	MACOX-002	1x2 ml	24 months
M	MACON-002	1x2 ml	24 months
L	MACOS-002	1x2 ml	24 months

ALPHA-1 GLYCOPROTEIN ACIDE/OSOMUCOÏD

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	GPTUR-B00/ARC	R1 1 x 37,5 ml R1 + 1 x 7,5 ml R2	97	36 months	serum	0.143 to 2.42 (g/l)
	GPTUR-H00/ARC	R2 2 x 75 ml R1 + 2 x 15 ml R2	443			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA-1 ANTITRYPSIN

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	ATTUR-B00/ARC	R1 1 x 45 ml R1 + 1 x 9 ml R2	122	24 months	serum or plasma	0.007 to 3.5 (g/l)
	ATTUR-H00/ARC	R2 2 x 75 ml R1 + 2 x 15 ml R2	443			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

ALPHA-1 MICROGLOBULIN (URINE)

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	A1DGX-B00/ARC	R1 1 x 15 ml R1 + 1 x 15 ml R2	150	17 months	urine	1,29 to 166,6 (mg/l)
	A1DGX-H00/ARC	R2 1 x 55 ml R1 + 1 x 55 ml R2	456			

CALIBRATORS

KIT 5 LEV	A1REK-000	5x1 ml	14 months
-----------	-----------	--------	-----------

CONTROLS

M	A1CON-002	1x2 ml	14 months
L	A1COS-002	1x2 ml	14 months

ALPHA-2 MACROGLOBULIN - ALPHA-2 MACROGLOBULIN URINARY

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
KIT	MGTUR-H00/ARC	R1 1 x 55 ml R1 + 1 x 11 ml R2	155	36 months	serum/plasma or urine	0.026 to 4.96 (g/l) serum 3.36 to 199 (mg/l) urine
	MGTUR-L00/ARC	R2 10 x 55 ml R1 + 10 x 11 ml R2	1550			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months
URINARY LOW	MGCU1-002	1x2 ml	24 months
URINARY HIGH	MGCU2-002	1x2 ml	24 months

ANTI STREPTOLYSIN O (ASO)							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	ASLTX-B00/ARC	1 x 70 ml R1 + 1 x 8 ml R2		296	20 months	serum or plasma	19.99 to 864.85 (IU/ml)
	ASLTX-H00/ARC	2 x 87 ml R1 + 2 x 10 ml R2		741			
CALIBRATORS							
H	ASREK-000	5x1 ml			18 months		
CONTROLS							
M	ASCON-002	1x2 ml			12 months		

APOLIPOPROTEIN A1							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	AATUR-B00/ARC	1 x 38 ml R1 + 1 x 10 ml R2		107	36 months	serum or plasma	0.029 to 4.0 (g/l)
	AATUR-H00/ARC	2 x 75 ml R1 + 2 x 20 ml R2		474			
CALIBRATORS							
KIT 5 LEV	ABREK-000	5x0.5 ml			18 months		
CONTROLS							
H	ABCOX-002	1x2 ml			18 months		
M	ABCON-002	1x2 ml			18 months		

APOLIPOPROTEIN B							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	ABTUR-B00/ARC	1 x 41 ml R1 + 1 x 8 ml R2		151	36 months	serum or plasma	0.019 to 4.0 (g/l)
	ABTUR-H00/ARC	2 x 67 ml R1 + 2 x 13 ml R2		543			
CALIBRATORS							
KIT 5 LEV	ABREK-000	5x0.5 ml			18 months		
CONTROLS							
H	ABCOX-002	1x2 ml			18 months		
M	ABCON-002	1x2 ml			18 months		

BETA 2 MICROGLOBULIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	B2DGX-B00/ARC	1 x 27 ml R1 + 1 x 12 ml R2		111	30 months	serum or plasma	0,132 to 95 (mg/l)
	B2DGX-H00/ARC	2 x 45 ml R1 + 2 x 20 ml R2		396			
	B2DGX-L00/ARC	10 x 45 ml R1 + 10 x 20 ml R2		1981			

CALIBRATORS

KIT 5 LEV	B2REK-000	5x1 ml	30 months
-----------	-----------	--------	-----------

CONTROLS

H	B2COX-001	1 x 1 ml	30 months
M	B2CON-003	1x3 ml	30 months
L	B2COS-003	1x3 ml	30 months

BETA 2 MICROGLOBULIN (URINE)							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	BUDGX-B00/ARC	1 x 23 ml R1 + 1 x 10 ml R2		117	30 months	Urine	0.055 to 9.3 (mg/l)
	BUDGX-H00/ARC	2 x 23 ml R1 + 2 x 10 ml R2		234			

CALIBRATORS

KIT 5 LEV	BUREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

M	B2CON-003	1x3 ml	30 months
L	B2COS-003	1x3 ml	30 months

CERULOPLASMIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	CETUR-B00/ARC	1 x 40 ml R1 + 1 x 8 ml R2		128	36 months	serum or plasma	9.48 to 1051.20 (mg/l)
	CETUR-H00/ARC	2 x 60 ml R1 + 2 x 12 ml R2		419			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

COMPLEMENT C3							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	C3TUR-B00/ARC	1 x 41 ml R1 + 1 x 8 ml R2		137	36 months	serum or plasma	0.015 to 4.31 (g/l)
	C3TUR-H00/ARC	2 x 56 ml R1 + 2 x 11 ml R2		405			
CALIBRATORS							
KIT 5 LEV	MPREK-000	5x1 ml			24 months		
CONTROLS							
H	MPCOX-002	1x2 ml			24 months		
M	MPCON-002	1x2 ml			24 months		
L	MPCOS-002	1x2 ml			24 months		

COMPLEMENT C4							
REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	C4TUR-B00/ARC	1 x 40 ml R1 + 1 x 8 ml R2		128	36 months	serum or plasma	0.002 to 0.825 (g/l)
	C4TUR-H00/ARC	2 x 60 ml R1 + 2 x 12 ml R2		419			
CALIBRATORS							
KIT 5 LEV	MPREK-000	5x1 ml			24 months		
CONTROLS							
H	MPCOX-002	1x2 ml			24 months		
M	MPCON-002	1x2 ml			24 months		
L	MPCOS-002	1x2 ml			24 months		

CRP GOLD LATEX							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	CPDGX-B00/ARC	1 x 15 ml R1 + 1 x 15 ml R2		105	24 months	serum or plasma	0.031 to 300 (mg/l)
	CPDGX-H00/ARC	1 x 55 ml R1 + 1 x 55 ml R2		456			
	CPDGX-L00/ARC	10 x 55 ml R1 + 10 x 55 ml R2		4561			
CALIBRATORS							
KIT 5 LEV	CPRLK-000	5x1 ml			24 months		
CONTROLS							
H	CPCOX-002	1x2 ml			24 months		
M	CPCON-002	1x2 ml			24 months		
L	CXCON-002	1x2 ml			24 months		

CRP XL WIDE RANGE

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	CPTXL-B00/ARC	1 x 55,5 ml R1 + 1 x 8 ml R2		140	36 months	serum or plasma	0.607 to 425 (mg/l)
	CPTXL-H00/ARC	2 x 85 ml R1 + 2 x 12 ml R2		459			
	CPTXL-L00/ARC	7 x 85 ml R1 + 7 x 12 ml R2		1607			

CALIBRATORS

KIT 5 LEV	CPREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	CPCOX-002	1x2 ml	24 months
M	CPCON-002	1x2 ml	24 months
L	CXCON-002	1x2 ml	24 months

CYSTATIN C - Colloidal Gold

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	CYCOL-B00/ARC	1 x 27 ml R1 + 1 x 11 ml R2		107	24 months	serum or plasma	0.15 to 9.355 (mg/l)
	CYCOL-H00/ARC	2 x 48 ml R1 + 2 x 20 ml R2		421			
	CYCOL-L00/ARC	10 x 48 ml R1 + 10 x 20 ml R2		2103			

CALIBRATORS

KIT 5 LEV	CYREK-000	5x1 ml	24 months
------------------	-----------	--------	-----------

CONTROLS

H	CYCOX-002	1x2 ml	24 months
M	CYCON-002	1x2 ml	24 months
L	CYCOS-002	1x2 ml	24 months

FERRITIN - Colloidal Gold

REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	FNCOL-B00/ARC	1 x 40 ml R1 + 1 x 8 ml R2		128	12 months	serum	6.344 to 500 (µg/l)
	FNCOL-H00/ARC	2 x 55 ml R1 + 2 x 11 ml R2		378			

CALIBRATORS

KIT 5 LEV	FNRGK-000	5x1 ml	18 months
------------------	-----------	--------	-----------

CONTROLS

M	FNCON-002	1x2 ml	18 months
----------	-----------	--------	-----------

HAPTOGLOBIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	HATUR-B00/ARC	1 x 37,5 ml R1 + 1 x 7,5 ml R2		97	36 months	serum or plasma	0.002 to 3.86 (g/l)
	HATUR-H00/ARC	2 x 75 ml R1 + 2 x 15 ml R2		443			
CALIBRATORS							
KIT 5 LEV	MPREK-000	5x1 ml			24 months		
CONTROLS							
H	MPCOX-002	1x2 ml			24 months		
M	MPCON-002	1x2 ml			24 months		
L	MPCOS-002	1x2 ml			24 months		

IMMUNOGLOBULINS A							
REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	IATUR-B00/ARC	1 x 38 ml R1 + 1 x 8 ml R2		144	36 months	serum	0.02 to 8.0 (g/l)
	IATUR-H00/ARC	2 x 54 ml R1 + 2 x 11 ml R2		425			
CALIBRATORS							
KIT 5 LEV	MPREK-000	5x1 ml			24 months		
CONTROLS							
H	MPCOX-002	1x2 ml			24 months		
M	MPCON-002	1x2 ml			24 months		
L	MPCOS-002	1x2 ml			24 months		

IMMUNOGLOBULINS G							
REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	IGTUR-B00/ARC	1 x 25 ml R1 + 1 x 9,5 ml R2		104	36 months	serum	0.125 to 35 (g/l)
	IGTUR-H00/ARC	2 x 50 ml R1 + 2 x 19 ml R2		461			
CALIBRATORS							
KIT 5 LEV	MPREK-000	5x1 ml			24 months		
CONTROLS							
H	MPCOX-002	1x2 ml			24 months		
M	MPCON-002	1x2 ml			24 months		
L	MPCOS-002	1x2 ml			24 months		

IMMUNOGLOBULINS G (CSF)							
REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	IGLCR-B00/ARC	1 x 18,5 ml R1 + 1 x 5 ml R2		100	36 months	cerebrospinal fluid	3,375 to 249,78 (mg/l)
	IGLCR-H00/ARC	1 x 74 ml R1 + 1 x 17 ml R2		464			
CALIBRATORS							
KIT 5 LEV	IGREK-000	5x1 ml			24 months		
CONTROLS							
H	IGCOX-002	1x2 ml			24 months		
M	IGCON-002	1x2 ml			24 months		
L	IGCOS-002	1x2 ml			24 months		

IMMUNOGLOBULINS M							
REAGENTS	Reference	Content		Tests*	Shelf Life*	Sample**	Analytical Range
		R1	R2				
KIT	IMTUR-B00/ARC	1 x 38 ml R1 + 1 x 8 ml R2		144	36 months	serum	0.038 to 7.0 (g/l)
	IMTUR-H00/ARC	2 x 49,5 ml R1 + 2 x 10,5 ml R2		402			
CALIBRATORS							
KIT 5 LEV	MPREK-000	5x1 ml			24 months		
CONTROLS							
H	MPCOX-002	1x2 ml			24 months		
M	MPCON-002	1x2 ml			24 months		
L	MPCOS-002	1x2 ml			24 months		

LIPOPROTEIN (A)							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	LPTUR2-B00/ARC	1 x 64 ml R1 + 1 x 8 ml R2		191	12 months	serum or plasma	0.0235 to 0.946 (g/l)
	LPTUR2-H00/ARC	2 x 80 ml R1 + 2 x 10 ml R2		503			
CALIBRATORS							
KIT 4 LEV	LPREK-000	4x0,5 ml			18 months		
CONTROLS							
M	LPCON-002	1x2 ml			24 months		

MICROALBUMIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	MATUR-B00/ARC	1 x 22,5 ml R1 + 1 x 6,1 ml R2		100	36 months	urine	1.86 to 400 (mg/l)
	MATUR-H00/ARC	1 x 74 ml R1 + 1 x 17,5 ml R2		361			
	MATUR-L00/ARC	7 x 75 ml R1 + 7 x 15 ml R2		2126			

CALIBRATORS

KIT 5 LEV	MAREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MACOX-002	1x2 ml	24 months
M	MACON-002	1x2 ml	24 months
L	MACOS-002	1x2 ml	24 months

PREALBUMIN							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	PATUR-H00/ARC	1 x 55 ml R1 + 1 x 11 ml R2		155	36 months	serum/plasma	6.53 to 1000 (mg/l)

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

RETINOL BINDING PROTEIN (RBP)							
REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	RBTUR-B00/ARC	1 x 31 ml R1 + 1 x 10 ml R2		138	36 months	serum or plasma	7,81 to 141 (mg/l)
	RBTUR-H00/ARC	2 x 31 ml R1 + 2 x 10 ml R2		277			
	RBTUR-L00/ARC	10 x 31 ml R1 + 10 x 10 ml R2		1383			

CALIBRATORS

KIT 5 LEV	RBREK-000	5x1 ml	18 months
-----------	-----------	--------	-----------

CONTROLS

H	RBCOX-002	1x2 ml	18 months
M	RBCON-002	1x2 ml	18 months
L	RBCOS-002	1x2 ml	18 months

RHEUMATOID FACTOR (RF) - Colloidal Gold

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	RFCOL-B00/ARC	1 x 30 ml R1 + 1 x 20 ml R2		122	18 months	serum and plasma	2.228 to 120 (kIU/l)
	RFCOL-H00/ARC	2 x 60 ml R1 + 2 x 40 ml R2		512			

CALIBRATORS

KIT 5 LEV	RFREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	RFCOX-002	1x2 ml	24 months
M	RFCON-002	1x2 ml	24 months
L	RFCOS-002	1x2 ml	24 months

TRANSFERRIN

REAGENTS	Reference	Content		Tests*	Shelf Life**	Sample***	Analytical Range****
		R1	R2				
KIT	TFTUR-B00/ARC	1 x 37,5 ml R1 + 1 x 10,5 ml R2		108	36 months	serum and plasma	0.076 to 8.0 (g/l)
	TFTUR-H00/ARC	2 x 72 ml R1 + 2 x 20 ml R2		449			
	TFTUR-L00/ARC	7 x 72 ml R1 + 7 x 20 ml R2		1570			

CALIBRATORS

KIT 5 LEV	MPREK-000	5x1 ml	24 months
-----------	-----------	--------	-----------

CONTROLS

H	MPCOX-002	1x2 ml	24 months
M	MPCON-002	1x2 ml	24 months
L	MPCOS-002	1x2 ml	24 months

Please contact us for more information on the following tests:

CALPROTECTIN FECAL

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
----------	-----------	---------	--------	--------------	-----------	----------------------

IMMUNOGLOBULIN E

REAGENTS	Reference	Content	Tests*	Shelf Life**	Sample***	Analytical Range****
----------	-----------	---------	--------	--------------	-----------	----------------------

Your expert in specific proteins



(*) Tests : Amount of real tests. May vary from a analyzer to another one.

() Shelf Life :** Diagam provides to distributors a minimum of **12 months**. The time period written in this document shows the maximum shelf life time after manufacturing.

(*) Sample :** Diagam may provide some applications for other samples, such as CSF, Urin, Veterinary samples etc.... Please contact your distributor or us directly at mail@diagam.com or dial **+32 68 55 14 82**.

(**)**
Analytical range : For general information purposes only. More details can be found in the package insert (instructions for use). Please contact your distributor or us directly at mail@diagam.com or dial **+32 68 55 14 82** in order to get the updated data.



Medical device for in vitro diagnosis.



DIAGAM SA, *Rue du Parc Industriel 40, 7822 Ghislenghien, BELGIUM*



The products fulfill the requirements of European directive 98/79 EC concerning the medical devices of in vitro diagnosis.

© Diagam. 2022. All rights reserved. All product names, registered trademarks, company names in this document remain the property of their respective owners. For general information purposes only. More details can be found in the package insert (instructions for use) and labelling. Not for the US market except for research use only.